




16th International congress on antiphospholipid antibodies task force report on antiphospholipid syndrome laboratory diagnostics and trends

Lupus
2023, Vol. 32(14) 1625–1636
© The Author(s) 2023



Article reuse guidelines:
sagepub.com/journals-permissions
DOI: 10.1177/09612033231211820
journals.sagepub.com/home/lup



Tatsuya Atsumi^{1,*}, Cecilia B Chighizola^{2,*}, Yuichiro Fujieda^{1,*} , Ian Mackie^{3,*}, Massimo Radin^{4,*} , Robert Roubey^{5,*} and Maria Laura Bertolaccini^{6,*} 

Abstract

Classification criteria for antiphospholipid syndrome (APS) require IgG or IgM isotypes of the anticardiolipin (aCL) antibodies, anti- β 2 glycoprotein I (anti- β 2GPI) antibodies, and/or the lupus anticoagulant (LA) to satisfy the laboratory disease definition. Over the past 20 years, non-criteria antiphospholipid antibodies (aPL) directed to other proteins of the coagulation cascade (i.e. prothrombin and/or phosphatidylserine–prothrombin complex) or to some domains of β 2GPI have been proposed. This task force concentrated and reviewed the literature on data including aPS/PT, antibodies to domain 4/5 of β 2GPI and the newly described antibodies to protein/HLA-DR complex. In addition, we discussed testing of LA in the ‘new’ oral anticoagulants’ era and the value of triple positivity in the risk assessment of aPL. The conclusions were presented at a special session during the 16th International Congress on aPL, Manchester, UK, September 2019.

Keywords

Phosphatidylserine-dependent antiprothrombin antibodies, antibodies to domains of β 2-glycoprotein-I, anti- β 2Glycoprotein-I antibodies

Date received: 13 March 2023; accepted: 17 October 2023

Introduction

Classification criteria for antiphospholipid syndrome (APS) require IgG or IgM isotypes of the anticardiolipin (aCL) antibodies, anti- β 2-glycoprotein I (anti- β 2GPI) antibodies, and/or the lupus anticoagulant (LA) to satisfy the laboratory criterion for disease definition.¹ Over the past years several ‘non-criteria’ antiphospholipid antibodies (aPL), directed to proteins of the coagulation cascade (i.e., prothrombin and/or phosphatidylserine–prothrombin complex) or to specific domains of β 2GPI have been focused upon.^{2,3}

The task force on antiphospholipid syndrome laboratory diagnostics and trends focused and reviewed the literature on data including aPS/PT, antibodies to domain 4/5 of β 2GPI and the newly described antibodies to protein/HLA-DR complex. In addition, we discussed testing of LA in the ‘new’ oral anticoagulants’ era and the value of triple positivity in the risk assessment of aPL. The conclusions were presented at a special session during the 16th International

¹Department of Rheumatology, Endocrinology and Nephrology, Faculty of Medicine, Hokkaido University, Sapporo, Japan

²Unit of Pediatric Rheumatology, ASST G. Pini - CTO, Department of Clinical Sciences and Community Health, University of Milan, Milan, Italy.

³Department of Haematology, Haemostasis Research Unit, University College London, London, UK

⁴Department of Clinical and Biological Sciences, and SCU Nephrology and Dialysis, S. Giovanni Bosco Hospital, Center of Research of Immunopathology and Rare Diseases - Coordinating Center of Piemonte and Valle d'Aosta Network for Rare Diseases, Turin, Italy

⁵Division of Rheumatology, Allergy & Immunology, Thurston Arthritis Research Center, University of North Carolina, Chapel Hill, NC, USA

⁶Academic Department of Vascular Surgery, School of Cardiovascular and Metabolic Medicine & Sciences, King's College London, London, UK

*All authors have equally contributed to this report

Corresponding author:

Maria Laura Bertolaccini, Department of Vascular Surgery, School of Cardiovascular and Metabolic Medicine & Sciences, King's College London, 1st Floor North Wing, St Thomas' Hospital, London SE1 7EH, UK. Email laura.bertolaccini@kcl.ac.uk

Congress on aPL, Manchester, UK, September 2019 and updated previous to this publication.

Phosphatidylserine-dependent Antiprothrombin Antibodies (aPS/PT)

Many reports show the clinical utility of phosphatidylserine-dependent antiprothrombin antibodies (aPS/PT) assay in the diagnosis of APS. The inclusion of aPS/PT antibodies as a laboratory criterion of APS has been previously considered by our task force and deemed unwarranted because of poor standardization of the available assays and because reproducibility of the strong correlations between aPS/PT and APS manifestations needed confirmation in larger studies.⁴

Anti-phosphatidylserine/prothrombin (aPS/PT) antibodies are a low affinity heterogeneous class of antibodies directed against a complex of negatively charged phospholipid, other than cardiolipin and prothrombin (PT). The methodology for their detection has improved over the years and currently aPS/PT antibodies are identified by enzyme immunosorbent assays using prothrombin in complex with phosphatidylserine in the presence of calcium. Calcium ion aids the binding of prothrombin to phosphatidylserine inducing major conformational changes to the prothrombin structure which, in turn, exposes cryptic or neo-epitopes that act as target for aPS/PT.⁵ The site of binding on the prothrombin molecule is still under research and it may be possible that different types of antibodies can recognize different sites of the prothrombin molecule.⁶ So far, it has been established that aPS/PT and aPT represent two different types of antibodies that can be present concurrently in some cases.⁷ Prothrombin 1 and fragment 1 and fragment 1+2 have been reported as potential antigens recognized by antiprothrombin antibodies, suggesting that the dominant epitopes are likely to be located near the phospholipid-binding site of the prothrombin molecule.⁵ A recent study by Chinnaraj and colleagues used prothrombin mutants and identified 2 subpopulations of aPS/PT, namely type I and type II, which engage fragment 1 of prothrombin at different epitopes.⁸

Testing for aPS/PT antibodies has been proposed as an additional tool to be considered when investigating a patient suspected of having APS, particularly in the absence of routine aPL positivity,^{2,9} or as a part of risk assessment strategies.¹⁰ aPS/PT represent in fact a stronger risk factor for thrombosis, both arterial and/or venous, than aPT¹¹ and in combination with LA and anti- β 2GPI offer the best diagnostic indication of APS.¹² This is supported by a recent study by Pengo and colleagues, that demonstrated that the search for aPS/PT antibodies along anti- β 2GPI antibodies, in patients positive for LA, might be useful to identify two distinct subgroups of patients at different risk of thromboembolic events.¹³

In addition, patients with triple positivity for LA, anti- β 2GPI and aPS/PT have been shown to be at a higher risk of

developing thromboembolic events, risk even higher than that seen for the 'classical' aCL, anti- β 2GPI, and LA triple positivity.^{12,14}

An early systematic review evaluated papers from 1988 to 2013 and assessed the correlation between aPT and aPS/PT antibodies and the risk of thrombosis.¹¹ Among 10 studies on aPS/PT, comprehensive of 1775 patients and 628 controls, eight of them confirmed the association with thrombosis, but only seven compared the ORs. When performed, multivariate analysis sustained aPS/PT correlation with thrombosis and venous events. The correlation was confirmed in a more recent systematic review¹⁵ that analyzed the studies on aPS/PT and their correlation with clinical manifestations of APS from 2012 to 2019. Briefly, the patient population included 1219 patients classified as APS according to Sidney criteria,¹ 285 patients with isolated persistently positive aPL and 1397 patients with a clinical suspicion of APS.

Twelve studies, including 1888 patients, analyzed the association between aPS/PT antibodies and thrombosis, observing a statistically significant association between aPS/PT IgG/IgM positivity and thrombotic events (mean OR 6.8 [95% CI 3.18–16.4], $p < 0.05$), confirmed when analyzing aPS/PT IgG (mean OR 6.7 [95% CI 3.04–21.6], $p < 0.05$) and aPS/PT IgM (mean OR 4.35 [95% CI 1.54–17.77], $p < 0.05$) separately. Seven studies, including 1388 patients, evaluated the association between aPS/PT antibodies and pregnancy morbidity. When pooled together, a statistically significant association between any pregnancy morbidity and aPS/PT IgG/IgM positivity (mean OR 10.6 [95% CI 3.54–35.38], $p < 0.05$), particularly aPS/PT IgG positivity (mean OR 6.7 [95% CI 3.04–21.6], $p < 0.05$) was found.

Overall, the current available data highlights the strong association between aPS/PT and the clinical manifestations of APS. With the available level of evidence, aPS/PT testing can be considered as a robust test applicable in the investigation of patients suspected of APS, also beyond the research settings.

Antibodies against domains of β 2GPI: a focus on domain 4/5

β 2 glycoprotein I (β 2GPI) provides the main antigen targeted by aPL, the diagnostic serum biomarkers and pathogenic effectors of APS. β 2GPI is believed to exert a relevant biological function, as suggested by the evolutionary conserved structure of this molecule. Although β 2GPI-deficient mice are apparently healthy, recent findings suggest that β 2GPI may serve as a bridge between the innate immune system and the coagulation cascade.¹⁶ Indeed, β 2GPI can engage lipopolysaccharide (LPS), being potentially deputed to remove endotoxin from the circulation.¹⁷ In addition, β 2GPI can

affect the activation of complement cascade: on one hand, it enhances the degradation of C3 by factor I, on the other it activates the lectin pathway by binding directly to the mannose-binding lectin. The latter interaction is mediated by the glycan content of β 2GPI, which is approximately 20%.^{18,19} In turn, the complement lectin pathway can up-regulate coagulation factors thus promoting the activation of coagulation.²⁰ Furthermore, β 2GPI directly modulates coagulation by exerting both a procoagulant and an anticoagulant action. The latter includes the prevention of platelet aggregation induced by ADP and von Willebrand factor, the inhibition of thrombin, factor Xa and tissue activator of plasminogen. However, procoagulant mechanisms prevail: β 2GPI inhibit procoagulant protein C, displace anticoagulant Annexin A5 and prevent the formation of thrombomodulin/thrombin complexes.¹⁶ Even the opposite interaction exists. Plasmin can inactivate β 2GPI by cleaving the molecule in a motif shared with plasminogen (named kringle domain); interestingly, such inactivation can be further augmented by heparin.¹⁶

β 2GPI is composed of five domains (D): D1 to D4 share a homologous structure which consists in 60 amino acids also found in complement control proteins, while D5 is aberrant. It contains an allosteric disulphide bond, which easily undergoes redox modifications, and a loop of lysine residues, which – due to its positively charge – is deputed to interaction with anionic phospholipids (PL), including cardiolipin (CL), as well as with coagulation and complement factors.¹⁶ More than 90% of circulating β 2GPI adopts a circular conformation, with D1 interacting with D5. β 2GPI opens to a J-shaped conformation upon binding to CL and LPS or following changes in pH and oxidative state. This is extremely relevant, as the opening of the molecule leads to the exposure of cryptic epitopes.²¹

From early days in APS research, efforts have focused on the characterization of the reactivity of anti- β 2GPI antibodies from APS patients against the different portions of the β 2GPI molecule. The first pioneer studies envisaged linear peptides, the more recent ones employ conformational epitopes and more sophisticated techniques. As detailed in Table 1, available studies documented a polyreactivity of anti- β 2GPI antibodies from APS patients against different domains of the molecule. However, it is now widely acknowledged that a cryptic and conformation-dependent structure in the N-terminal D1 of the molecule provides the most relevant epitope involved in β 2GPI/anti- β 2GPI antibody binding.^{22–31} Consistently, antibodies against D1 have been extensively documented to exert a pathogenic role in mediating both thrombotic and obstetric complications of the syndrome.^{28,32–34} Anti-D1 antibodies are the prevalent autoantibody subsets among APS patients, being detected across the different cohorts in 40%–88% of subjects with thrombotic APS and 17%–84% of women with obstetric complications.³⁵ Such wide range in the

Table 1. Reactivity against the different domains of β 2 glycoprotein I of antiphospholipid antibodies from patients with antiphospholipid syndrome.

Author, year (REF)	D1	D2	D3	D4	D5
Hunt, 1994 ⁸⁴					X
Wang, 1995 ⁸⁵				X	
Igarashi, 1996 ²²	X	X	X	X	
Yang, 1998 ⁸⁶					X
George, 1998 ⁸⁷				X	
Iverson, 1998 ²³	X				
Blank, 1999 ⁸⁸			X	X	
Reddel, 2000 ²⁴	X				
McNeeley, 2001 ²⁵	X				
Iverson, 2002 ²⁶	X				
Shoenfeld, 2003 ²⁷	X	X	X	X	X
Kawahara, 2005 ⁸⁹				X	
De Laat, 2005 ²⁸	X				
Ioannou, 2007 ²⁹	X				
De Laat, 2011 ³⁰	X				
De Moerloose, 2017 ³¹	X				
Serrano, 2019 ⁹⁰			X	X	X

D: domain.

positivity rates of anti-D1 antibodies might be ascribed to the composition of study populations and to the detection methods.

Besides D1, other studies have focused on the characterization of the reactivity of antibodies against D4 and D5 of β 2GPI.^{36–39} The available studies used a research ELISA kit from INOVA Diagnostics (San Diego, CA, USA), which employs recombinant D4 and D5 bound to the polystyrene microwell plate under conditions preserving the native state. These studies are concordant in describing a low prevalence of anti-D4/5 antibodies (below 30%) in patients with APS (Table 2). Furthermore, the reports reject any association between anti-D4/5 antibodies and clinical manifestations of APS, including vascular events and obstetric manifestations (Table 2). Notably, among patients with aPL positivity associated with other systemic autoimmune rheumatic disease (SARD), the positivity rate of anti-D4/5 antibodies is similar (ranging between 16.7% and 35.3%).^{36,37} Likewise, in these studies, anti-D4/5 antibody reactivity is not associated with thrombotic complications.³⁶ Interestingly, subjects with high-risk aPL profile (namely, two or three positive aPL tests) present lower anti-D4/5 titres and a lower anti-D4/5 positivity rate when compared to those with single aPL positivity.^{38,39} Consistently, patients with anti-D4/5 antibodies tend to display a single positivity for anti- β 2GPI antibodies with negative aCL and LA.^{39–43} Anti-D4/5 antibodies do not react in the β 2GPI-dependent aCL assay, with a magnitude effect of binding inhibition dependent on the β 2GPI concentration, the hypothesis being that anti-D4/5 antibodies do not react against

Table 2. Antibodies against domains 4 and 5 of $\beta 2$ glycoprotein I and clinical manifestations of antiphospholipid syndrome.

Author, year ^{REF}	Isotype	N	Study population	Thrombosis		Pregnancy morbidity		Controls	
				Positivity rate	Association	Positivity rate	Association	N	Positivity rate
Despierres, 2014 ³⁶	IgA	59	SLE	35.3%	No association	N.R.	N.R.	N.R.	N.R.
Andreoli, 2015 ³⁷	IgG	87	anti- $\beta 2$ GPI + PAPS	19.6%	No association	16.1%	No association	30	40%
Pengo, 2015 ³⁸	IgG	65	anti- $\beta 2$ GPI + subjects	N.R.	No association	N.R.	N.R.	40	N.R.
Chighizola, 2018 ³⁹	IgG	108	anti- $\beta 2$ GPI + PAPS	29.1%	No association	29.5%	No association	27	37%

N: number; N.R.: not reported; SLE: systemic lupus erythematosus; APS: antiphospholipid syndrome; PAPS: primary anti-phospholipid syndrome.

$\beta 2$ GPI complexed with CL due to the spatial proximity of the two binding sites.^{41–43}

The above discussed data clearly suggests that anti-D4/5 antibodies could be useful in assessing risk in individuals with anti- $\beta 2$ GPI antibodies with data suggesting that their presence imply a lower risk. The other side of the coin relates to the so-called asymptomatic aPL carriers or subjects with clinical manifestations not related to aPL positivity. Indeed, in these populations, anti- $\beta 2$ GPI antibodies preferentially recognize D4 and D5, as documented in many studies. As detailed in Table 2, aPL carriers (in other words, individuals who are persistently positive for aPL but do not develop any manifestation of the syndrome) present a higher positivity rates and higher titres of anti-D4/5 compared to APS patients.^{20–22} Antibodies purified from subjects with isolated positivity for anti-D4/5 antibodies have been recently shown to react against D5 by means of recombinant domains in the immunoenzymatic assay.⁴³ Other investigators characterized the domain reactivity of anti- $\beta 2$ GPI antibodies from subjects with non-APS condition, localizing the epitope alternately in D4 or D5. Using domain-deleted mutants, Iverson and colleagues observed that anti- $\beta 2$ GPI IgA from 29 patients with atherosclerotic syndrome preferentially target D4.²¹ In another study, anti- $\beta 2$ GPI antibodies of IgG and/or IgM isotypes were detected in 39% of 176 lepromatous patients. The authors performed in vitro experiments using $\beta 2$ GPI deleted mutants and 8C3, a monoclonal antibody that binds to D1, and concluded that anti- $\beta 2$ GPI antibodies in subjects with leprosy react against D5.⁴² Paediatric subjects are an additional population of interest, since anti- $\beta 2$ GPI antibodies are frequently detected in healthy children. It is possible that the *de novo* production of these autoantibodies follows the ubiquitous environmental diffusion of $\beta 2$ GPI. In a cohort of 93 children with different allergic diseases, a high frequency of anti- $\beta 2$ GPI IgG was found in those with atopic dermatitis (42%). The observation that proteolytic cleavage of PL-binding site in the C-terminal loop in D5 abolished antibody binding to

$\beta 2$ GPI suggested that the antigenic site resides in close vicinity.⁴¹ A subsequent study found an even higher positivity rate of anti- $\beta 2$ GPI in 33 children with atopic dermatitis (54.5%), with a polarization towards D4/5 reactivity (33%). Anti-D4/5 antibodies have been detected in 21 of 57 one-year-old healthy children born to mothers with SARD (36.8%), whereas anti-D1 antibodies tested positive in nine children (15.7%).⁴⁴

The observation that anti-D4/5 antibodies do not react with $\beta 2$ GPI complexed with CL might imply that such autoantibody subpopulation does not react against cell-bound $\beta 2$ GPI and thus are not able to elicit a pathogenic potential.²⁷ As a matter of fact, when the pathogenic potential of antibodies purified from subjects with isolated positivity for anti-D4/5 antibodies in mediating vascular occlusion was explored in LPS-treated Wistar rats, these antibodies failed to promote thrombosis. This was evaluated as vessel occlusion in mesenteric vessels using intravital microscopy and was significantly different from what was observed with anti-D1 IgG.⁴³ In vitro experiments showed that anti-D5 antibodies bind significantly less to D5 after incubation with $\beta 2$ GPI at increasing concentrations, suggesting its interaction with $\beta 2$ GPI in soluble form. These results led to formulate the hypothesis that anti-D5 might antagonize the procoagulant activity of anti-D1 antibodies by competing for $\beta 2$ GPI binding. As approximately 15% of APS patients display positivity for both anti-D1 and anti-D5 IgG,³⁷ these observations might exert important clinical implications. Indeed, due to the interaction of anti-D5 IgG with soluble $\beta 2$ GPI, it is tempting to speculate that anti-D5 IgG might prevent $\beta 2$ GPI binding to target cells, thus antagonizing the procoagulant activity of anti-D1 antibodies. The magnitude of this competitive effect might depend on the relative anti-D1 and anti-D5 antibody levels. This intriguing hypothesis has been partially explored in patients by calculating the ratio between anti-D1 and anti-D4/5 antibodies (the so-called ‘anti-D1:anti-D4/5 ratio’). In two studies, the relevance of this tool has been evaluated in anti-

β 2GPI IgG positive patients subgrouped upon clinical presentation: thrombotic and/or obstetric APS, SARD, and asymptomatic aPL positivity. Interestingly, employing two research ELISAs, it was evinced that an anti-D1:anti-D4/5 ratio above 1.5 is predictive of systemic autoimmunity (APS and SARD).³⁷ Other authors detected anti-D1 and anti-D4/5 antibodies by a novel line immunoassay; an anti-D1:anti-D4/5 ratio above 4.6 could significantly distinguish APS patients from subjects with SARD.⁴⁵ In another study including exclusively patients with primary APS, anti-D1:anti-D4/5 ratio above 2.1 yielded an odds ratio of 2.7 for APS diagnosis, with a sensitivity of 62% and a specificity of 63%.³⁹

Overall, current evidence rejects any association between D4/5 epitope specificity and both thrombotic and obstetric manifestations of APS. Consistently, available data negate a role for anti-D4/5 antibodies in the pathological processes of APS. However, the hypothesis of a protective role of anti-D4/5 antibodies warrants further confirmation. Available evidence comes from a single group, and some data are controversial. When rats were treated with three different autoantibody preparations (anti-D1 antibodies, anti-D4/5 antibodies, IgG from normal healthy subjects) no difference in vascular deposition of β 2GPI could be evidenced at immunofluorescence of mesenteric tissue, a finding that conflicts with the hypothesis that anti-D4/5, by binding to β 2GPI, might prevent its endothelial localization. At present, anti-D4/5 antibodies might be regarded as a second-line test to be reserved to anti- β 2GPI antibody positive subjects to refine the process of stratification of the hazard of future clinical events. Future works should aim at clarifying the clinical significance of anti- β 2GPI antibodies not reacting against D1 and D4/5, which is the case in approximately one fourth of APS subjects.

Antibodies against protein/HLA-DR complex

Anti- β 2GPI/Human Leukocyte Antigen (HLA)-DR antibodies has been reported as a new aPL test based on a novel mechanism of antigen presentation.⁴⁶ The specific APS susceptible HLA-DR alleles influence the cell-surface expression of β 2GPI/HLA-DR complexes. β 2GPI/HLA-DR complexes are considered as major target antigens for autoantibodies in patients with APS.

Genome-wide analysis has confirmed that major histocompatibility complex (MHC) class II loci showed strong association with susceptibility to many autoimmune diseases.⁴⁷ MHC class II molecules are primarily expressed by professional antigen presentation cells, such as dendritic cells, macrophages, B cells. MHC class II α - and β -chains assemble in the endoplasmic reticulum (ER) and form a complex with the invariant chain.⁴⁸ The invariant chain-

MHC class II heterotrimer is transported through the Golgi to the MHC class II compartment, either directly and/or via the plasma membrane. Consequently, peptides derived from the endocytic compartment are presented as antigen to CD4+ T cells. In general, misfolded proteins localized in the ER are degraded promptly inside cells. The novel theory, however, is that misfolded proteins are rescued from protein degradation in the ER2 and transported to the cell surface by MHC class II molecules as 'neo-self-antigens'. Jiang et al.⁴⁶ found that misfolded proteins can associate with MHC class II molecules in the ER instead of the invariant chain, depending on allelic polymorphisms in MHC class II genes. In addition, 'neo-self-antigens' activate the antigen-specific B cells,⁴⁹ suggesting that these complexes of misfolded proteins with MHC class II molecules are considered major targets of autoantibodies in multiple autoimmune diseases.⁵⁰

aPL bind to some part of β 2GPI/HLA-DR complexes which are transported to the cell surface without degradation. aPL may recognise an epitope on the β 2GPI/HLA-DR complexes as well as phospholipid binding β 2GPI. Anti- β 2GPI/HLA-DR antibodies are detected by a 'cell-based assay' using 293T cells co-transfected with β 2GPI and HLA-DR.⁵¹ Since almost half of the APS patients were single positive for anti- β 2GPI/HLA-DR antibodies, while negative for anti- β 2GPI and aCL, β 2GPI/HLA-DR complexes may expose unique epitopes that are not present on plate-bound β 2GPI or β 2GPI/CL complexes. Interestingly, the expression levels of β 2GPI on cell surface were different between HLA-DR7 and HLA-DR8. The comprehensive analysis of HLA-DR alleles revealed that, in addition to HLA-DR7, HLA-DR4 was transporting high levels of β 2GPI to the cell surface and, conversely, several other HLA-DR alleles were transporting very little β 2GPI. HLA-DR7 or HLA-DR4 has been reported as APS susceptibility alleles,⁵² suggesting that specific HLA-DR alleles might be associated with expressing 'neo-self-antigens' as well as susceptibility allele to APS. Unlike T cell receptors, HLA-DR allele itself does not affect the autoantibody binding to β 2GPI/HLA-DR complexes. Sera from 120 patients with APS were examined and over 80% of APS patients showed anti- β 2GPI/HLA-DR7 antibodies.⁵¹ Considering these results, using β 2GPI/HLA-DR7 complex as an assay antigen is enough to analyze the prevalence of anti- β 2GPI/HLA-DR antibodies in patients with APS.

Anti- β 2GPI/HLA-DR antibodies are also associated with obstetric complication in patients with APS. In fact, β 2GPI/HLA-DR complexes were found in uterine decidual tissues from patients with APS.⁵¹ Recently, by conducting prospective, multicentre, cross-sectional study, Tanimura et al.⁵³ reported that women with unexplained recurrent pregnancy loss (RPL) had anti- β 2GPI/HLA-DR antibodies. RPL was defined as the loss of ≥ 2 pregnancies and β 2GPI/HLA-DR7 complex was used as an antigen to detect anti- β 2GPI/HLA-DR antibodies in these women. Of the

227 women with RPL, 45 (19.8%) tested positive for aPL such as LA, aCL, or anti- β 2GPI and 52 (22.9%) tested positive for anti- β 2GPI/HLA-DR antibodies. Of these, 35 (67.3%) were single-positive. Among the women with unexplained RPL, 24 (19.8%) were positive for anti- β 2GPI/HLA-DR antibodies. Of the 112 women who did not meet criteria for APS, 21 (18.8%) were positive for anti- β 2GPI/HLA-DR antibodies. These results suggest that anti- β 2GPI/HLA-DR antibodies may be useful in classifying women with unexplained RPL into obstetric APS.

Anti- β 2GPI/HLA-DR antibodies have a potential as disease-specific antibodies in APS. However, the presence of β 2GPI/HLA-DR complexes in the human body remains to be investigated. β 2GPI is produced not only by hepatocytes, but also by endothelial cells and placental villous tissue.^{54,55} Non-immune cells strongly express MHC class II molecules in response to stimulation from interferon γ (IFN γ), but do not express costimulatory molecules required for the induction of T cell responses, such as CD80 or CD86. Stimulated non-immune cells can activate B cell directly. β 2GPI/HLA-DR complexes can be expressed on the cell surface during infection or inflammation and could stimulate B cells to produce autoantibodies. This aberrant HLA expression in APS would explain the pathogenicity of autoimmune disease, leading to novel therapeutic approaches.

Lupus anticoagulant and new anticoagulants

The ISTH survey on LA testing⁵⁶ showed that 70.3% of participants thought that LA testing should not be performed in patients receiving direct oral anticoagulants (DOAC). Among the participants, 17% suggested that testing could be performed in the trough period, 11% after pre-treatment of the sample with commercial DOAC adsorbent or antidote preparations, while 2.7% felt that LA testing may be undertaken in some circumstances in patients on DOACs during the peak concentration period. This indicates that there is considerable uncertainty about what action to take in patients receiving DOACs.

There are numerous publications demonstrating that DOACs interfere with LA testing, causing false positive or false negative results^{57–66} and there are particular problems with direct thrombin inhibitors. There are a variety of potential ways around the problems:

1. Wait until the patient has stopped DOAC or switch to a different anticoagulant
2. Collect sample in the trough period
3. Use a LA test that is not affected by DOAC
4. Neutralise/remove DOAC

Stopping the anticoagulant or waiting until the patient has finished their course of anticoagulant treatment are not usually practical. LA detection could influence the drug choice, intensity of anticoagulation and duration of treatment. Stopping anticoagulation or switching to a different anticoagulant (bridging therapy) could put the patient at risk of adverse events. Low molecular weight heparin is less likely to affect LA tests (most commercial reagents contain heparin neutralisers effective up to about 1 IU/mL). There is also a potential for patient confusion about any altered dosing or change of medication type if the anticoagulant is changed to facilitate blood tests.

Testing for LA during the trough period would appear to be a suitable approach, since anticoagulant levels would be too low to have much influence on the tests. However, false positive dilute Russel's viper venom time (DRVVT) results have been reported in trough even at very low (<50ng/mL) rivaroxaban levels.^{61,65} It is also difficult to be certain that the patient is in the trough period, as there may have been errors or variability in the time of taking medication. One answer might be to assay the DOAC, but the sensitivity varies between assays and the minimum concentration of each DOAC type that affects LA tests is not known (and different LA tests vary in their sensitivity to DOACs). External quality assurance studies have shown a small degree of inter-laboratory variability in DOAC assays (10%–12%), but the greatest differences between reagents were observed for rivaroxaban, especially at concentrations above 100ng/mL, where the cv was 10%–15% and there was a difference of about 20ng/mL between methods.⁶⁷ Many laboratories have set up their assays to measure samples at peak DOAC levels and the lower detection limit of their method may be only 50ng/mL. A recent publication⁶⁸ studied 60 venous thromboembolism patients receiving DOAC (30 rivaroxaban, 30 apixaban) during the trough period and found mean levels of 23ng/mL (range <18–68) and 42ng/mL (19–99) respectively. 93% of rivaroxaban and 40% apixaban patients had a false positive DRVVT due to a marked effect of the DOAC on the screening part of the test. Similarly, 40 and 30% of patients had false positive silica clotting time and 17 and 20% showed false positive APTT ratio comparing LA sensitive and insensitive reagents. In addition, LA results did not necessarily correlate with the drug level, even when patients appeared to be truly in the trough period.

Dabigatran appears to affect most LA tests, since it is a direct thrombin inhibitor and all well recognised LA diagnostic tests employ a coagulation end point method. For direct factor Xa inhibitors, the Taipan/Ecarin venom clotting time ratio appears to be suitable,^{58,59,69,70} but the reagents are currently only available from one supplier (Diagnostic Reagents Ltd, Thame, UK) and the method has not been validated for apixaban. The method requires further standardisation and wider assessment before it can be adopted into routine clinical practice (ISTH guidelines).

Drugs have been developed that are therapeutic antidotes and neutralise DOACS. Adexanet alpha has been used *in vitro* but does not fully neutralise rivaroxaban and therefore does not normalise the clotting time,⁷¹ but *in vitro* addition of idarucizumab (a therapeutic humanized antibody fragment) neutralized dabigatran.⁷² The regular use of these therapeutic agents *in vitro* in clinical diagnostic laboratories is likely to be difficult and cost prohibitive. It would be difficult to know whether neutralization/removal had occurred and the relevant DOAC would probably have to be assayed after sample treatment. It is also unknown whether residual neutralizing agent can interfere in the DOAC assay.

A more successful approach appears to be with the use of activated charcoal products to remove DOAC, although this is not without problems. There are two commercial branded products, DOAC-STOP (Haematex Research) and DOAC-Remove (5-Diagnostics) that have been widely studied.^{72–76} DOAC-STOP was shown to reduce DOAC levels to <30ng/mL in 60%–100% of VTE patient samples and known APS patient samples remained LA positive after treatment. However, some studies have reported 5%–15% false LA positives, which may be partly due to incomplete DOAC removal, particularly if drug levels are >350ng/mL. DOAC-Remove is generally effective, but one study showed that DOAC was incompletely removed in 4% of samples. An unbranded activated carbon product (Norit Carbomix, Kela Pharma) has also been shown to be effective, but a few false positive LA results were observed. The authors suggested that the normal plasma used for mixing studies and cut-off determination should also be treated.⁷⁷

In summary, the various activated charcoal products are not completely problem free. Some of the studies have used DOAC spiked samples rather than performing *ex vivo* studies and even then, the number of patient samples has been limited. DOAC plasma levels were not always assayed after sample treatment to ensure adequate removal. Further studies are underway with larger amounts of removal agent which may resolve some of the problems. Unless the DOAC removal methods are sanctioned by LA test reagent/analyser manufacturers, the technique may be deemed as manipulation and create regulatory and accreditation problems, meaning that full, local validation exercises will have to be undertaken before adoption into routine clinical laboratory practice.

Assessing ‘Triple Positivity’

Antiphospholipid antibody (aPL) ‘triple positivity’ is defined as positivity in the three most used clinical aPL laboratory tests: (1) aCL IgG or IgM, medium or high level, (2) anti-β2GPI IgG or IgM, medium or high level, (3) LA performed in accordance with international standards. It is well-established and broadly accepted that aPL ‘triple

positivity’ is a high-risk phenotype for the clinical manifestations of the antiphospholipid syndrome (APS).

To understand the nature of ‘triple positivity’ it is important to recognize that it is defined based on clinical tests and not on antibody specificity. Historically this has been a source of confusion. Figure 1 schematically shows the relationship among major clinical laboratory tests, the antibodies detected in each test, and the clinical associations of these antibodies. Typical aCL assays can detect antibodies of various specificities. Although purified cardiolipin is the intended antigen, bovine serum is typically used as the blocking agent/sample diluent and cardiolipin-binding proteins in bovine serum (notably, β2GPI) can also serve as target antigens. Thus, aCL assays can detect antibodies binding directly to cardiolipin (seen in patients with syphilis, certain other infections, and occasionally in normal individuals), antibodies to bovine β2GPI (associated with APS), and possibly antibodies to other cardiolipin-binding proteins in bovine serum (without known clinical associations). Anti-β2GPI assays detect antibodies to purified human β2GPI (associated with APS). Research has shown that antibodies responsible for lupus anticoagulant activity (in the setting of APS) can be directed against β2GPI or prothrombin. Available data also suggest that antibodies to D1 of β2GPI have LA activity²⁸ while the value of antibodies to D4/5 is still being investigated. In contrast to the relatively lower affinity anti-prothrombin antibodies in APS patients, very high affinity anti-prothrombin antibodies, which are quite rare, can be associated with hypoprothrombinemia and a bleeding disorder. Anti-prothrombin and aPS/PT immunoassays which detect anti-prothrombin antibodies are emerging.

aCL and anti-β2GPI tests detect largely overlapping subsets of antibodies to β2GPI. In APS patients the concordance of these two tests is high. By definition, ‘triple positive’ patients are concordant. There are several reasons why there may be discordant results in these two tests. Firstly, antibodies recognizing cardiolipin in aCL assays are not detected in anti-β2GPI tests. Secondly, although anti-β2GPI in most APS patients are reactive with both human and bovine β2GPI there are rare patients with antibodies specific for human, but not bovine, β2GPI. Thirdly, some discordance may be due to the epitope specificity of anti-β2GPI antibodies and the differential display of epitopes in the two assays.

LA tests are coagulation-based assays, not immunoassays. These tests detect certain antibodies based on the antibodies’ interference with a clotting reaction. A considerable body of research demonstrates that anti-β2GPI or anti-prothrombin antibodies exert LA activity if these antibodies can form high avidity cross-linked phospholipid-bound antigen/antibody complexes in *in vitro* coagulation assays thereby inhibiting these coagulation reactions by decreasing the available anionic phospholipid surface. There are two major factors that appear to determine whether anti-β2GPI and/or anti-prothrombin antibodies have LA activity. The first is




Figure 1. Relationship of aPL clinical tests, antibodies detected, and associated clinical manifestations.

Table 3. ‘aPL Load’, LA positivity and ‘triple positivity’.

A	LA	N	Average ‘aPL Load’	
	Negative	44	63.1 ± 78.4	
	Positive	281	134.1 ± 123.0	

B	‘aPL Load’ (quintiles)	N	LA+	% LA+
	Q1 (0–25)	64	46	72%
	Q2 (26–50)	66	52	79%
	Q3 (51–117)	65	59	91%
	Q4 (118–220)	66	64	97%
	Q5 (221–665)	64	60	94%
	All subjects	325	281	86%

C	‘Triple positivity’	N	Average ‘aPL Load’	
	Not ‘triple positive’	245	72.2 ± 64.3	
	‘Triple positive’	80	284.6 ± 110.3	

$p < 0.0001$, Mann-Whitney test.
 X^2 for trend = 21.23, $p < 0.0001$.
 $p < 0.0001$, Mann-Whitney test.

epitope specificity. For example, it is likely that antibodies to D1 of β2GPI can cross-link phospholipid-bound β2GPI since D1 is exposed when the protein is bound to a membrane. In contrast, antibodies to D4/5 of β2GPI may not be able to cross-link bound β2GPI because these domains are not exposed or less exposed on the membrane-bound molecule. The second key factor, and one that is sometimes overlooked, is antibody concentration. Simply put, in contrast to the exquisite analytical sensitivity of ELISAs and other immunoassays, antibody inhibition of an in vitro coagulation reaction requires high antibody concentration. This has been clearly demonstrated with both monoclonal and polyclonal antibodies.^{78–81} For monoclonal antibodies, LA activity

requires 100- to 1000-fold higher concentrations compared to detection in ELISAs. Using polyclonal IgG purified from 2 APS patients with high levels of IgG anti-β2GPI, we observed that LA activity required 5 to 10 times the IgG concentration that was strongly positive in an anti-β2GPI ELISA. These observations suggest that LA activity and, therefore, ‘triple positivity’ are proxies for high antibody titer. Additionally, positive LA tests can reflect the combined effect of all relevant antibodies in a given specimen (aCL, anti-β2GPI, anti-prothrombin, all isotypes of each specificity) and do not indicate a single antibody profile. A positive LA can be due to anti-β2GPI in some patients, anti-prothrombin antibodies in other patients, and perhaps to both specificities in still others. In a cohort of 254 patients from the Antiphospholipid Syndrome Collaborative Registry, LA activity was associated with higher levels of IgG aCL and IgM aCL, and this association was more pronounced when IgG and IgM levels were combined.⁸¹

More recently the association of antibody levels, LA activity, and ‘triple positivity’ was evaluated in a group of patients enrolled in the APS ACTION Registry Database. Complete core lab data (aCL, anti-β2GPI, and interpretable LA results) were available from 325 patients. Most of these patients (281, 86%) were LA positive; 80 patients (25%) were ‘triple positive’. An ‘aPL load’ was crudely calculated as the sum of four test results (IgG aCL, IgM aCL, IgG anti-β2GPI, IgM anti-β2GPI). Note that this calculation does not include aPS/PT, anti-PT, IgA aCL, or IgA anti-β2GPI. As shown in Table 3, ‘aPL Load’ was strongly associated with LA positivity and ‘triple positivity’.

It has long been recognized that the risks for thrombosis, pregnancy morbidity and mortality, and other clinical manifestations of APS are associated with higher aPL titers. Higher risk is also associated with LA positivity and with

‘triple positivity’. It is likely that LA activity requires both epitope specificities that allow for the cross-linking of membrane-bound antigens, for example, antibodies to D1 of β 2GPI, as well as high levels of such antibodies. Importantly, the factors that contribute to LA activity in vitro appear to play an important role in the pathogenicity of aPL in vivo. For example, dimerized D5 of β 2GPI (which mimics membrane-bound β 2GPI cross-linked by antibodies to D1) has been shown to upregulate monocyte procoagulant activity mimicking the effects of anti- β 2GPI.⁸²

Scoring strategies that combine and weight aPL tests, such as the antiphospholipid score⁸³ and the Global APS Score¹⁰ (the latter of which also considers other risk factors) may prove useful in applying the concept of aPL ‘load’ to clinical evaluation.

Declaration of conflicting interests

The author(s) declared no potential conflicts of interest with respect to the research, authorship, and/or publication of this article.

Funding

The author(s) disclosed receipt of the following financial support for the research, authorship, and/or publication of this article: MLB is supported by the King’s BHF Centre for Award Excellence RE18/2/34213.

ORCID iDs

Yuichiro Fujieda  <https://orcid.org/0000-0003-4705-341X>

Massimo Radin  <https://orcid.org/0000-0003-1941-2606>

Maria Laura Bertolaccini  <https://orcid.org/0000-0002-8623-139X>

References

- Miyakis S, Lockshin MD, Atsumi T, et al. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost* 2006; 4(2): 295–306.
- Zohoury N, Bertolaccini ML, Rodriguez-Garcia JL, et al. Closing the serological gap in the antiphospholipid syndrome: the value of “Non-criteria” antiphospholipid antibodies. *J Rheumatol* 2017; 44: 1597–1602. DOI: [10.3899/jrheum.170044](https://doi.org/10.3899/jrheum.170044).
- Bertolaccini ML and Sanna G. The clinical relevance of noncriteria antiphospholipid antibodies. *Semin Thromb Haemost* 2017; 44: 453–457. DOI: [10.1055/s-0037-1601328](https://doi.org/10.1055/s-0037-1601328).
- Bertolaccini ML, Amengual O, Artim-Eser B, et al. Clinical and prognostic significance of non-criteria antiphospholipid antibody tests. In: Erkan D and Lockshin MD (eds) *Antiphospholipid Syndrome: Current Research Highlights and Clinical Insights*. Cham: Springer International Publishing, 2017, pp. 171–187.
- Bertolaccini ML. Antibodies to prothrombin. *Lupus* 2012; 21(7): 729–731. (In eng). DOI: [10.1177/0961203312443299](https://doi.org/10.1177/0961203312443299).
- Pozzi N and Di Cera E. Prothrombin structure: unanticipated features and opportunities. *Expert Rev Proteomics* 2014; 11(6): 653–655. DOI: [10.1586/14789450.2014.971763](https://doi.org/10.1586/14789450.2014.971763).
- Bertolaccini ML, Atsumi T, Koike T, et al. Antiprothrombin antibodies detected in two different assay systems. Prevalence and clinical significance in systemic lupus erythematosus. *Thromb Haemost* 2005; 93(2): 289–297.
- Chinnaraj M, Planer W, Pengo V, et al. Discovery and characterization of 2 novel subpopulations of aPS/PT antibodies in patients at high risk of thrombosis. *Blood Advances* 2019; 3(11): 1738–1749. DOI: [10.1182/bloodadvances.2019030932](https://doi.org/10.1182/bloodadvances.2019030932).
- Rodriguez-Garcia V, Ioannou Y, Fernandez-Nebro A, et al. Examining the prevalence of non-criteria anti-phospholipid antibodies in patients with anti-phospholipid syndrome: a systematic review. *Rheumatology* 2015; 54(11): 2042–2050. DOI: [10.1093/rheumatology/kev226](https://doi.org/10.1093/rheumatology/kev226).
- Sciascia S, Sanna G, Murru V, et al. GAPSS: the Global Anti-Phospholipid Syndrome Score. *Rheumatology* 2013; 52(8): 1397–1403. DOI: [10.1093/rheumatology/kes388](https://doi.org/10.1093/rheumatology/kes388).
- Sciascia S, Sanna G, Murru V, et al. Anti-prothrombin (aPT) and anti-phosphatidylserine/prothrombin (aPS/PT) antibodies and the risk of thrombosis in the antiphospholipid syndrome. A systematic review. *Thromb Haemost* 2014; 111(2): 354–364. DOI: [10.1160/TH13-06-0509](https://doi.org/10.1160/TH13-06-0509).
- Sciascia S, Murru V, Sanna G, et al. Clinical accuracy for diagnosis of antiphospholipid syndrome in systemic lupus erythematosus: evaluation of 23 possible combinations of antiphospholipid antibody specificities. *J Thromb Haemost* 2012; 10(12): 2512–2518. (In eng). DOI: [10.1111/jth.12014](https://doi.org/10.1111/jth.12014).
- Pengo V, Del Ross T, Ruffatti A, et al. Lupus anticoagulant identifies two distinct groups of patients with different antibody patterns. *Thromb Res* 2018; 172: 172–178. DOI: [10.1016/j.thromres.2018.11.003](https://doi.org/10.1016/j.thromres.2018.11.003).
- Tonello M, Mattia E, Favaro M, et al. IgG phosphatidylserine/prothrombin antibodies as a risk factor of thrombosis in antiphospholipid antibody carriers. *Thromb Res* 2019; 177: 157–160. DOI: [10.1016/j.thromres.2019.03.006](https://doi.org/10.1016/j.thromres.2019.03.006).
- Radin M, Foddai SG, Cecchi I, et al. Antiphosphatidylserine/prothrombin antibodies: an update on their association with clinical manifestations of antiphospholipid syndrome. *Thromb Haemost* 2020; 120(4): 592–598. DOI: [10.1055/s-0040-1705115](https://doi.org/10.1055/s-0040-1705115).
- McDonnell T, Wincup C, Buchholz I, et al. The role of beta-2-glycoprotein I in health and disease associating structure with function: more than just APS. *Blood Rev* 2020; 39: 100610. DOI: [10.1016/j.blre.2019.100610](https://doi.org/10.1016/j.blre.2019.100610).
- Agar C, de Groot PG, Morgelin M, et al. beta(2)-glycoprotein I: a novel component of innate immunity. *Blood* 2011; 117(25): 6939–6947. DOI: [10.1182/blood-2010-12-325951](https://doi.org/10.1182/blood-2010-12-325951).
- Gambino R, Ruiu G, Pagano G, et al. Study of the glycosylation of apolipoprotein H. *Chem Phys Lipids* 1999; 103(1-2): 161–174. DOI: [10.1016/s0009-3084\(99\)00108-5](https://doi.org/10.1016/s0009-3084(99)00108-5).

19. Ruben E, Planer W, Chinnaraj M, et al. The J-elongated conformation of beta2-glycoprotein I predominates in solution: implications for our understanding of antiphospholipid syndrome. *J Biol Chem* 2020; 295(31): 10794–10806. DOI: [10.1074/jbc.RA120.013939](https://doi.org/10.1074/jbc.RA120.013939).
20. Durigutto P, Macor P, Pozzi N, et al. Complement activation and thrombin generation by MBL bound to beta2-glycoprotein I. *J Immunol* 2020; 205(5): 1385–1392. DOI: [10.4049/jimmunol.2000570](https://doi.org/10.4049/jimmunol.2000570).
21. Chighizola CB, Gerosa M and Meroni PL. New tests to detect antiphospholipid antibodies: anti-domain I beta-2-glycoprotein-I antibodies. *Curr Rheumatol Rep* 2014; 16(2): 402. DOI: [10.1007/s11926-013-0402-7](https://doi.org/10.1007/s11926-013-0402-7).
22. Igarashi M, Matsuura E, Igarashi Y, et al. Human beta2-glycoprotein I as an anticardiolipin cofactor determined using mutants expressed by a baculovirus system. *Blood* 1996; 87(8): 3262–3270.
23. Iverson GM, Victoria EJ and Marquis DM. Anti-beta2 glycoprotein I (beta2GPI) autoantibodies recognize an epitope on the first domain of beta2GPI. *Proc Natl Acad Sci U S A* 1998; 95(26): 15542–15546.
24. Reddel SW, Wang YX, Sheng YH, et al. Epitope studies with anti-beta 2-glycoprotein I antibodies from autoantibody and immunized sources. *J Autoimmun* 2000; 15(2): 91–96. DOI: [10.1006/jaut.2000.0427](https://doi.org/10.1006/jaut.2000.0427).
25. McNeeley PA, Dlott JS, Furie RA, et al. Beta2-glycoprotein I-dependent anticardiolipin antibodies preferentially bind the amino terminal domain of beta2-glycoprotein I. *Thromb Haemost* 2001; 86(2): 590–595.
26. Iverson GM, Reddel S, Victoria EJ, et al. Use of single point mutations in domain I of beta 2-glycoprotein I to determine fine antigenic specificity of antiphospholipid autoantibodies. *J Immunol* 2002; 169(12): 7097–7103.
27. Shoenfeld Y, Krause I, Kvapil F, et al. Prevalence and clinical correlations of antibodies against six beta2-glycoprotein-I-related peptides in the antiphospholipid syndrome. *J Clin Immunol* 2003; 23(5): 377–383.
28. de Laat B, Derksen RH, Urbanus RT, et al. IgG antibodies that recognize epitope Gly40-Arg43 in domain I of beta 2-glycoprotein I cause LAC, and their presence correlates strongly with thrombosis. *Blood* 2005; 105(4): 1540–1545. DOI: [10.1182/blood-2004-09-3387](https://doi.org/10.1182/blood-2004-09-3387) [pii].
29. Ioannou Y, Pericleous C, Giles I, et al. Binding of antiphospholipid antibodies to discontinuous epitopes on domain I of human beta(2)-glycoprotein I: mutation studies including residues R39 to R43. *Arthritis Rheum* 2007; 56(1): 280–290. (In eng). DOI: [10.1002/art.22306](https://doi.org/10.1002/art.22306).
30. de Laat B, van Berkel M, Urbanus RT, et al. Immune responses against domain I of beta(2)-glycoprotein I are driven by conformational changes: domain I of beta(2)-glycoprotein I harbors a cryptic immunogenic epitope. *Arthritis Rheum* 2011; 63(12): 3960–3968. DOI: [10.1002/art.30633](https://doi.org/10.1002/art.30633).
31. de Moerloose P, Fickentscher C, Boehlen F, et al. Patient-derived anti-beta2GPI antibodies recognize a peptide motif pattern and not a specific sequence of residues. *Haematologica* 2017; 102(8): 1324–1332. DOI: [10.3324/haematol.2017.170381](https://doi.org/10.3324/haematol.2017.170381).
32. Agostinis C, Durigutto P, Sblattero D, et al. A non complement-fixing antibody to beta2 glycoprotein I as a novel therapy to control abortions and thrombosis in antiphospholipid syndrome. *Blood*. 2014; 123(22): 3478–3487. doi:[10.1182/blood-2013-11-537704](https://doi.org/10.1182/blood-2013-11-537704), 24642748.
33. Zandman-Goddard G, Pierangeli SS, Gertel S, et al. Tolerogenic dendritic cells specific for beta2-glycoprotein-I Domain-I, attenuate experimental antiphospholipid syndrome. *J Autoimmun* 2014; 54: 72–80. DOI: [10.1016/j.jaut.2014.06.001](https://doi.org/10.1016/j.jaut.2014.06.001).
34. Pericleous C, Ruiz-Limon P, Romay-Penabad Z, et al. Proof-of-concept study demonstrating the pathogenicity of affinity-purified IgG antibodies directed to domain I of beta2-glycoprotein I in a mouse model of anti-phospholipid antibody-induced thrombosis. *Rheumatology* 2015; 54(4): 722–727. DOI: [10.1093/rheumatology/keu360](https://doi.org/10.1093/rheumatology/keu360).
35. de Jesus GR, Benson AE, Chighizola CB, et al. 16th International congress on antiphospholipid antibodies task force report on obstetric antiphospholipid syndrome. *Lupus* 2020; 29(12): 1601–1615. DOI: [10.1177/0961203320954520](https://doi.org/10.1177/0961203320954520).
36. Despierres L, Beziane A, Kaplanski G, et al. Contribution of anti-beta2glycoprotein I IgA antibodies to the diagnosis of anti-phospholipid syndrome: potential interest of target domains to discriminate thrombotic and non-thrombotic patients. *Rheumatology* 2014; 53(7): 1215–1218. DOI: [10.1093/rheumatology/keu003](https://doi.org/10.1093/rheumatology/keu003).
37. Andreoli L, Chighizola CB, Nalli C, et al. Clinical Characterization of Antiphospholipid Syndrome by Detection of IgG Antibodies Against beta2 -Glycoprotein I Domain 1 and Domain 4/5: Ratio of Anti-Domain 1 to Anti-Domain 4/5 As a Useful New Biomarker for Antiphospholipid Syndrome. *Arthritis Rheumatol* 2015; 67(8): 2196–2204. DOI: [10.1002/art.39187](https://doi.org/10.1002/art.39187).
38. Pengo V, Ruffatti A, Tonello M, et al. Antibodies to Domain 4/5 (Dm4/5) of beta2-Glycoprotein 1 (beta2GPI) in different antiphospholipid (aPL) antibody profiles. *Thromb Res* 2015; 136(1): 161–163. DOI: [10.1016/j.thromres.2015.04.031](https://doi.org/10.1016/j.thromres.2015.04.031).
39. Chighizola CB, Pregnotato F, Andreoli L, et al. Beyond thrombosis: Anti-beta2GPI domain 1 antibodies identify late pregnancy morbidity in anti-phospholipid syndrome. *J Autoimmun* 2018; 90: 76–83. DOI: [10.1016/j.jaut.2018.02.002](https://doi.org/10.1016/j.jaut.2018.02.002).
40. Pengo V, Ruffatti A, Tonello M, et al. Antiphospholipid syndrome: antibodies to Domain 1 of beta2-glycoprotein 1 correctly classify patients at risk. *J Thromb Haemost* 2015; 13(5): 782–787. DOI: [10.1111/jth.12865](https://doi.org/10.1111/jth.12865).
41. Ambrozic A, Avicin T, Ichikawa K, et al. Anti-beta(2)-glycoprotein I antibodies in children with atopic dermatitis. *Int Immunol* 2002; 14(7): 823–830.
42. Arvieux J, Renaudineau Y, Mane I, et al. Distinguishing features of anti-beta2 glycoprotein I antibodies between patients with leprosy and the antiphospholipid syndrome. *Thromb Haemost* 2002; 87(4): 599–605.

43. Durigutto P, Grossi C, Borghi MO, et al. New insight into antiphospholipid syndrome: antibodies to beta2glycoprotein I-domain 5 fail to induce thrombi in rats. *Haematologica* 2019; 104(4): 819–826. DOI: [10.3324/haematol.2018.198119](https://doi.org/10.3324/haematol.2018.198119).
44. Andreoli L, Nalli C, Motta M, et al. Anti-beta(2)-glycoprotein I IgG antibodies from 1-year-old healthy children born to mothers with systemic autoimmune diseases preferentially target domain 4/5: might it be the reason for their 'innocent' profile? *Ann Rheum Dis* 2011; 70(2): 380–383. DOI: [10.1136/ard.2010.137281](https://doi.org/10.1136/ard.2010.137281).
45. Nalli C, Somma V, Andreoli L, et al. Anti-phospholipid IgG antibodies detected by line immunoassay differentiate patients with anti-phospholipid syndrome and other autoimmune diseases. *Auto Immun Highlights* 2018; 9(1): 6. DOI: [10.1007/s13317-018-0106-0](https://doi.org/10.1007/s13317-018-0106-0).
46. Jiang Y, Arase N, Kohyama M, et al. Transport of misfolded endoplasmic reticulum proteins to the cell surface by MHC class II molecules. *Int Immunol* 2013; 25(4): 235–246. DOI: [10.1093/intimm/dxs155](https://doi.org/10.1093/intimm/dxs155).
47. Matzaraki V, Kumar V, Wijmenga C, et al. The MHC locus and genetic susceptibility to autoimmune and infectious diseases. *Genome Biol* 2017; 18(1): 76. DOI: [10.1186/s13059-017-1207-1](https://doi.org/10.1186/s13059-017-1207-1).
48. Neefjes J, Jongasma ML, Paul P, et al. Towards a systems understanding of MHC class I and MHC class II antigen presentation. *Nat Rev Immunol* 2011; 11(12): 823–836. DOI: [10.1038/nri3084](https://doi.org/10.1038/nri3084).
49. Arase H. Rheumatoid rescue of misfolded cellular proteins by MHC Class II Molecules: a new hypothesis for autoimmune diseases. *Adv Immunol* 2016; 129: 1–23. DOI: [10.1016/bs.ai.2015.09.005](https://doi.org/10.1016/bs.ai.2015.09.005).
50. Jin H, Arase N, Hirayasu K, et al. Autoantibodies to IgG/HLA class II complexes are associated with rheumatoid arthritis susceptibility. *Proc Natl Acad Sci U S A* 2014; 111(10): 3787–3792. DOI: [10.1073/pnas.1401105111](https://doi.org/10.1073/pnas.1401105111).
51. Tanimura K, Jin H, Suenaga T, et al. beta2-Glycoprotein I/HLA class II complexes are novel autoantigens in antiphospholipid syndrome. *Blood* 2015; 125(18): 2835–2844. DOI: [10.1182/blood-2014-08-593624](https://doi.org/10.1182/blood-2014-08-593624).
52. Sebastiani GD, Minisola G and Galeazzi M. HLA class II alleles and genetic predisposition to the antiphospholipid syndrome. *Autoimmun Rev* 2003; 2(6): 387–394. DOI: [10.1016/s1568-9972\(03\)00068-5](https://doi.org/10.1016/s1568-9972(03)00068-5).
53. Tanimura K, Saito S, Nakatsuka M, et al. The beta 2 -glycoprotein I/HLA-DR complex as a major autoantibody target in obstetric antiphospholipid syndrome. *Arthritis Rheumatol* 2020; 72(11): 1882–1891. DOI: [10.1002/art.41410](https://doi.org/10.1002/art.41410).
54. Caronti B, Calderaro C, Alessandri C, et al. Beta2-glycoprotein I (beta2-GPI) mRNA is expressed by several cell types involved in anti-phospholipid syndrome-related tissue damage. *Clin Exp Immunol* 1999; 115(1): 214–219. DOI: [10.1046/j.1365-2249.1999.00770.x](https://doi.org/10.1046/j.1365-2249.1999.00770.x).
55. Chamley LW, Allen JL and Johnson PM. Synthesis of beta2 glycoprotein 1 by the human placenta. *Placenta* 1997; 18(5-6): 403–410. DOI: [10.1016/s0143-4004\(97\)80040-9](https://doi.org/10.1016/s0143-4004(97)80040-9).
56. Cohen H, Mackie IJ, Devreese KMJ, et al. Standardization Committee for Lupus Anticoagulant/Antiphospholipid A. Clinical and laboratory practice for lupus anticoagulant testing: An International Society of Thrombosis and Haemostasis Scientific and Standardization Committee survey. *J Thromb Haemost* 2019; 17(10): 1715–1732. DOI: [10.1111/jth.14560](https://doi.org/10.1111/jth.14560).
57. Merriman E, Kaplan Z, Butler J, et al. Rivaroxaban and false positive lupus anticoagulant testing. *Thromb Haemost* 2011; 105(2): 385–386. DOI: [10.1160/TH10-08-0511](https://doi.org/10.1160/TH10-08-0511).
58. van Os GM, de Laat B, Kamphuisen PW, et al. Detection of lupus anticoagulant in the presence of rivaroxaban using Taipan snake venom time. *J Thromb Haemost* 2011; 9(8): 1657–1659. DOI: [10.1111/j.1538-7836.2011.04395.x](https://doi.org/10.1111/j.1538-7836.2011.04395.x).
59. Arachchillage DR, Mackie IJ, Efthymiou M, et al. Interactions between rivaroxaban and antiphospholipid antibodies in thrombotic antiphospholipid syndrome. *J Thromb Haemost* 2015; 13(7): 1264–1273. DOI: [10.1111/jth.12917](https://doi.org/10.1111/jth.12917).
60. Bonar R, Favaloro EJ, Mohammed S, et al. The effect of the direct factor Xa inhibitors apixaban and rivaroxaban on haemostasis tests: a comprehensive assessment using in vitro and ex vivo samples. *Pathology* 2016; 48(1): 60–71. DOI: [10.1016/j.pathol.2015.11.025](https://doi.org/10.1016/j.pathol.2015.11.025).
61. Ratzinger F, Lang M, Belik S, et al. Lupus-anticoagulant testing at NOAC trough levels. *Thromb Haemost* 2016; 116(2): 235–240. DOI: [10.1160/TH16-02-0081](https://doi.org/10.1160/TH16-02-0081).
62. Hoxha A, Banzato A, Ruffatti A, et al. Detection of lupus anticoagulant in the era of direct oral anticoagulants. *Autoimmun Rev* 2017; 16(2): 173–178. DOI: [10.1016/j.autrev.2016.12.010](https://doi.org/10.1016/j.autrev.2016.12.010).
63. Mani H, Hesse C, Stratmann G, et al. Ex vivo effects of low-dose rivaroxaban on specific coagulation assays and coagulation factor activities in patients under real life conditions. *Thromb Haemost* 2013; 109(1): 127–136. DOI: [10.1160/TH12-04-0228](https://doi.org/10.1160/TH12-04-0228).
64. Martinuzzo ME, Barrera LH, Ma D, et al. Frequent false-positive results of lupus anticoagulant tests in plasmas of patients receiving the new oral anticoagulants and enoxaparin. *Int J Lab Hematol* 2014; 36(2): 144–150. DOI: [10.1111/ijlh.12138](https://doi.org/10.1111/ijlh.12138).
65. Martinuzzo ME, Forastiero R, Duboscq C, et al. False-positive lupus anticoagulant results by DRVVT in the presence of rivaroxaban even at low plasma concentrations. *Int J Lab Hematol* 2018; 40(5): e99–e101. DOI: [10.1111/ijlh.12865](https://doi.org/10.1111/ijlh.12865).
66. Antovic A, Norberg EM, Berndtsson M, et al. Effects of direct oral anticoagulants on lupus anticoagulant assays in a real-life setting. *Thromb Haemost* 2017; 117(9): 1700–1704. DOI: [10.1160/TH17-03-0204](https://doi.org/10.1160/TH17-03-0204).
67. Hollestelle MJ and Meijer P. International external quality assessment for measurements of direct oral anticoagulants: results and recommendations. *Br J Haematol* 2020; 188(3): 460–464. DOI: [10.1111/bjh.16179](https://doi.org/10.1111/bjh.16179).
68. Gay J, Duchemin J, Imarazene M, et al. Lupus anticoagulant diagnosis in patients receiving direct oral FXa inhibitors at trough levels: A real-life study. *Int J Lab Hematol* 2019; 41(6): 738–744. DOI: [10.1111/ijlh.13101](https://doi.org/10.1111/ijlh.13101).

69. Pouplard C, Vayne C, Berthomet C, et al. The Taipan snake venom time can be used to detect lupus anticoagulant in patients treated by rivaroxaban. *Int J Lab Hematol* 2017; 39(3): e60–e63. DOI: [10.1111/ijlh.12611](https://doi.org/10.1111/ijlh.12611).
70. Moore GW, Culhane AP, Maloney JC, et al. Taipan snake venom time coupled with ecarin time enhances lupus anticoagulant detection in nonanticoagulated patients. *Blood Coagul Fibrinolysis* 2016; 27(4): 477–480. DOI: [10.1097/MBC.0000000000000486](https://doi.org/10.1097/MBC.0000000000000486).
71. Favaloro EJ, Gilmore G, Arunachalam S, et al. Neutralising rivaroxaban induced interference in laboratory testing for lupus anticoagulant (LA): A comparative study using DOAC Stop and andexanet alfa. *Thromb Res* 2019; 180: 10–19. DOI: [10.1016/j.thromres.2019.05.013](https://doi.org/10.1016/j.thromres.2019.05.013).
72. Jacquemin M, Toelen J, Feyen L, et al. The adsorption of dabigatran is as efficient as addition of idarucizumab to neutralize the drug in routine coagulation assays. *Int J Lab Hematol* 2018; 40(4): 442–447. DOI: [10.1111/ijlh.12807](https://doi.org/10.1111/ijlh.12807).
73. Exner T, Michalopoulos N, Pearce J, et al. Simple method for removing DOACs from plasma samples. *Thromb Res* 2018; 163: 117–122. DOI: [10.1016/j.thromres.2018.01.047](https://doi.org/10.1016/j.thromres.2018.01.047).
74. Favresse J, Lardinois B, Sabor L, et al. Evaluation of the DOAC-Stop® Procedure to Overcome the Effect of DOACs on Several Thrombophilia Screening Tests. *TH Open* 2018; 2(2): e202–e209.
75. Platton S and Hunt C. Influence of DOAC Stop on coagulation assays in samples from patients on rivaroxaban or apixaban. *Int J Lab Hematol* 2019; 41(2): 227–233. DOI: [10.1111/ijlh.12950](https://doi.org/10.1111/ijlh.12950).
76. Zabczyk M, Kopytek M, Natorska J, et al. The effect of DOAC-Stop on lupus anticoagulant testing in plasma samples of venous thromboembolism patients receiving direct oral anticoagulants. *Clin Chem Lab Med / FESCC* 2019; 57(9): 1374–1381. DOI: [10.1515/cclm-2018-1197](https://doi.org/10.1515/cclm-2018-1197).
77. Frans G, Meeus P and Bailleul E. Resolving DOAC interference on aPTT, PT, and lupus anticoagulant testing by the use of activated carbon. *J Thromb Haemost* 2019; 17(8): 1354–1362. DOI: [10.1111/jth.14488](https://doi.org/10.1111/jth.14488).
78. Jennings I, Mackie I, Arnout J, et al. Lupus anticoagulant testing using plasma spiked with monoclonal antibodies: performance in the UK NEQAS proficiency testing programme. *J Thromb Haemost* 2004; 2(12): 2178–2184. DOI: [10.1111/j.1538-7836.2004.01028.x](https://doi.org/10.1111/j.1538-7836.2004.01028.x).
79. Takeya H, Mori T, Gabazza EC, et al. Anti-beta2-glycoprotein I (beta2GPI) monoclonal antibodies with lupus anticoagulant-like activity enhance the beta2GPI binding to phospholipids. *J Clin Invest* 1997; 99(9): 2260–2268.
80. Roubey RA, Pratt CW, Buyon JP, et al. Lupus anticoagulant activity of autoimmune antiphospholipid antibodies is dependent upon beta 2-glycoprotein I. *J Clin Invest* 1992; 90(3): 1100–1104.
81. Roubey RA. Risky business: the interpretation, use, and abuse of antiphospholipid antibody tests in clinical practice. *Lupus* 2010; 19(4): 440–445. DOI: [10.1177/0961203310361356](https://doi.org/10.1177/0961203310361356).
82. Kolyada A, Barrios DA and Beglova N. Dimerized domain V of beta2-glycoprotein I is sufficient to upregulate procoagulant activity in PMA-treated U937 monocytes and require intact residues in two phospholipid-binding loops. *Antibodies* 2017; 6(2): 8. DOI: [10.3390/antib6020008](https://doi.org/10.3390/antib6020008).
83. Otomo K, Atsumi T, Amengual O, et al. Efficacy of the antiphospholipid score for the diagnosis of antiphospholipid syndrome and its predictive value for thrombotic events. *Arthritis Rheum* 2012; 64(2): 504–512. (Research Support, Non-U.S. Gov't). DOI: [10.1002/art.33340](https://doi.org/10.1002/art.33340).
84. Hunt J and Krilis S. The fifth domain of beta 2-glycoprotein I contains a phospholipid binding site (Cys281-Cys288) and a region recognized by anticardiolipin antibodies. *J Immunol* 1994; 152: 653–659.
85. Wang MX, Kandiah DA, Ichikawa K, et al. Epitope specificity of monoclonal anti-beta 2-glycoprotein I antibodies derived from patients with the antiphospholipid syndrome. *J Immunol* 1995; 155(3): 1629–1636.
86. Yang CD, Chen SL, Shen N, et al. Detection of anti-recombinant beta 2-glycoprotein I and anti-recombinant beta 2-glycoprotein I fifth-domain antibodies in sera from patients with systemic lupus erythematosus. *Rheumatol Int* 1998; 18(1): 5–10. DOI: [10.1007/s002960050046](https://doi.org/10.1007/s002960050046).
87. George J, Gilburd B, Hojnik M, et al. Target recognition of beta2-glycoprotein I (beta2GPI)-dependent anticardiolipin antibodies: evidence for involvement of the fourth domain of beta2GPI in antibody binding. *J Immunol* 1998; 160(8): 3917–3923.
88. Blank M, Waisman A, Mozes E, et al. Characteristics and pathogenic role of anti-beta2-glycoprotein I single-chain Fv domains: induction of experimental antiphospholipid syndrome. *Int Immunol* 1999; 11(12): 1917–1926.
89. Kasahara H, Matsuura E, Kaihara K, et al. Antigenic structures recognized by anti-beta2-glycoprotein I auto-antibodies. *Int Immunol* 2005; 17(12): 1533–1542. DOI: [10.1093/intimm/dxh330](https://doi.org/10.1093/intimm/dxh330).
90. Serrano M, Martinez-Flores JA, Norman GL, et al. The IgA Isotype of anti-beta2 glycoprotein I antibodies recognizes epitopes in domains 3, 4, and 5 that are located in a lateral zone of the molecule (L-Shaped). *Front Immunol* 2019; 10: 1031. DOI: [10.3389/fimmu.2019.01031](https://doi.org/10.3389/fimmu.2019.01031).