

AperTO - Archivio Istituzionale Open Access dell'Università di Torino

Conventional and narrow bore short capillary columns with cyclodextrin derivatives as chiral selectors to speed-up enantioselective gas chromatography and enantioselective gas chromatography-mass spectrometry analyses

This is the author's manuscript

Original Citation:

Availability:

This version is available <http://hdl.handle.net/2318/132061> since

Terms of use:

Open Access

Anyone can freely access the full text of works made available as "Open Access". Works made available under a Creative Commons license can be used according to the terms and conditions of said license. Use of all other works requires consent of the right holder (author or publisher) if not exempted from copyright protection by the applicable law.

(Article begins on next page)



UNIVERSITÀ DEGLI STUDI DI TORINO

This Accepted Author Manuscript (AAM) is copyrighted and published by Elsevier. It is posted here by agreement between Elsevier and the University of Turin. Changes resulting from the publishing process - such as editing, corrections, structural formatting, and other quality control mechanisms - may not be reflected in this version of the text. The definitive version of the text was subsequently published in [*Journal of Chromatography A*, Volume: 1217, Issue: 7, Pages: 1106-1113, date: FEB 12 2010, doi:10.1016/j.chroma.2008.10.013].

You may download, copy and otherwise use the AAM for non-commercial purposes provided that your license is limited by the following restrictions:

- (1) You may use this AAM for non-commercial purposes only under the terms of the CC-BY-NC-ND license.
- (2) The integrity of the work and identification of the author, copyright owner, and publisher must be preserved in any copy.
- (3) You must attribute this AAM in the following format: Creative Commons BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/deed.en>), doi:10.1016/j.chroma.2008.10.013

1 **Conventional and narrow bore short capillary columns with cyclodextrin derivatives as chiral**
2 **selectors to speed-up enantioselective gas chromatography and enantioselective gas**
3 **chromatography mass spectrometry analyses**

4

5 **Carlo Bicchi, Erica Liberto, Cecilia Cagliero, Chiara Cordero, Barbara Sgorbini, Patrizia**
6 **Rubiolo***

7 Dipartimento di Scienza e Tecnologia del Farmaco - Università degli Studi di Torino

8 Via Pietro Giuria, 9 10125 Torino – Italia

9

10 **Address for correspondence:**

11 Prof. Dr. Patrizia Rubiolo

12 Dipartimento di Scienza e Tecnologia del Farmaco, Facoltà di Farmacia,

13 Università degli Studi di Torino, Via Pietro Giuria 9, Turin 10125, Italy;

14 Fax: +39 011 6707687; e-mail: patrizia.rubiolo@unito.it

15

16 Abstract

17 The analysis of complex real-world samples of vegetable origin requires rapid and accurate routine
18 methods, enabling laboratories to increase sample throughput and productivity while reducing
19 analysis costs. This study examines shortening enantioselective GC (ES-GC) analyses following the
20 approaches used in fast-GC. ES-GC separations are conditioned by a weak enantiomer-CD host-
21 guest interaction and the separation is thermodynamically driven and strongly influenced by
22 temperature. As a consequence, fast temperature rates can interfere with enantiomeric
23 discrimination; thus the use of short and/or narrow bore columns is a possible approach to speeding-
24 up ES-GC analyses. The performance of ES-GC with a conventional inner diameter (I.D.) column
25 (25 m length x 0.25 mm I.D., 0.15 μm and 0.25 μm d_f) coated with 30% of 2,3-di-O-ethyl-6-O-tert-
26 butyldimethylsilyl- β -cyclodextrin in PS-086 is compared to those of conventional I.D. short column
27 (5 m length x 0.25 mm I.D., 0.15 μm d_f) and of different length narrow bore columns (1, 2, 5 and 10
28 m long x 0.10 mm I.D., 0.10 μm d_f) in analysing pesticide standard, racemate standards and real-
29 world-samples in the flavour and fragrance field.

30 Short conventional I.D. columns gave shorter analysis time and comparable or lower resolutions
31 with the racemate standards, depending mainly on analyte volatility.

32 Narrow-bore columns were tested under different analysis conditions; they provided shorter
33 analysis time and resolutions comparable to those of conventional I.D. ES columns. The narrow-
34 bore columns offering the most effective compromise between separation efficiency and analysis
35 time are the 5 and 2 m columns; in combination with mass spectrometry as detector, applied to
36 lavender and bergamot essential oil analyses, these reduced analysis time by a factor of at least
37 three while separation of chiral markers remained unaltered.

38

39 **Keywords:** fast enantioselective GC-MS, cyclodextrin derivatives, linear retention indices, mass
40 spectra, flavours, fragrances.

41

42 1. Introduction

43 Cyclodextrin derivatives (CDs) have gained an increasing acceptance as chiral stationary
44 phases (CSPs) in GC separation of enantiomers because of their wide enantioselectivity and ability
45 to separate underivatized enantiomers of different volatilities. When used in enantioselective gas
46 chromatography (ES-GC), CDs are generally diluted in apolar or moderately polar polysiloxanes to
47 obtain highly efficient capillary GC columns, as proposed by Schurig *et al.* [1-3]. The success of a
48 separation technique depends on its suitability for use in the routine analysis of real-world samples,
49 whose requirements include the ability to identify characteristic markers in a sample and short
50 analysis time. In general, a correct enantiomeric excess (EE) or enantiomeric ratio (ER)
51 determination of optically active markers in complex real-world samples requires a two-
52 dimensional approach, since enantioselective analysis can make the resulting chromatographic
53 profile more complex, because enantiomer separation may cause peak doubling, increasing the risk
54 of peak overlap. Two complementary but distinct approaches are those most commonly adopted:
55 the first and most popular one is based on a second dimension *in separation* with heart-cut GC-GC
56 [4-7] or comprehensive two-dimensional GC (GCxGC) [8-11]; the second approach includes a
57 second dimension *in identification* using mass spectrometry (MS) (or very rarely FT-IR) as ES-GC
58 detector. MS is known not to be a selective chiral probe making enantiomer MS spectra
59 indistinguishable. These spectra can therefore successfully be exploited to determine EE or ER
60 correctly in both the total and the extract ion modes although they can only be used to locate the
61 enantiomers in the chromatogram since their linear retention indices (I^T) [12] are necessary to
62 identify them unequivocally [13].

63

64 One of the approaches used to satisfy the increasing demand for analytical controls is to
65 adopt fast chromatographic techniques, which enable high sample throughput and laboratory
66 productivity, while at the same time reducing analysis costs. fast-GC and fast-GC-MS have long
67 been established for routine analysis [14-19]; these methods involve contemporary contributions of

67 several parameters including column length and inner diameter (I.D.), stationary phase film
68 thickness (d_f), temperature programme and linear carrier gas velocity. A limit of ES-GC in routine
69 quality control is the long analysis time since the separation depends on a small difference in the
70 energy of the association between each enantiomer and the CD chiral selector and requires a very
71 high chromatographic efficiency. The ES-GC separation of enantiomers by CD derivatives is
72 known to be based on fast kinetics and entirely governed by thermodynamics [20, 21], and therefore
73 strongly controlled by temperature. In contrast to what happens with fast conventional GC, the
74 temperature rates that can be used to speed up an ES-GC separation with CDs as CSPs are therefore
75 rather low. Fast-ES-GC can thus mainly be achieved by operating on column length, inner diameter
76 and/or flow-rate and in particular with short narrow bore (NB) columns. These columns not only
77 increase analysis speed and analyte detectability due to peak sharpening [22], but also reduce
78 enantiomer elution temperature, resulting in a gain of enantioselectivity that compensates (in full or
79 in part) for the loss of efficiency (N). Short columns were applied in ES-GC with CDs since the
80 beginning of the 1990s enabling enantiomer separations even in a few seconds [23-26]. After an in-
81 depth study of ES-GC, Schurig and Czesla [27] concluded that short conventional 0.25 mm I.D.
82 columns should be used for fast-ES-GC because of their good loadability, integration characteristic,
83 use of conventional instrumentation and lower consumption of carrier gas.

84 However, conventional I.D. short columns can be applied in monodimensional ES-GC only
85 when optically active compounds must be determined in low complexity samples and/or when a
86 low number of enantiomers must be analysed simultaneously. With medium-to-high complexity
87 samples, a highly-efficient separation system combined with single- or multiple-ion monitoring-MS
88 detection is necessary to determine EE and/or ER correctly. This article examines how to optimise
89 the parameters conditioning the speeding up of ES-GC-MS analysis of medium-to-high- and
90 medium-to-low-volatility analytes; it also looks at the application of fast-ES-GC-qMS to routine
91 analysis of real world samples in the field of flavours and fragrances.

92

93 **2. Experimental**

94 *2.1 Samples*

95 Pure standards of racemic linalool, linalyl propionate, γ -lactones (C6, C7, C8, C11, C12, C14,
96 C15), δ -octalactone, α -HCH, trichlorfon, *trans*-chlordane, and heptachlor, together with the
97 components of the laboratory-prepared chiral test mixture consisting of limonene, 2-octanol,
98 camphor, isobornyl acetate, linalyl acetate, 2-methyl-(3Z)-hexenyl butyrate, menthol,
99 hydroxycitronellal, γ -decalactone and δ -decalactone racemates [28], were from the collection of
100 standards in the authors' laboratory or, if unavailable there, were obtained from Sigma-Aldrich
101 (Milan, Italy). All standard compounds were dissolved in cyclohexane at a concentration of 100
102 ppm each. Solvents were all HPLC grade from Riedel-de Haen (Seelze, Germany). Lavender
103 (*Lavandula angustifolia* P. Mill.) and bergamot (*Citrus bergamia* Risso et Poiteau) essential oils
104 (e.o.), obtained by hydrodistillation following the method described in the European Pharmacopoeia
105 (6th edition) [29], were diluted 1:200 in cyclohexane before analysis.

106

107 *2.2 Instrumental set-up*

108 A Shimadzu QP2010 GC-MS system provided with Shimadzu GCMS Solution 2.51 software
109 (Shimadzu, Milan, Italy) was used.

110 Columns: GC analyses were carried out on the set of columns reported in table 1 and coated with
111 2,3-di-O-ethyl-6-O-*tert*-butyldimethylsilyl- β -CD (2,3DE6TBDMS- β -CD) [30] as CSP diluted at 30
112 % in PS086. All columns were from MEGA (Legnano, Italy).

113 GC-MS conditions: injection mode: split; split ratio: 1: 20 for conventional I.D. column and split
114 ratio 1:100 for narrow bore (NB) columns; injection volume: 1 μ l. Temperatures: injector: 220°C,
115 transfer line: 230°C; ion source: 200°C; carrier gas: He. The MS was operated in electron impact
116 ionization mode (EI) at 70 eV, the scan rate was 1111 u/s and a mass range of 35–350 m/z suitable
117 to cover the full fragmentation pattern of most e.o. components was chosen. Flavour and fragrance

118 compounds were analysed starting from 50°C while pesticides from 100°C. Temperature
119 programmes and flow rates are listed in table 1. Resolution was calculated with the following
120 equation: $R = 1.18 (t_{R(2)} - t_{R(1)}) / (w_{0.5(1)} + w_{0.5(2)})$
121
122

123 3. Results and discussion

124 This section is divided into two main parts: the first comprises an evaluation of the effect of
125 temperature and flow rates, stationary phase film thickness, inner diameter and column length on
126 speeding up and performance of ES-GC or ES-GC-MS with CD derivatives as CSPs, and the second
127 concerns the application of these techniques to the analysis of real-world samples in the e.o. field.
128 Some pesticide standard racemates were also analysed to extend the range of structures and
129 characteristics of the analytes investigated. If not specified otherwise, all results are compared to
130 those obtained with a 25 m long, 0.25 mm I.D., 0.25 μm d_f taken as reference column (hereafter
131 called “reference column”).
132

133 3.1 Effect of column length, inner diameter, film thickness and chromatographic conditions on the 134 speed of an enantioselective GC analysis.

135 One of the conditions for a correct ES-GC analysis is that the enantiomers of the chiral
136 compounds investigated are baseline separated must to afford a correct EE (or ER) determination.
137 The excess of efficiency with several chiral compounds of the columns coated with the last
138 generation of CD derivatives can be exploited to shorten analysis time by acting on temperature and
139 flow rates, stationary phase film thickness, inner diameter and column length. This approach has
140 here first been applied to the analysis of a group of standards of racemates with different polarity
141 and volatility (table 2) with a set of columns coated with 2,3-DE-6TBDMS- β -CD (30% in PS-086)
142 as CSP [30] that includes the reference column, two conventional 0.25 mm I.D. columns 5 m long x
143 0.15 μm d_f and 25 m x 0.15 μm and four 0.10 mm narrow bore (NB) columns 1, 2, 5, and 10 m long
144 and 0.10 μm d_f . The influence of the different columns and analysis conditions on the separation
145 performance of the chromatographic systems was evaluated by the separation measure, S , [31],
146 calculated by analysing a chiral test mixture used routinely in the authors’ laboratory [28]. The
147 adopted analysis conditions were chosen to achieve the simultaneous baseline separation of as many
148 chiral components as possible for a given real-world sample. Analyses were therefore operated in
149 temperature programmed mode and not in isothermal mode, as it is usual for theoretical studies
150 where the separation of each chiral analyte must be optimized.

151 The first experiments concerned the influence of length and film thickness on the performance of
152 conventional I.D. columns. All experiments were run at a heating rate of 2 °C/min, without
153 optimising elution temperatures. Table 2 reports retention time (t_R) of the second eluted enantiomer,
154 elution temperature of both enantiomers (T_{e1-2}), resolution (R) and % analysis time reduction
155 calculated on the second eluting enantiomer of the chiral standard compounds investigated, when
156 analysed with conventional I.D. columns of different lengths (25 and 5 m) and film thicknesses
157 (0.25 and 0.15 μm). Figure 1 reports the ES-GC separation of α -HCH with the three 0.25 I.D.
158 columns in question. The results show that a thinner CSP film thickness produces a different
159 reduction of analysis time, ranging from 37% for limonene to 13% for γ -pentadecalactone, with
160 increased resolution for all analytes. These results can be explained by the decreased elution
161 temperature, –between 13°C –and 20°C lower, which improves the separation capability of the
162 chromatographic system.

163 Under the same conditions, the effect of column length is even more evident. Compared to the
164 25 m - 0.15 μm d_f , analysis times with the 5 m - 0.15 μm d_f column decreased from about 69% for
165 limonene to about 26% for γ -pentadecalactone, with a loss of resolution ranging from about 65%
166 for γ -heptalactone to about 10% for heptachlor. As expected, the decreased efficiency due to
167 column shortening had a bigger influence on the resolution of the better-separated enantiomers.
168 Moreover, the percent reduction of analysis time decreased with analyte volatility, as shown with

169 the homologous series of γ -lactones, where it ranges from about -60% with γ -hexalactone to -26%
170 with γ -pentadecalactone. The combined effect of thinner CSP film thickness and column shortening
171 produced a reduction in analysis time ranging from about 80% for limonene to 35% for γ -
172 pentadecalactone, with a decrease in elution temperature ranging between 28°C and 55°C. In real
173 terms, the limonene analysis time decreased from about 17 to 3.5 min, from about 69 to 44 min for
174 γ -pentadecalactone and from about 28 to 6 min for α -HCH; for the last two analytes, the resolution
175 remained almost unaltered.

176 Column length strongly influences the separation capability of the ES-GC system: the
177 separation measure S calculated on the chiral test mixture for the reference column decreased by
178 more than 50% with the 5 m x 0.15 μm d_f column (from 942 to 491) while S was comparable (i.e.
179 972) when film thickness was reduced to 0.15 μm with the 25 m column (table 3).

180 This loss of separation capability (and efficiency) of the chromatographic system can affect
181 routine analysis of real-world samples, but can be overcome by increasing efficiency with columns
182 with narrower I.D. [14-19] and/or selectivity by applying MS detection in extract ion mode. With
183 the same phase ratios, narrow-bore columns are known to give higher efficiency compared to
184 conventional I.D. columns, thus enabling shorter columns to be used. The most widely used I.D. is
185 0.10 mm, because the resulting columns enable reasonable split ratios and/or dilutions that avoid
186 column overloading, while reducing analysis time by a factor varying from 4 to 10. The same value
187 of enantiomer resolution of 1.5 was adopted as the limit of analysis acceptance. Two racemates with
188 different volatility and polarity (linalool and γ -pentadecalactone) and the chiral test mixture [28]
189 were analysed with four 0.10 mm I.D. columns of different lengths (10, 5, 2, 1 m) but coated with
190 the same CSP film thickness (0.10 μm of 30% 2,3-DE-6TBDMS- β -CD in PS-086) and under
191 different temperature and flow rates. All analysis were started at the same temperature (50°C). The
192 ES-GC conditions are reported in table 1. Figure 2 shows the ES-GC patterns of linalool analysed
193 with the four NB columns under investigation, at 1 mL flow rate at 10°C/min. Figure 3 shows the
194 relationships of resolution (A) and retention time (B) *versus* flow and temperature rates of γ -
195 pentadecalactone. In the diagram A, the limit of accepted resolution (1.5) is shown as a full line.
196 The highest resolution with linalool was obtained in all cases with the 5 m narrow bore column,
197 followed by that 2 m long. The low values for the 10 m column are not surprising and may be due
198 to the long linalool residence time in the column, which increases elution temperature. Linalool
199 resolution in all conditions was above 1.5; thus the 1 m column could be used successfully, with a
200 reduction in the analysis time of about 95% compared to the reference column (24.3 min *vs.* 1.1
201 min) and of 80% compared to the 10 m NB column (5.3 min *vs.* 1.1 min) when quicker ES-GC
202 conditions are applied. Similar considerations can be made for γ -pentadecalactone, for which the 5
203 and the 2 m NB columns provided base-line separation up to 5°C/min in contrast to the 10 and 1m
204 NB columns, with which 1.5 resolution could only be achieved at 2°C/min. With the 5 and 2 m NB
205 columns, the reduction in analysis time under the quickest conditions (1.0 mL/min – 5°C/min) was
206 about 65% (69.1 min *vs.* 24.7 min) compared to the reference column and 56% compared to the 10
207 m NB column base-line separation (56.8 min *vs.* 24.7 min).

208 The most effective NB columns are therefore the 2 m and the 5 m long columns, the latter
209 probably offering the best compromise between analyte residence time and system separation
210 efficiency. These considerations were confirmed when the separation capability of the four NB
211 columns was measured on the chiral test mixture through the separation measure S . Table 3 clearly
212 shows that S for the 5 m NB column is higher than that of 10 m column and more than double the
213 value for the corresponding conventional I.D. 5m column, even at 10°C/min. The low separation
214 capability of the 10 and 1 m NB columns is due to different factors: the efficiency of the 1 m NB
215 column is too low and is not compensated for by temperature decrease, while for the 10 m NB
216 column, efficiency is high but the high elution temperature interferes with effective chiral
217 discrimination

218 The reliability of fast ES-GC-MS with NB columns *versus* EE and ER determination was also
219 evaluated, by analysing linalool with the four NB columns investigated under all applied analysis

220 conditions. Table 4 shows % areas of each enantiomer on the total areas of the *R*- and *S*-linalool
221 enantiomers and related SD and RSD% compared to those obtained with the reference column
222

223 3.2 Fast ES-GC/MS analysis of real-world samples with narrow bore columns

224 The second part of this section deals with the analysis of real-world samples. Two
225 commercially important essential oils of medium complexity, i.e. lavender and bergamot, were
226 analysed by ES-GC-MS.

227 Lavender essential oil contains several optically active components: α - and β -pinene, camphene, β -
228 phellandrene, limonene, 1-octen-3-ol, camphor, linalool, borneol, linalyl acetate, terpinen-4-ol,
229 lavandulol, α -terpineol and lavandulyl acetate. A sample of lavender essential oil was first analysed
230 with the four NB column at different temperature rates (2, 3.5, 5 and 10°C/min) and at constant He
231 flow of 1 mL/min. Figure 4 gives the ES-GC-MS profiles of the lavender e.o. analysed with the
232 reference column. Table 5 reports enantiomer resolutions of optically active components, analysis
233 time, % analysis time reduction and separation measure *S*, when lavender e.o. was analysed with the
234 10, 5 and 2 m NB columns, and compares them to those of the reference column. The 1 m NB
235 column was not used, because the sample complexity was too high for its low efficiency. In this
236 case too, the results were considered to be acceptable when all optically active components had a
237 resolution above the limit of 1.5. These results also show that the 10m NB column can only be used
238 at 2°C/min, because at higher rates terpinen-4-ol, linalyl acetate and borneol are not base-line
239 separated. Figure 5a-b reports the ES-GC-MS profiles of the lavender e.o. analysed with the 5 m NB
240 column at 5 and 10°C/min. The highest applicable temperature rate was 5°C/min, because at
241 10°C/min terpinen-4-ol resolution decreased to 1.2, as is clear from its ion profiles at *m/z* 71, 111,
242 154, in Figure 5c-d. Under these conditions, the analysis time was reduced from about 43 min with
243 the reference column to about 14 minutes.

244 The 2 m NB column gave an analysis time similar to that of the 5 m NB column (about 13 min at
245 5°C/min) and comparable enantiomer resolutions (see table 5). When one of the enantiomers of the
246 markers co-eluted with other components, its resolution was calculated from their extracted ion MS
247 profiles, selecting suitable diagnostic ions. Figure 6 reports the 2 m NB column ES-GC-MS profiles
248 of the lavender e.o. together with those of extracted ion profiles of limonene (*m/z* 68, 107) linalool
249 (*m/z* 80, 93) and borneol (*m/z* 95, 110) whose enantiomers in TIC mode coelute with those of other
250 components. The use of shorter columns (2 and 5 m) results in lower values of the separation
251 measure *S*, that however are still sufficient for a reliable separation of all lavender e.o. optically
252 active markers. Table 6 reports the % area of each enantiomer calculated vs. the total area of the two
253 enantiomers, and gives the relative standard deviation (RSD%) of three optically active components
254 with different volatility, polarity and abundance (β -pinene, linalyl acetate and α -terpineol), analysed
255 with the 5 m NB column at different temperature rates and with the 10, 5, and 2 m NB columns at
256 the most effective temperature rate. These results show a high repeatability of enantiomer % areas
257 under all conditions. % Areas of each enantiomer were chosen instead of EE or ER because it better
258 shows repeatability under all conditions applied. The high RSD% of (*S*)-linalyl acetate is due to the
259 noise interference with area determination because it is in trace amount.

260 The bergamot e.o. was analysed with the same NB columns and under the same conditions as those
261 adopted for lavender e.o.. Five markers were selected, i.e. β -pinene, limonene, linalool, linalyl
262 acetate and α -terpineol. This e.o. was chosen partly because it contains a highly predominant
263 component (limonene) in its turn with a very high ER (EE), together with a set of minor markers
264 important to define its quality. This situation is quite common in the field of real-world samples of
265 natural origin. Narrow bore column capacity is lower than that of conventional I.D. columns,
266 therefore the simultaneous analysis of all markers produces either (*R*)-limonene peak distortion or,
267 with a suitably diluted sample, missing of minor components. The tailing of (*R*)-limonene produces
268 a peak asymmetry of 5.4 with the reference column and of 3.9 and 6.8 for the 2 and 5 m NB
269 columns, the late giving a reduction of *calculated* resolution below 1.5, although the two limonene
270 enantiomers are more than baseline separated, as is clear from figure 7. A similar situation, though

271 with lower relative abundance, also occurred with linalool in lavender e.o.. In this case, the
272 limonene resolution was calculated on the basis of the peak asymmetry value of (R)-enantiomer
273 (see table 7). A series of analyses on the bergamot e.o. at different dilutions (from 1: 200 to 1:5000)
274 also showed the linearity of (R)-limonene signal under the applied conditions (data not reported).
275 Figure 7 gives the ES-GC-MS profiles of a sample of bergamot e.o. analysed with the reference
276 column, and with the 2 and 5 m NB columns, together with the extract ion MS profiles of limonene
277 (m/z 68) on the reference 2 m NB and 5m NB columns, required because of the coelution of
278 ocimene and 2a limonene enantiomer on all columns, and of α -terpinene and 2b enantiomer on the
279 2 m NB column. Table 7 reports the enantiomer resolutions of five bergamot e.o. optically active
280 markers, their analysis time, % analysis time reduction and *S* after analysis with the 10, 5, and 2m
281 NB columns, analysed under different analytical conditions. In this case the limiting marker is
282 linalyl acetate, whose enantiomers can only be separated at the base line at 2°C/min with the 10 m
283 NB and up to 5°C/min with the 5 m NB column, while with the 2 m NB column they are separated
284 in all the conditions. With the latter column the extract ion approach has only to be applied to
285 limonene to obtain a correct EE or ER.
286 These results confirms those obtained with lavender e.o. analysis and show that 2 and 5 m NB
287 columns can successfully be used to speed up bergamot e.o. analysis, affording a very marked
288 reduction in analysis time: about 65% with the 5 m NB column at 5°C/min, and about 80% with the
289 2 m NB column at 10°C/min. The 2 and 5 m NB columns, at 5°C/min, gave separation measures (*S*)
290 that were reasonably lower than those of the reference column (max 35%), thus enabling a
291 temperature rate of 10°C/min to be applied to the 2 m NB column, with a further reduction in
292 analysis time.

293

294 **4. Conclusions**

295 The results show that short conventional and narrow bore columns, combined with CD
296 derivatives as chiral selectors, and MS as detector, can successfully be used for fast ES-GC analysis
297 of real-world samples. These results have been obtained by adopting: 1) a resolution limit of 1.5 to
298 enable a correct EE and ER determination; 2) mass spectrometry as detector to operate in both TIC
299 and extract ion modes, to avoid problems due to co-elution with other components in the EE and ER
300 determination. This assumption is correct because MS is not a chiral probe and the mass spectra of
301 two enantiomers are indistinguishable; 3) temperature rates up to 10°C/min, in consideration that
302 the enantioselective separations with CD as chiral selector are strongly conditioned by temperature,
303 and 4) short conventional and narrow bore columns that not only contribute to shortening analysis
304 time, but may also increase enantioselectivity and (at least partially) compensate for the loss of
305 efficiency due to reduced length, which lowers elution temperature. The combination of these
306 approaches enabled us to decrease analysis time by as much as 85% compared to that with
307 conventional I.D. columns, while keeping a resolution suitable for a correct EE and ER
308 determination, and to increase the number of samples analysable in a routine laboratory very
309 markedly. These results also show that new CD derivatives with high enantioselectivity are still
310 necessary not only to increase the number of optically active compounds separated or improve their
311 resolution, but also to provide further speeding-up of ES-GC-(MS) routine analysis.
312 Further studies are under way to examine the influence of selectivity tuning with different CD
313 derivatives as chiral selectors and diluting phases.

314

315 **Acknowledgements**

316 This research was carried out within the project entitled: "Sviluppo di metodologie
317 innovative per l'analisi di prodotti agroalimentari" (FIRB Cod.: RBIP06SXMR_002) of the
318 *Ministero dell'Istruzione, dell'Università e della Ricerca* (MIUR) (Italy).

319

320 **References**

321 [1] V. Schurig, H.-P. Nowotny, J. Chromatogr. 441 (1988) 155

322 [2] H.-P. Nowotny, D. Schmalzing, D. Wistuba, V. Schurig, J. High Resolut. Chromatogr. 12
323 (1989) 383.

324 [3]. V. Schurig, M. Jung, D. Schmalzing, M. Schleimner, J. Duvekot, J.C. Buyten, J. A. Peene, P.
325 Mussche, J. High Resolut. Chromatogr. 13 (1990) 470.

326 [4] G. Schomburg, H. Husmann, E. Hubinger, W.A. Koenig, J. High Resolut. Chromatogr.
327 Chromatogr. Commun.Resolut. Chromatogr 7 (1984) 404.

328 [5] A. Bernreuther, N. Christoph, P. Schreier, J.Chromatogr. 481 (1989) 363.

329 [6] A. Mosandl, U. Hener, U. Hagenauer-Hener, A. Kustermann, J. High Resolut. Chromatogr. 12
330 (1989) 532.

331 [7] L. Mondello, M. Catalfamo, G. Dugo, P. Dugo, J.Chromatogr. Sci. 36 (1998) 201.

332 [8] M. Jung, D. Schmalzing, V. Schurig, J.Chromatogr. 552 (1991) 43.

333 [9] V. Schurig, H.-P. Nowotny, Angew.Chem. Int. Ed. Engl., 29 (1990) 939.

334 [10] R. Shellie, P. Marriott, C. Cornwell, Anal. Chem. 74 (2002) 5426.

335 [11] R. Shellie, L. Mondello, G. Dugo, P. Marriot, Flavour Fragr. J. 19 (2004) 582.

336 [12] H. Van den Dool, P. D. Kratz, J. Chromatogr. 11 (1963) 463.

337 [13] E. Liberto, C. Cagliero, B. Sgorbini, C. Bicchi, D. Sciarrone, B. D'Acampora Zellner, L.
338 Mondello, P. Rubiolo, J. Chromatogr. A, 1195 (2008) 117.

339 [14] C.A. Cramers, H.-G. Janssen, M.M.Van Deursen, P.A. Leclercq, J. Chromatogr. A 856 (1999)
340 315.

341 [15] F.David, D.R.Gere, F. Scanlan, P. Sandra, J. Chromatogr. A 842 (1999) 309.

342 [16] K. Maštovská, S. J. Lehotay, J. Chromatogr. A 1000 (2003) 153 and references cited therein.

343 [17] M. Kirkchner, E. Matisová, S. Hrouzková, J. De Zeeuw, J. Chromatogr. A 1090 (2005) 126.

344 [18] M. Dömötörova, M. Kirkchner, E. Matisová, J. De Zeeuw, J. Sep. Sci. 29 (2006) 1051.

345 [19] P. Rubiolo, E. Liberto, B. Sgorbini, R. Russo, J.-L. Veuthey, C. Bicchi, J. Sep. Sci. 31 (2008)
346 1074.

347 [20] V. Schurig, J.Chromatogr. A 906 (2001) 275 and references cited therein.

348 [21] P. A. Levkin, V. Schurig, J.Chromatogr. A 1184 (2008) 309 and references cited therein.

349 [22] I. Hardt, W. A. Koenig, J. Microcol. Sep., 5 (1993) 35.

350 [23] M. Lindstroem, J. High Resolut. Chromatogr. 14 (1991) 765.

351 [24] C. Bicchi, G. Artuffo, A. D'Amato, G.M. Nano, A. Galli, M. Galli, presented at the 2nd
352 International Symposium on Chiral Discrimination, Rome, 27-31 May .

353 [25] C. Bicchi, G. Artuffo, A. D'Amato, A. Galli and M. Galli, Chirality 4 (1992) 125.

354 [26] H. Grosenick, V. Schurig, J. Constante, A. Collet, Tetrahedron Asymmetr. 6 (1995) 87.

355 [27] V. Schurig, H. Czesla, Enantiomer, 6 (2001) 107.

356 [28] C. Bicchi, G. Artuffo, A. D'Amato, V. Manzin, J. High Resolut. Chromatogr. 15 (1992) 710.

357 [29] European Pharmacopoeia 6th Edn. 2008.

358 [30] C. Bicchi, A. D'Amato, V. Manzin, A. Galli, M. Galli, J.Chromatogr. A 742 (1996) 161.

359 [31] L.M. Blumberg, M. Klee, J.Chromatogr. A 933 (2001) 1.

360

361 Captions to figures

362

363 Figure 1: TIC ES-GC-MS profiles of α -HCH analysed with different 0.25 mm I.D. columns.

364

365 Figure 2: ES-GC-MS profile of linalool analysed with the four NB columns under investigation at 1
366 mL/min flow rate and 10°C/min temperature rate.

367

368

369 Figure 3: Resolution *versus* flow and temperature rates (A) and retention time *versus* flow and
370 temperature rates (B) of γ -pentadecalactone on different columns

371

372 Figure 4: ES-GC-MS profile of the lavender essential oil analysed with the reference column (25 m
373 x 0.25 mm I.D., 0.25 μ m d_f). For analysis conditions see text. Peak identification: 1: α pinene, 2:
374 camphene, 3: β -pinene, 4: β -phellandrene, 5: limonene, 6: 1-octen-3-ol, 7: camphor, 8: linalool, 9:
375 borneol, 10: linalyl acetate, 11: 4-terpineol, 12: lavandulol, 13: α -terpineol, 14: lavandulyl acetate;
376 a: (*R*) enantiomer, b: (*S*) enantiomer. 6: not separated; 7a and 12a: not detected;

377

378 Figure 5: ES-GC-MS profiles of the lavender e.o. analysed with the 5 m NB column at 5 (a) and
379 10°C/min (b). Extract ion profiles of terpinen-4-ol (71,111,154 m/z) at 5 (c) and 10°C/min (d). For
380 analysis conditions see text. For peak identification, see caption to Figure 4.

381

382 Figure 6: 2m NB columns ES-GC-MS profiles of the lavender e.o. and extracted ion profiles of
383 limonene, linalool and borneol. For analysis conditions see text. For peak identification, see caption
384 to Figure4.

385

386 Figure 7: ES-GC-MS profiles of bergamot e.o. (top) and extract ion profiles of limonene (bottom)
387 obtained with different columns and under different analysis conditions. Peak identification: 1: β
388 pinene, 2: limonene, 3: linalool, 4: linalyl acetate, 5: α terpineol; a: (*R*) enantiomer, b: (*S*)
389 enantiomer. %areas on total limonene area: (*S*)-limonene 2%, (*R*)-limonene 98% .

Table 1: Columns and analysis conditions adopted for the present study

Conventional (internal diameter 0.25 mm)				Narrow bore (internal diameter 0.10 mm)				
<i>Columns</i>		<i>Analysis conditions</i>		<i>Columns</i>		<i>Analysis conditions</i>		
Lenght (m)	Film thickness (μm)	Flow rate (mL/min)	Temperature rate ($^{\circ}\text{C}/\text{min}$)	Lenght (m)	Film thickness (μm)	Flow rate (mL/min)	Temperature rate ($^{\circ}\text{C}/\text{min}$)	
25	0.25	1	2	10	0.10	0.4	2	
25	0.15			5	0.10		0.7	3.5
5	0.15			2	0.10		1	5
				1	0.10		10	

Table 3: Separation measure S calculated with the chiral test mixture under all conditions applied and for all columns investigated.

	S			
	2°C/min	3.5°C/min	5°C/min	10°C/min
Conv 25 m 0.25 μm	942			
Conv 25 m 0.15 μm	972			
Conv 5 m	491			
NB 10 m	477	442	395	347
NB 5 m	826	775	754	686
NB 2 m	727	705	697	562
NB 1 m	216	277	278	266

Table 4: % Areas of each enantiomer on the total areas of the (R)- and (S)-linalool enantiomers and related SD and RSD%, analysed with the 10, 5, 2, and 1 m NB columns under all conditions applied, and with the reference column.

Column	% Area	Mean	SD	RSD%
	(R)/(S)-linalool			
Reference	(R)	49.0		
	(S)	51.0		
10 m (all anal cond.)	(R)	48.9	0.4	0.8
	(S)	51.1	0.4	0.7
5 m (all anal cond.)	(R)	49.0	0.4	0.8
	(S)	51.0	0.4	0.8
2 m (all anal cond.)	(R)	49.5	0.5	1.1
	(S)	50.5	0.5	1.1
1 m (all anal cond.)	(R)	48.8	0.4	0.7
	(S)	51.2	0.4	0.7
All columns and conditions	(R)	49.1	0.3	0.6
	(S)	50.9	0.3	0.6

Table 6: % Area of each enantiomer calculated vs. the total area of the two enantiomers of β -pinene, linalyl acetate and α -terpineol in lavender e.o., and related RSD%, analysed with the 5 m NB column at different temperature rates and with the 10, 5, and 2 m NB columns at the most effective temperature rate.

Column	%Area				Mean	RSD%
	5 m					
Temperature rate (°C/min)	2	3.5	5	10		
(+)- β -Pinene	56.3	54.7	54.3	54.5	55.0	1.7
(-)- β -Pinene	43.7	45.3	45.7	45.5	45.0	2.1
(<i>R</i>)-Linalyl acetate	99.5	99.7	99.6	99.6	99.6	0.1
(<i>S</i>)-Linalyl acetate	0.5	0.4	0.4	0.4	0.4	12.7
(<i>S</i>)- α -Terpineol	28.4	29.0	28.8	29.6	28.9	1.8
(<i>R</i>)- α -Terpineol	71.7	71.0	71.2	70.4	71.1	0.8
Column	Ref.	10 m	5 m	2 m	Mean	RSD%
Temperature rate (°C/min)	2	2	5	5		
(+)- β -Pinene	55.6	55.0	54.3	55.2	55.1	1.0
(-)- β -Pinene	44.4	45.0	45.7	44.8	45.0	1.2
(<i>R</i>)-Linalyl acetate	99.4	99.6	99.6	99.4	99.5	0.1
(<i>S</i>)-Linalyl acetate	0.6	0.4	0.4	0.6	0.5	18.0
(<i>S</i>)- α -Terpineol	27.5	28.7	28.8	28.2	28.3	2.0
(<i>R</i>)- α -Terpineol	72.5	71.4	71.2	71.8	71.7	0.8

Table 2: Retention time (R_t) of the second enantiomer, elution temperature of both enantiomers (T_{el1-2}), resolution (R) and % analysis time reduction (ATR) calculated on the second eluting enantiomer of the chiral standard compounds investigated.

Compound	<i>25 m x 0.25 μm d_f</i>			<i>25 m x 0.15 μm d_f</i>			<i>5 m x 0.15 μm d_f</i>			ATR <i>25 m 0.15 μm vs 25 m 0.25 μm</i>	ATR <i>5 m 0.15 μm vs 25 m 0.15 μm</i>	ATR <i>5 m 0.15 μm vs 25 m 0.25 μm</i>
	t_R (min) II	T_{el} (°C) I-II	R	t_R (min) II	T_{el} (°C) I-II	R	t_R (min) II	T_{el} (°C) I-II	R			
Limonene	18.67	85-87	8.2	11.72	72-73	8.4	3.66	57-58	4.1	37.2%	68.7%	80.4%
Linalool	25.15	99-100	7.9	18.77	86-87	7.8	7.94	65-66	4.2	25.4%	57.8%	68.4%
Linalyl acetate	27.80	105-106	3.0	20.60	90-91	3.7	8.76	67-68	2.5	25.9%	57.5%	68.5%
Linalyl propionate	31.78	113-114	1.0	24.15	98-98	1.2	11.16	72-72	1.0	24.9%	53.8%	64.9%
γ -Hexalactone	28.30	103-106	11.7	21.99	90-93	13.6	8.97	66-69	5.2	22.3%	59.2%	68.3%
γ -Heptalactone	32.50	112-115	11.9	26.04	99-102	13.9	13.08	73-76	5.2	20.6%	49.8%	62.7%
γ -Octalactone	36.73	121-123	10.0	30.05	107-110	11.6	16.69	80-83	4.8	18.2%	44.4%	56.7%
γ -Decalactone	46.06	141-142	5.7	38.96	126-128	7.2	24.88	98-100	3.8	15.4%	36.1%	46.0%
γ -Undecalactone	50.85	151-152	4.9	43.52	136-137	6.1	29.12	106-108	3.4	14.4%	34.2%	42.7%
γ -Dodecalactone	55.60	160-161	4.0	48.05	145-146	4.8	33.34	115-117	3.0	13.6%	33.1%	40.0%
γ -Tetradecalactone	65.41	179-180	2.4	56.82	163-164	3.2	41.34	131-133	2.7	13.1%	27.2%	36.8%
γ -Pentadecalactone	69.45	188-189	1.9	60.65	171-172	2.5	45.01	139-140	1.9	12.7%	25.7%	35.2%
δ -Octalactone	36.98	123-124	3.7	30.22	109-110	3.4	16.84	83-84	2.2	18.3%	44.3%	54.4%
α -HCH	29.00	156-158	5.9	21.09	140-142	8.8	7.60	114-115	5.3	27.3%	64.0%	73.8%
Trichlorfon	26.27	149-152	14.4	20.15	138-140	12.2	8.56	116-117	5.8	23.3%	57.5%	67.4%
<i>trans</i> -Chlordane	46.19	192-192	0.7	36.56	172-173	1.0	18.70	137-137	1.2	20.8%	48.9%	69.5%
Heptachlor	36.48	173-173	0.6	26.83	153-154	1.0	10.96	122-122	0.9	26.5%	59.2%	70.0%

Table 5: Enantiomer resolutions of optically active components, analysis time, % analysis time reduction and separation measure *S* obtained for a lavender e.o. analysed with the 10, 5 and 2 m NB columns in comparison to those of the reference column (Ref.).

	Ref.	10 m NB				5 m NB				2 m NB				
Temperature rate (°C/min)	2	2	3.5	5	10	2	3.5	5	10	2	3.5	5	10	
Analysis time (min)	40.67	29.35	19.60	14.77	8.67	26.81	17.77	13.60	8.01	28.83	16.57	12.73	7.56	
% Analysis time reduction		27.8	51.8	63.7	78.7	34.1	56.3	66.6	80.3	29.1	59.2	68.7	81.4	
<i>S</i>	1077	407	374	333	286	840	808	744	652	759	701	688	582	
	<i>I</i> ^{T*}	Resolution												
α-Pinene (1)	(<i>R</i>)921/(<i>S</i>)923	1.2	NS	NS	NS	NS	NS	NS	NS	NS	NS	NS	NS	NS
Camphene (2)	(-)-917/(+)-923	6.8	2.2	1.9	1.8	1.7	6.7	5.6	5.6	4.9	3.9	3.5	3.4	3.5
β-Pinene (3)	(+)-944/(-)-955	5.0	1.8	1.5	1.4	1.2	5.3	4.8	4.4	3.5	3.3	3.3	2.9	2.6
β-Phellandrene (4)	(-)-1049/(+)-1060	6.1	2.1	1.8	1.7	1.5	6.2	5.3	5.1	4.5	3.1	2.8	2.6	2.4
Limonene (5)	(<i>S</i>)1056/(<i>R</i>)1072	9.1	2.8	2.5	1.9	1.9	9.0	8.2	8.0	5.0	5.4	4.2	3.8	2.9
1-Octen-3-ol (6)	1126	NS	NS	NS	NS	NS	NS	NS	NS	NS	NS	NS	NS	NS
Camphor (7)	(<i>S</i>)1133/(<i>R</i>)1141	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E
Linalool (8)	(<i>R</i>)1174/(<i>S</i>)1189	6.1	2.6	2.4	2.3	1.9	3.1	3.3	3.2	2.5	3.1	3.4	3.0	3.0
Borneol (9)	(<i>S</i>)1192/(<i>R</i>)1200	3.0	1.6	1.2	1.2	1.0	2.9	2.5	2.0	1.9	2.1	2.0	1.7	1.6
Linalyl acetate (10)	(<i>R</i>)1231/(<i>S</i>)1237	3.2	2.6	2.5	1.6	NS	3.8	3.0	2.8	2.4	4.1	3.8	3.1	1.9
Terpinen-4-ol (11)	(<i>S</i>)1248/(<i>R</i>)1253	2.2	1.6	1.0	NS	NS	2.4	2.1	2.0	1.2	2.6	2.3	2.0	1.3
Lavandulol (12)	(<i>S</i>)1250/(<i>R</i>)1273	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E
α-Terpineol (13)	(<i>S</i>)1296/(<i>R</i>)1309	6.0	3.2	2.8	2.5	1.9	6.9	5.6	5.5	4.2	6.5	5.6	4.7	4.0
Lavandulyl acetate (14)	(<i>R</i>)1259/(<i>S</i>)1263	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E	1E

NS: not baseline separated

1E: only (*R*) enantiomer found

*obtained on the reference column

Table 7: Enantiomer resolutions of five bergamot e.o. optically active components, analysis time, % analysis time reduction and separation measure *S* analysed with the 10, 5 and 2 m NB columns in comparison to those of the reference column (Ref.).

		Ref.	10 m NB				5 m NB				2 m NB			
Temperature rate (°C/min)		2	2	3.5	5	10	2	3.5	5	10	2	3.5	5	10
Analysis time (min)		41.81	30.62	20.04	15.24	8.90	28.11	18.47	14.07	8.23	26.11	17.28	13.20	7.76
% Analysis time reduction			26.8	52.1	63.5	78.7	32.8	55.8	66.3	80.3	37.5	58.7	68.4	81.4
<i>S</i>		1148	431	408	388	322	973	892	791	760	876	805	750	669
	<i>I</i> ^{T**}	Resolution												
β-Pinene (1)	(+)-944/(-)-955	4.0	1.2	1.1	1.1	1.0	3.1	2.9	2.8	2.7	2.3	2.3	2.2	2.0
Limonene (2)	(<i>S</i>)-1056/(<i>R</i>)-1072	4.8	1.2*	1.0*	1.0*	0.9*	3.8*	3.1*	2.8*	2.8*	2.2*	2.4*	2.2*	2.2*
Linalool (3)	(<i>R</i>)-1174/(<i>S</i>)-1189	6.7	3.0	2.8	2.6	1.7	5.2	5.1	4.6	4.3	4.8	4.2	4.2	3.6
Linalyl acetate (4)	(<i>R</i>)-1231/(<i>S</i>)-1237	3.0	1.5	0.9	NS	NS	2.6	2.4	1.9	0.9	2.8	2.7	2.5	1.9
α-Terpineol (5)	(<i>S</i>)-1296/(<i>R</i>)-1309	5.0	3.8	2.8	2.6	1.6	7.3	5.2	5.1	4.3	7.8	5.8	5.2	3.1

NS: not baseline separated

* calculated through peak asymmetry

**obtained on the reference column

Figure 1




Figure 2




Figure 3






Figure 4




Figure 5




Figure 6







Figure 7

