

Supplementary Table 1. Comparison between Spadaro's and manufacturer's methods in wine matrix.

Clean-up method	Validation levels ($\mu\text{g/L}$)	Recovery (%) ($n = 3$)	RSD (%) ($n = 3$)	LOD ($\mu\text{g/L}$) ($n = 3$)	LOQ ($\mu\text{g/L}$) ($n = 3$)
Method of Spadaro <i>et al.</i>	2	84.03	19.22		
	10	102.24	15.51	0.14	0.48
	20	93.07	10.53		
Method of manufacturer	2	82.31	17.84		
	10	99.15	12.26	0.12	0.41
	20	93.96	13.95		