SCIENTIFIC DATA (11011)

OPEN Data Descriptor: RNA-Seq profiling in peripheral blood mononuclear cells of amyotrophic lateral sclerosis patients and controls

Received: 5 July 2018 Accepted: 14 December 2018 Published: 5 February 2019

Susanna Zucca^{1,*}, Stella Gagliardi^{1,*}, Cecilia Pandini^{1,2}, Luca Diamanti^{3,4}, Matteo Bordoni^{1,3}, Daisy Sproviero¹, Maddalena Arigoni⁵, Martina Olivero⁵, Orietta Pansarasa¹, Mauro Ceroni^{3,4}, Raffaele Calogero^{5,†} & Cristina Cereda^{1,†}

Coding and long non-coding RNA (IncRNA) metabolism is now revealing its crucial role in Amyotrophic Lateral Sclerosis (ALS) pathogenesis. In this work, we present a dataset obtained via Illumina RNA-seq analysis on Peripheral Blood Mononuclear Cells (PBMCs) from sporadic and mutated ALS patients (mutations in FUS, TARDBP, SOD1 and VCP genes) and healthy controls. This dataset allows the wholetranscriptome characterization of PBMCs content, both in terms of coding and non-coding RNAs, in order to compare the disease state to the healthy controls, both for sporadic patients and for mutated patients. Our dataset is a starting point for the omni-comprehensive analysis of coding and IncRNAs, from an easy to withdraw, manage and store tissue that shows to be a suitable model for RNA profiling in ALS.

Design Type(s)	parallel group design • individual genetic characteristics comparison design
Measurement Type(s)	transcription profiling assay
Technology Type(s)	RNA sequencing
Factor Type(s)	diagnosis • protein_altering_variant
Sample Characteristic(s)	Homo sapiens • peripheral blood mononuclear cell

¹Genomic and post-Genomic Center, IRCCS Mondino Foundation, Pavia, Italy. ²Department of Biology and Biotechnology "L. Spallanzani", University of Pavia, Pavia, Italy. 3Department of Brain and Behavioral Sciences, University of Pavia, Pavia, Italy. ⁴General Neurology, IRCCS Mondino Foundation, Pavia, Italy. ⁵Department of Molecular Biotechnology and Health Sciences, Bioinformatics and Genomics Unit, University of Turin, Torino, Italy. *These authors contributed equally to this work. †These authors jointly supervised this work. Correspondence and requests for materials should be addressed to S.Z. (email: susanna.zucca@mondino.it)

Background & Summary

Deep sequencing technologies allow high-throughput massive-parallel RNA sequencing, permitting an extensive characterization of transcriptomic profile of cell populations and tissues. RNA-seq analyses are the gold standard for in depth characterization of global changes in gene expression levels between different conditions (e.g.: disease and healthy subjects or treated and untreated samples)¹.

Such comprehensive analysis are currently bringing the light on the expression of largely unexplored RNA classes. While coding RNA involvement in a wide plethora of disorders has been subject of intense research, the more recently revealed class of lncRNAs is only starting to be characterized. LncRNAs are RNA transcripts greater than 200 nucleotides that lack an open reading frame and therefore do not encode proteins². The heterogeneity of lncRNAs, both in biological type and function, is a major obstacle that makes it difficult to extract information about lncRNAs function and to enrich our knowledge³. While coding genes are widely annotated, high-quality catalogues of lncRNAs and information on the tissues in which they are expressed are recently being constructed⁴⁻⁶. Recent efforts are directed to characterize this largely unexplored functional component of the genome. Deep sequencing is a precious instrument to obtain information suitable to investigate both coding and lncRNAs simultaneously in the same sample, and this comprehensive RNA analysis can elucidate fundamental mechanisms in gene expression regulation at post-transcriptional level^{7,8}.

There is mounting evidence that altered RNA metabolism, both involving coding and non-coding RNAs, plays an important role in ALS pathogenesis^{9–13}. ALS is an adult-onset, progressive and fatal neurodegenerative disease, caused by the selective loss of both upper and lower motor neurons in the cerebral cortex, brainstem and spinal cord. The pathogenesis of the disease is still unknown. Moreover, the investigation of ALS mechanism is challenging due to the difficulty to obtain a sufficient amount of biological material from ALS patients: it is indeed impossible to access biological samples from the most affected areas, such as spinal cord. PBMCs have been shown to be a convincing and realistic model for ALS, since many pathways, typically located in neurons, are also activated in these cells^{8,14,15}.

Furthermore, low ALS incidence (ranging between 1 and 2.5/100'000)¹⁶ and low percentage of familial ALS cases (5-10%) impose the requirement of data sharing and data aggregation among different hospitals and laboratories in order to obtain a sufficient sample size. Detailed curation of these datasets is paramount for accurate interpretation, widespread dissemination, and repurposing of data.

In this paper, we present a dataset suitable for the analysis for both coding RNAs and lncRNAs in PBMCs of 15 sporadic ALS (sALS) patients, 9 ALS patients with mutation in classical ALS-genes (*SOD1*, *FUS*, *TARDBP* and *VCP*) and 7 healthy controls. More in detail, these data are suitable for antisense lncRNAs detection and characterization because of the chemistry of the protocol properly design to maintain strandness information. Strand-specific RNA-Seq permits a more complete understanding of the transcriptome, enhancing resolution for sense/antisense profiling¹⁷.

Part of this dataset was used to perform an analysis of differentially expressed transcripts (coding RNAs and lncRNAs)⁸ aiming at investigating the importance of sense and antisense RNA regulation in central and peripheral systems in ALS disease.

To our knowledge, no dataset of RNA expression in PBMCs from ALS patients is publicly available. Our dataset is a starting point for the omni-comprehensive analysis of RNA classes, from an easy to withdraw, manage and store tissue that shows to be a suitable model for RNA profiling in ALS⁸.

It was shown that in RNA-seq studies, increasing sample size is the best choice to enhance statistical power¹⁸. For these reasons, the re-usability of this dataset is related to the possibility to expand studies with smaller sample sizes, especially for mutated patients, to identify new sense and antisense transcripts with altered expression in both sporadic and mutated ALS patients compared to healthy controls and to provide an open-angle point of view for a broad-spectrum characterization.

Methods

Isolation of human PBMCs from ALS patients and healthy controls

Subject enlistment. 24 ALS patients and 7 age- and sex-matched healthy controls (CTRL) were recruited after obtaining written informed consent (Table 1). ALS patients underwent clinical and neurologic examination at IRCCS National Neurological Institute "C. Mondino" (Pavia, Italy). All patients were diagnosed with ALS as defined by El Escorial criteria. All ALS patients were analysed for causative mutations in SOD1, TARDBP, FUS, C9orf72, ANG and VCP genes: 15 patients resulted to be sporadic (sALS) and are indicated as ALS-s1, ALS-s2 and so on. Three patients were mutated in SOD1 gene (SOD1-m1, SOD1-m2 and SOD1-m3), three in FUS gene (FUS-m1, FUS-m2 and FUS-m3), two in TARDBP gene (TARDBP-m1 and TARDBP-m2) and one in VCP gene (VCP-m1). The seven control subjects (indicated as CTRL1-CTRL7) have been recruited at the Transfusional Service and Centre of Transplantation Immunology, Foundation San Matteo, IRCCS (Pavia, Italy). All details are reported in Table 1.

The study protocol to obtain PBMCs from patients and controls was approved by the Ethical Committee of the National Neurological Institute "C. Mondino", IRCCS (Pavia, Italy). Before being enrolled, the subjects participating in the study signed an informed consent form (Protocol n°375/04 – version 07/01/2004). All experiments were performed in accordance with relevant guidelines and regulations. These methods and following methods are expanded versions of descriptions in our related work.

CTRLI GSE106443 M 48 na na CTRL2 GSE10643 M 60 na na CTRL3 GSE10643 M 68 na na CTRL4 GSE115259 M 38 na na CTRL5 GSE115259 F 40 na na CTRL6 GSE115259 F 40 na na CTRL6 GSE115259 F 40 na na CTRL7 GSE10543 M 66 65 na ALS-41 GSE10643 F 61 0 na ALS-32 GSE10643 F 61 0 na ALS-34 GSE10643 F 61 0 na ALS-35 GSE10643 F 55 33 na ALS-49 GSE10643 F 56 55 na ALS-49 GSE10643 F 58 3 <	Sample	GEO	Sex	Age	Age of onset	Mutation
CREATS GSE106443 M 68 na na CTRL4 GSE115259 M 38 na na CTRL5 GSE115259 F 36 na na CTRL6 GSE115259 F 40 na na CTRL7 GSE115259 F 40 na na ALS-41 GSE106443 M 66 65 na ALS-32 GSE106443 M 63 61 na ALS-33 GSE106443 F 61 60 na ALS-34 GSE106443 F 70 68 na ALS-45 GSE106443 F 55 55 na ALS-45 GSE106443 F 56 55 na ALS-48 GSE106443 F 83 73 na ALS-49 GSE106443 F 83 73 na ALS-410 GSE106443 F 88 <td< td=""><td>CTRL1</td><td>GSE106443</td><td>М</td><td>48</td><td>na</td><td>na</td></td<>	CTRL1	GSE106443	М	48	na	na
CTRL4 GSE115259 M 38 ma na CTRL5 GSE115259 F 36 ma na CTRL6 GSE115259 F 40 na na CTRL7 GSE116259 F 40 na na ALS-81 GSE106443 M 66 65 na ALS-82 GSE106443 M 63 61 na ALS-83 GSE106443 F 61 60 na ALS-84 GSE106443 F 70 68 na ALS-85 GSE106443 F 55 53 na ALS-86 GSE106443 F 55 53 na ALS-87 GSE106443 F 83 73 na ALS-89 GSE106443 F 88 58 na ALS-810 GSE106443 F 88 58 na ALS-81 GSE106443 F 86 <td< td=""><td>CTRL2</td><td>GSE106443</td><td>М</td><td>60</td><td>na</td><td>na</td></td<>	CTRL2	GSE106443	М	60	na	na
CREL5 GSE115259 F 36 na na CTRL6 GSE115259 F 40 na na CTRL7 GSE115259 F 40 na na ALS-s1 GSE106443 M 66 65 na ALS-s2 GSE106443 M 63 61 na ALS-s3 GSE106443 F 61 60 na ALS-s4 GSE106443 F 70 68 na ALS-s6 GSE106443 F 55 53 na ALS-s6 GSE106443 F 55 53 na ALS-s7 GSE106443 F 56 55 na ALS-s8 GSE106443 F 83 73 na ALS-s8 GSE106443 F 58 58 na ALS-s10 GSE106443 F 58 58 na ALS-s11 GSE106443 M 66 <	CTRL3	GSE106443	М	68	na	na
CTRL6 GSE115259 F 40 na na CTRL7 GSE115259 F 40 na na ALS-s1 GSE106443 M 66 65 na ALS-s2 GSE106443 M 63 61 na ALS-s3 GSE106443 F 61 60 na ALS-s4 GSE106443 F 70 68 na ALS-s5 GSE106443 F 55 53 na ALS-s6 GSE106443 F 55 53 na ALS-s7 GSE106443 F 56 55 na ALS-s8 GSE106443 F 56 55 na ALS-s9 GSE106443 F 58 58 na ALS-s9 GSE106443 F 58 58 na ALS-s10 GSE106443 F 58 58 na ALS-s11 GSE106443 M 68	CTRL4	GSE115259	М	38	na	na
CTRL7 GSE115259 F 40 na na na ALS-s1 GSE106443 M 66 65 na na ALS-s2 GSE106443 M 63 61 na ALS-s3 GSE106443 F 61 60 na ALS-s3 GSE106443 F 70 68 na ALS-s5 GSE106443 F 70 68 na ALS-s5 GSE106443 F 70 68 na ALS-s5 GSE106443 F 70 68 na ALS-s6 GSE106443 F 70 68 na ALS-s6 GSE106443 F 70 70 68 na ALS-s6 GSE106443 F 70 70 68 na ALS-s6 GSE106443 F 70 70 70 70 70 70 70 70 70 70 70 70 70	CTRL5	GSE115259	F	36	na	na
ALS-s1 GSE106443 M 66 66 65 na ALS-s2 GSE106443 F 61 60 na ALS-s3 GSE106443 F 61 60 na ALS-s4 GSE106443 F 70 68 na ALS-s5 GSE106443 F 70 68 na ALS-s6 GSE106443 F 55 55 na ALS-s7 GSE106443 F 55 55 na ALS-s8 GSE06443 F 56 55 na ALS-s8 GSE06443 F 56 55 na ALS-s9 GSE106443 F 58 83 73 na ALS-s9 GSE106443 F 58 83 73 na ALS-s10 GSE106443 F 58 83 73 na ALS-s10 GSE106443 F 58 83 73 na ALS-s11 GSE106443 F 58 83 73 na ALS-s11 GSE106443 F 58 83 73 na ALS-s12 GSE105443 F 58 85 8 na ALS-s13 GSE105443 F 58 85 8 na ALS-s14 GSE105443 F 58 86 65 na ALS-s15 GSE115259 M 86 68 65 na ALS-s16 GSE115259 F 68 68 65 na ALS-s15 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 69 67 na ALS-s16 GSE115259 F 69 69 67 na ALS-s16 GSE106443 M 49 48 FUS R521C FUS-m1 GSE106443 F 50 49 FUS P160R FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE105259 F 50 49 FUS P160R FUS-m3 GSE105259 F 74 73 SOD1 G93D FUS-m3 GSE105259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP ASZT TARDBP-m2 GSE106443 F 52 50 TARDBP ASZT	CTRL6	GSE115259	F	40	na	na
ALS-s2	CTRL7	GSE115259	F	40	na	na
ALS-33 GSE106443 F 61 60 na ALS-34 GSE106443 F 70 68 na ALS-35 GSE106443 M 56 55 na ALS-36 GSE106443 F 55 53 na ALS-37 GSE106443 F 56 55 na ALS-38 GSE106443 F 83 73 na ALS-39 GSE106443 F 58 58 na ALS-310 GSE106443 F 58 58 na ALS-311 GSE106443 M 66 60 na ALS-312 GSE115259 M 86 84 na ALS-313 GSE115259 F 68 65 na ALS-316 GSE115259 F 69 67 na ALS-316 GSE106443 F 50 49 FUS P160R FUS-m2 GSE106443 F 50<	ALS-s1	GSE106443	М	66	65	na
ALS-84 GSE106443 F 70 68 na ALS-85 GSE106443 M 56 55 na ALS-86 GSE106443 F 55 53 na ALS-87 GSE106443 F 56 55 na ALS-88 GSE106443 F 83 73 na ALS-89 GSE106443 M 66 60 na ALS-810 GSE106443 F 58 58 na ALS-811 GSE106443 M 68 65 na ALS-812 GSE105259 M 86 84 na ALS-813 GSE15259 F 68 65 na ALS-815 GSE15259 F 69 67 na ALS-816 GSE15259 F 69 69 67 na ALS-816 GSE166443 F 70 70 NA ALS-816 GSE166443 F 70 70 NA ALS-816 GSE166443 F 70 NA ALS-816 NA ALS-81	ALS-s2	GSE106443	М	63	61	na
ALS-s5 GSE106443 M 56 55 na ALS-s6 GSE106443 F 55 53 na ALS-s7 GSE106443 F 56 55 na ALS-s8 GSE106443 F 83 73 na ALS-s9 GSE106443 M 66 60 na ALS-s10 GSE106443 F 58 58 na ALS-s10 GSE106443 M 68 65 na ALS-s11 GSE106443 M 68 65 na ALS-s12 GSE115259 M 86 84 na ALS-s13 GSE115259 F 68 65 na ALS-s15 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS R521C FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A3S2T TARDBP-m2 GSE106443 F 52 50 TARDBP A3S2T	ALS-s3	GSE106443	F	61	60	na
ALS-s6 GSE106443 F 55 53 na ALS-s7 GSE106443 F 56 55 na ALS-s8 GSE106443 F 83 73 na ALS-s9 GSE106443 M 66 60 na ALS-s10 GSE106443 F 58 58 na ALS-s11 GSE106443 M 68 65 na ALS-s12 GSE115259 M 86 84 na ALS-s13 GSE115259 F 68 65 na ALS-s15 GSE115259 F 69 67 na FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 F 74 73 SOD1 G147S SOD1-m2 GSE106443 <td< td=""><td>ALS-s4</td><td>GSE106443</td><td>F</td><td>70</td><td>68</td><td>na</td></td<>	ALS-s4	GSE106443	F	70	68	na
ALS-87 GSE106443 F 56 55 na ALS-88 GSE106443 F 83 73 na ALS-99 GSE106443 M 66 60 na ALS-10 GSE106443 F 58 58 na ALS-11 GSE106443 M 68 65 na ALS-11 GSE106443 M 68 65 na ALS-12 GSE115259 M 86 84 na ALS-13 GSE115259 F 68 65 na ALS-15 GSE115259 F 69 67 na ALS-16 GSE115259 F 69 67 na ALS-16 GSE115259 F 69 67 pa ALS-16 GSE115259 F 69 67 pa ALS-16 GSE115259 F 69 67 pa ALS-16 GSE115259 F 69 69 67 pa ALS-16 GSE115259 F 69 69 67 pa ALS-16 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP GS57D	ALS-s5	GSE106443	М	56	55	na
ALS-88 GSE106443 F 83 73 na ALS-99 GSE106443 M 66 60 na ALS-10 GSE106443 F 58 58 na ALS-11 GSE106443 M 68 65 na ALS-12 GSE115259 M 86 86 65 na ALS-13 GSE115259 F 68 65 na ALS-15 GSE115259 F 69 67 na ALS-16 GSE115259 F 69 67 na ALS-17 NO NA ALS-18 GSE115259 F 69 67 NA ALS-18 GSE115259 F 69 67 NA ALS-19 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m1 GSE106443 M 68 66 TARDBP G357D	ALS-s6	GSE106443	F	55	53	na
ALS-s9 GSE106443 M 66 60 na ALS-s10 GSE106443 F 58 58 na ALS-s11 GSE106443 M 668 65 na ALS-s12 GSE115259 M 86 84 na ALS-s13 GSE115259 F 68 65 na ALS-s15 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na ALS-s16 GSE106443 M 49 48 FUS R521C FUS-m1 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 F 74 73 SOD1 L106F SOD1-m2 GSE106443 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m1 GSE106443 F 52 50 TARDBP G357D	ALS-s7	GSE106443	F	56	55	na
ALS-s10 GSE106443 F 58 58 na ALS-s11 GSE106443 M 68 65 na ALS-s12 GSE115259 M 86 84 na ALS-s13 GSE115259 F 68 65 na ALS-s15 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 F 74 73 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE106443 F 55 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 55 50 TARDBP A382T TARDBP-m2 GSE106443 F 55 50 TARDBP G357D	ALS-s8	GSE106443	F	83	73	na
ALS-s11 GSE106443 M 68 65 na ALS-s12 GSE115259 M 86 84 na ALS-s13 GSE115259 F 68 65 na ALS-s15 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na EUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R EUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	ALS-s9	GSE106443	М	66	60	na
ALS-s12 GSE115259 M 86 84 na ALS-s13 GSE115259 F 68 65 na ALS-s15 GSE115259 M 71 70 na ALS-s16 GSE115259 F 69 67 na ALS-s16 GSE115259 F 69 67 na FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS G302A FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	ALS-s10	GSE106443	F	58	58	na
ALS-s13 GSE115259 F 68 65 na ALS-s15 GSE115259 M 71 70 na ALS-s16 GSE115259 F 69 67 na FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	ALS-s11	GSE106443	М	68	65	na
ALS-s15 GSE115259 M 71 70 na ALS-s16 GSE115259 F 69 67 na FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	ALS-s12	GSE115259	М	86	84	na
ALS-s16 GSE115259 F 69 67 na FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 M 68 66 TARDBP G357D	ALS-s13	GSE115259	F	68	65	na
FUS-m1 GSE106443 M 49 48 FUS R521C FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	ALS-s15	GSE115259	М	71	70	na
FUS-m2 GSE106443 F 50 49 FUS P160R FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	ALS-s16	GSE115259	F	69	67	na
FUS-m3 GSE115259 M 57 56 FUS G302A SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	FUS-m1	GSE106443	M	49	48	FUS R521C
SOD1-m1 GSE106443 M 50 49 SOD1 L106F SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	FUS-m2	GSE106443	F	50	49	FUS P160R
SOD1-m2 GSE106443 F 74 73 SOD1 G147S SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	FUS-m3	GSE115259	М	57	56	FUS G302A
SOD1-m3 GSE115259 F 58 57 SOD1 G93D TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	SOD1-m1	GSE106443	M	50	49	SOD1 L106F
TARDBP-m1 GSE106443 F 52 50 TARDBP A382T TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	SOD1-m2	GSE106443	F	74	73	SOD1 G147S
TARDBP-m2 GSE106443 M 68 66 TARDBP G357D	SOD1-m3	GSE115259	F	58	57	SOD1 G93D
	TARDBP-m1	GSE106443	F	52	50	TARDBP A382T
VCP-m1 GSE115259 M 51 50 VCP R191Q	TARDBP-m2	GSE106443	М	68	66	TARDBP G357D
	VCP-m1	GSE115259	М	51	50	VCP R191Q

Table 1. RNA-seq profiling to evaluate differential gene expression. All samples had the same source (blood sample) and were processed with the following steps: PBMCs isolation, RNA extraction and RNA-seq. All samples included in this study were Italian. All ALS patients had spinal onset and sample collection was performed after 6 months from exordium.

Isolation of human PBMCs. PBMCs were prepared by centrifugation. Peripheral blood was layered (density = 1.077) and centrifuged at 950 g for 30 min. After isolation on a Ficoll-Histopaque layer (Sigma Aldrich, Italy), cell viability was assayed by a trypan blue exclusion test and cells were then used for RNA extraction.

RNA extraction. Samples were homogenized and total RNA was isolated by Trizol reagent (Life Science Technologies, Italy) following the manufacturer's protocol. RNAs were quantified using a Nanodrop ND-100 Spectrophotometer (Nanodrop Technologies, Wilmington, USA) and a 2100 Bioanalyzer (Agilent RNA 6000 Nano Kit, Waldbronn, Germany); RNAs with a 260:280 ratio of ≥1.5 and an RNA integrity number of ≥8 were deep sequenced (Supplementary Figure 1).

RNA-seq library preparation, sequencing and analysis. Sequencing libraries were prepared with the Illumina TruSeq Stranded RNA Library Prep, version 2, Protocol D, using 500-ng total RNA (Illumina, USA). The qualities of the libraries were assessed by 2100 Bioanalyzer with a DNA1000 assay. Libraries were quantified by qPCR using the KAPA Library Quantification kit for Illumina sequencing platforms (KAPA Biosystems); RNA processing has been carried out using Illumina NextSeq 500 Sequencing, using 8 samples for each run mixing samples and controls in each flowcell to avoid not manageable batch-effects.

Sample	Number of input reads	Average input read length	Uniquely mapped reads %	Number of reads mapped to multiple loci	% of reads mapped to multiple loci	Seq. batch
CTRL1	9.85E + 07	2 × 75	22.68%	5.55E + 07	56.38%	1
CTRL2	8.51E + 07	2 × 75	25.65%	5.49E + 07	64.53%	1
CTRL3	9.00E + 07	2 × 75	26.99%	5.53E + 07	61.47%	2
CTRL4	2.00E + 07	2 × 75	6.50%	1.23E + 07	61.53%	2
CTRL5	2.28E + 08	2 × 75	11.02%	1.71E + 08	74.93%	3
CTRL6	1.05E + 08	2 × 150	4.17%	2.88E + 07	27.54%	4
CTRL7	5.53E + 07	2 × 150	2.96%	1.77E + 07	31.95%	4
ALS-s1	8.30E + 07	2 × 75	26.69%	3.96E + 07	47.72%	1
ALS-s2	9.25E + 07	2 × 75	27.52%	5.73E + 07	61.92%	1
ALS-s3	9.39E + 07	2 × 75	33.88%	5.23E + 07	55.73%	1
ALS-s4	9.10E + 07	2 × 75	37.71%	2.56E + 07	28.17%	1
ALS-s5	9.18E + 07	2 × 75	31.32%	5.51E + 07	60.00%	2
ALS-s6	8.74E + 07	2 × 75	32.43%	5.29E + 07	60.55%	2
ALS-s7	1.03E + 08	2 × 75	17.58%	7.60E + 07	73.79%	2
ALS-s8	7.73E + 07	2 × 75	28.06%	4.71E + 07	60.97%	2
ALS-s9	7.32E + 07	2 × 75	24.14%	3.69E + 07	50.35%	3
ALS-s10	9.12E + 07	2 × 75	14.76%	5.79E + 07	63.49%	3
ALS-s11	7.83E + 07	2 × 75	34.72%	4.48E + 07	57.26%	3
ALS-s12	4.11E + 07	2 × 150	7.12%	2.12E + 07	51.64%	4
ALS-s13	7.17E + 07	2 × 150	17.87%	3.76E + 07	52.36%	4
ALS-s15	5.36E + 07	2 × 150	10.09%	3.52E + 07	65.63%	4
ALS-s16	8.36E + 07	2 × 150	7.68%	5.53E + 07	66.15%	4
FUS-m1	9.63E + 07	2 × 150	22.49%	6.84E + 07	71.03%	4
FUS-m2	9.82E + 07	2 × 75	13.62%	7.33E + 07	74.65%	2
FUS-m3	1.65E + 06	2 × 150	10.43%	1.28E + 07	81.24%	4
SOD1-m1	1.14E + 08	2 × 75	28.38%	6.86E + 07	60.15%	3
SOD1-m2	8.50E + 07	2 × 75	31.92%	4.88E + 07	57.46%	1
SOD1-m3	3.51E + 07	2 × 75	22.36%	2.33E + 07	66.42%	3
TARDBP- m1	1.16E + 08	2 × 75	31.21%	6.21E + 07	53.46%	3
TARDBP- m2	7.94E + 07	2 × 75	29.40%	3.50E + 07	44.14%	3
VCP-m1	2.98E + 07	2 × 75	13.46%	2.07E + 07	69.66%	1

Table 2. RNA-seq read statistics. All RNA samples used to perform RNA-Seq analysis had RIN = 8 and 260/280 ratio 1.5, as recommended in manufacturer's protocol.

FastQ files were generated via llumina bcl2fastq2 (Version 2.17.1.14 - http://support.illumina.com/downloads/bcl2fastq-conversion-software-v217.html) starting from .bcl files produced by Illumina NextSeq sequencer (see Table 2).

Quality validation and read mapping

Quality of individual sequences were evaluated using FastQC software (see Code Availability 1) after adapter trimming with cutadapt software. Per base sequence quality plots, showing the Phred quality score distribution among all sample reads are shown in Figure 1. Gene and transcript intensities were computed using STAR/RSEM software¹⁹ using the "--strandness forward" option (see Code Availability 2 and 3 and Figure 2). Human genome reference used for the alignment was GRCh38 (Gencode release 27): in a rapidly evolving field like the one of non coding RNA analysis it is indeed fundamental to use up-to-date reference versions containing all the available information about annotated coding and non coding RNAs.

Mapping results are summarized in Table 2. Percentage of uniquely mapped reads is 21.8% on average, with standard deviation 10.35%. Remaining reads belong to ribosomal RNA because of the non perfect ribosomal RNA depletion. Since 10 to 25 M reads are suggested for RNA-Seq experiments²⁰, this dataset is suitable for differential expression analysis, since about 60% of samples (21/36) have at least 10 M reads. Furthermore, biological replicates are the most effective strategy to improve statistical power




Figure 1. Quality assessment FASTQ data. Quality assessment of raw FASTQ sequence data for paired end and left (sample_name_R1) and right reads (sample_name_R2). Box and whisker plots demonstrate the distribution of per base quality for each left and right read position read for each of the analyzed samples. Mean value is indicated by the blue line and the yellow box represents the interquartile range (25–75%) with the lower and upper whiskers represent the 10 and 90% points, respectively. Plots were generated by FastQC program (see Code Availability 1).




Figure 2. Experimental overiview and evaluation of sample variance. (a) The flowchart represents RNA-Seq workflow and data analysis. (b) An estimate of the dispersion parameter for each gene is shown. (c) Principal component analysis results. (d) Heatmap showing the sample-to-sample distance. It was obtained with DeSeq2 package on regularized-logarithm transformed counts. Color code is reported above the heatmap.

in such experiments. For this reason, we decided to share our dataset with the scientific community, to have the possibility of integrating our data with similar data obtained in different laboratories.

Expressed transcripts per sample were evaluated imposing a minimum threshold of 5 counts per gene to consider it as expressed. Figure 3 shows the amount of coding and non coding transcripts detected in each sample. With the exception of 4 samples which have a low number of detected transcripts (probably due to a lower number of available reads), 7000 to 14000 expressed coding genes and 500 to 3000 non coding genes were detected per sample. This is in accordance with the currently available knowledge about non coding transcripts that result to be globally less ex-pressed than coding ones within the cell²¹.

Downstream analysis

Differential expression analysis for all transcripts (coding and non coding RNAs) was performed with the R package DESeq2²² (see Code Availability 4). This tool was selected because of its superior performance in identifying isoforms differential expression²³.




Figure 3. Detected transcripts per sample. Number of detected transcripts per sample for coding (dark grey bars) and lncRNAs (light grey bars), separately. Transcripts were considered if covered by at least 5 reads. Dark grey bars represent coding RNAs while light grey bars represent lncRNAs.

Genes were considered differentially expressed and retained for further analysis with $llog_2$ (disease sample/healthy control) $l \ge 1$ and a FDR ≤ 0.1 . Analysis workflow is summarized in Figure 2.

For each sample, we evaluated the number of detected transcripts for coding and lncRNAs, separately. Transcripts were considered if covered by at least 5 reads²⁴. Results are summarized in Figure 3.

Code availability

The following software and versions were used for quality control and data analysis as described in the main text:

- 1. FastQC v0.11.6 was used for quality assessment of raw FASTQ sequencing data: http://www.bioinformatics.babraham.ac.uk/projects/fastqc/
- STAR 2.6 was used for reads mapping to the human GRCh38 genome assembly: https://github.com/ alexdobin/STAR/releases
- 3. RSEM v1.3.0, was used for gene expression quantification: http://bioinf.wehi.edu.au/featureCounts/
- DESeq2 v1.20.0 was used for differential gene expression analysis: http://bioconductor.org/packages/ DESeq2

Data Records

Raw FASTQ files for the RNA-seq libraries were deposited to the NCBI Sequence Read Archive (SRA), and have been assigned BioProject accession PRJNA416880, Table 1, (Data Citation 1) and BioProject accession PRJNA474387, Table 1 (Data Citation 2). Raw counts data were deposited to the NCBI Gene Expression Omnibus (GEO) with accession number GSE106443, Table 1 (Data Citation 3) and GSE115259, Table 1 (Data Citation 4). First dataset has SRA SRP123453, second dataset has SRA SRP149638.

Technical Validation

RNA integrity assessment

RNAs were quantified using a Nanodrop ND-100 Spectrophotometer (Nanodrop Technologies, Wilmington, USA) and a 2100 Bioanalyzer (Agilent RNA 6000 Nano Kit, Waldbronn, Germany); RNAs with a 260:280 ratio of \geq 1.5 and an RNA integrity number of \geq 8 were subjected to deep sequencing.

RNA-seq data qyality assessment

We applied FastQC v0.11.5 software to verify that per base quality scores are suitable for downstream analysis (Figure 1). Mapping percentage have been computed and are reported in Table 2. PCA plot, distance matrix and dispersion estimates are also shown in Figure 2.

References

- 1. Conesa, A. et al. A survey of best practices for RNA-seq data analysis. Genome Biol. 1, 13 (2016).
- 2. Lee, J. T. Epigenetic regulation by long noncoding RNAs. Science 338, 1435-1439 (2012).
- 3. Lina, M. A., Bajic, V. B. & Zhang, Z. On the classification of long non-coding RNAs. RNA Biol. 10(6), 924-993 (2013).
- 4. Harrow, J. et al. GENCODE: the reference human genome annotation for The ENCODE Project. Genome Res. 22, 1760–1774 (2012).
- 5. Aken, B. et al. The Ensembl gene annotation system. Database 2016, 1–19 (2016).
- Derrien, T. et al. The GENCODE v7 catalog of human long noncoding RNAs: analysis of their gene structure, evolution, and expression. Genome Res. 22, 1775–1789 (2012).
- Cabili, M. N. et al. Integrative annotation of human large intergenic noncoding RNAs reveals global properties and specific subclasses. Gene Dev 25, 1915–1927 (2011).
- 8. Gagliardi, S. et al. Long non-coding and coding RNAs characterization in Peripheral Blood Mononuclear Cells and Spinal Cord from Amyotrophic Lateral Sclerosis patients. Sci. Rep 8, 2378 (2018).

- Nishimoto, Y. et al. The long non-coding RNA nuclear-enriched abundant transcript 1_2 induces paraspeckle formation in the motor neuron during the early phase of amyotrophic lateral sclerosis. Mol. Brain 6, 31 (2013).
- Gagliardi, S., Milani, P., Sardone, V., Pansarasa, O. & Cereda, C. From Transcriptome to Noncoding RNAs: Implications in ALS Mechanism. Neurol. Res. Int 2012, 275–728 (2012).
- 11. Wan, P., Su, W. & Zhuo, Y. The Role of Long Noncoding RNAs in Neurodegenerative Diseases. *Mol. Neurobiol.* **54**(3), 2012–2021 (2017).
- 12. Strong, M. J. The evidence for altered RNA metabolism in amyotrophic lateral sclerosis (ALS). J. Neurol. Sci. 288, 1-12 (2010).
- 13. Milani, P. et al. Posttranscriptional regulation of SOD1 gene expression under oxidative stress: Potential role of ELAV proteins in sporadic ALS. Neurobiol. Dis. 60, 51–60 (2013).
- Pansarasa, O. et al. Lymphoblastoid cell lines as a model to understand amyotrophic lateral sclerosis disease mechanisms. Dis. Model Mec. 11, dmm031625 (2018).
- 15. Gagliardi, S. et al. SOD1 mRNA expression in sporadic amyotrophic lateral sclerosis. Neurobiol. Dis. 39, 198-203 (2010).
- 16. Logroscino, G. et al. Incidence of amyotrophic lateral sclerosis in southern Italy: a population based study. J. Neurol. Neurosur. Ps 76, 1094–1098 (2005).
- 17. Mills, J. et al. Strand-specific RNA-seq provides greater resolution of transcriptome profiling. Curr. Genomics 14, 173-181 (2013).
- 18. Ching, T., Huang, S. & Garmire, L. X. Power analysis and sample size estimation for RNA-Seq differential expression. RNA 20, 1684–1696 (2014).
- Li, B. & Dewey, C. N. RSEM: accurate transcript quantification from RNA-Seq data with or without a reference genome. BMC Bioinformatics 12, 323 (2011).
- 20. Liu, Y. et al. RNA-seq differential expression studies: more sequence or more replication? Bioinformatics 30, 301-304 (2014).
- 21. Pelechano, V. & Lars, M. S. Gene regulation by antisense transcription. Nat. Rev. Genet. 14, 880 (2013).
- 22. Love, M. I., Huber, W. & Anders, S. Moderated estimation of fold change and dispersion for RNA-seq data with DESeq 2. Genome Biol. 15, 550 (2014).
- 23. Carrara, M. et al. Alternative splicing detection workflow needs a careful combination of sample prep and bioinformatics analysis. BMC Bioinformatics 16, S2 (2015).
- 24. Tarazona, S., García-Alcalde, F., Dopazo, J., Ferrer, A. & Conesa, A. Differential expression in RNA-seq: a matter of depth. *Gen. Res* 21(12), 2213–2223 (2011).

Data Citations

- 1. NCBI Sequence Read Archive SRP123453 (2017).
- 2. NCBI Sequence Read Archive SRP149638 (2018).
- 3. Gagliardi, S. et al. Gene Expression Omnibus GSE106443 (2017).
- 4. Zucca, S. et al. Gene Expression Omnibus GSE115259 (2018).

Acknowledgements

This work was supported by grants from the Italian Agency for Research on ALS-ARiSLA (Pilot project "LNCinALS"2015-2016 and Full project "GRANULOPATHY"2015-2017). This work was also supported by "Fondazione Regionale per la Ricerca Biomedica" (FRRB 2015-0023).

Author Contributions

S.Z. and S.G. wrote the manuscript. S.G. and C.P. performed the experiments. S.Z., M.A. and M.O. performed bioinformatic analysis. L.D. and M.C. participated to patients and controls recruitment. M.B. participated to experimental plan. D.S. and O.P. reviewed the manuscript. R.C. and C.C. supervised this work. All authors reviewed and accepted the final version of this manuscript.

Additional Information

Supplementary information accompanies this paper at http://www.nature.com/sdata

Competing interests: The authors declare no competing interests.

How to cite this article: Zucca, S. *et al.* RNA-Seq profiling in Peripheral Blood Mononuclear Cells of Amyotrophic Lateral Sclerosis patients and controls. *Sci. Data.* 6:190006 https://doi.org/10.1038/sdata.2019.6 (2019).

Publisher's note: Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Open Access This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons license, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons license and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this license, visit http://creativecommons.org/licenses/by/4.0/

The Creative Commons Public Domain Dedication waiver http://creativecommons.org/publicdomain/zero/1.0/ applies to the metadata files made available in this article.

© The Author(s) 2019