

AperTO - Archivio Istituzionale Open Access dell'Università di Torino

Toward a definition of blueprint of virgin olive oil by comprehensive two-dimensional gas chromatography

This is the author's manuscript

Original Citation:

Availability:

This version is available <http://hdl.handle.net/2318/144439> since 2016-12-01T13:15:47Z

Published version:

DOI:10.1016/j.chroma.2014.01.067

Terms of use:

Open Access

Anyone can freely access the full text of works made available as "Open Access". Works made available under a Creative Commons license can be used according to the terms and conditions of said license. Use of all other works requires consent of the right holder (author or publisher) if not exempted from copyright protection by the applicable law.

(Article begins on next page)



UNIVERSITÀ DEGLI STUDI DI TORINO

This Accepted Author Manuscript (AAM) is copyrighted and published by Elsevier. It is posted here by agreement between Elsevier and the University of Turin. Changes resulting from the publishing process - such as editing, corrections, structural formatting, and other quality control mechanisms - may not be reflected in this version of the text. The definitive version of the text was subsequently published in [*Journal of Chromatography A*, Volume 1334, 21 March 2014, Pages 101–111, DOI: 10.1016/j.chroma.2014.01.067].

You may download, copy and otherwise use the AAM for non-commercial purposes provided that your license is limited by the following restrictions:

- (1) You may use this AAM for non-commercial purposes only under the terms of the CC-BY-NC-ND license.
- (2) The integrity of the work and identification of the author, copyright owner, and publisher must be preserved in any copy.
- (3) You must attribute this AAM in the following format: Creative Commons BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/deed.en>), <http://dx.doi.org/10.1016/j.chroma.2014.01.067>

**TOWARD A DEFINITION OF BLUEPRINT OF VIRGIN OLIVE OIL BY COMPREHENSIVE TWO-DIMENSIONAL
GAS CHROMATOGRAPHY**

Giorgia Purcaro^{1,*}, Chiara Cordero^{2,*}, Erica Liberto², Carlo Bicchi² and Lanfranco S. Conte¹

¹ Dipartimento di Scienze degli Alimenti, Università di Udine, via Sondrio 2A, I-33100 Udine, Italy

² Dipartimento di Scienza e Tecnologia del Farmaco, Università di Torino, via Pietro Giuria 9, I-10125 Torino, Italy

* Authors for correspondence:

Dr. Giorgia Purcaro - Dipartimento di Scienze degli Alimenti, Università di Udine, via Sondrio 2A, I-33100 Udine, Italy

e-mail: giopurcaro@gmail.com or giorgia.purcaro@uniud.it; phone: +39 0432 558170; fax: +39 0432 558130

Dr. Chiara Cordero - Dipartimento di Scienza e Tecnologia del Farmaco, Università di Torino, Via P. Giuria 9, I-10125 Torino, Italy

e-mail: chiara.cordero@unito.it ; phone: +39 011 6707662; fax: +39 011 2367662

Abstract

This study investigates the applicability of an iterative approach aimed at defining a chemical blueprint of virgin olive oil volatiles to be correlated to the product sensory quality. The investigation strategy proposed allows to fully exploit the informative content of a comprehensive multidimensional gas chromatography (GC × GC) coupled to a mass spectrometry (MS) data set. Olive oil samples (19), including 5 reference standards, obtained from the International Olive Oil Council, and commercial samples, were submitted to a sensory evaluation by a Panel test, before being analysed in two laboratories using different instrumentation, column set, and software elaboration packages in view of a cross-validation of the entire methodology. A first classification of samples based on untargeted *peak features* information, was obtained on raw data from two different column combinations (apolar × polar and polar × apolar) by applying unsupervised Multivariate Analysis (i.e., Principal Component Analysis - PCA). However, to improve effectiveness and specificity of this classification, *peak features* were reliably identified (261 compounds), on the basis of the MS spectrum and linear retention index matching, and subjected to successive pair-wise comparisons based on 2D patterns, which revealed peculiar distribution of chemicals correlated with samples sensory classification. The most informative compounds were thus identified and collected in a “blueprint” of specific defects (or combination of defects) successively adopted to discriminate Extra Virgin from defected oils (i.e. *lampante* oil) with the aid of a supervised approach, i.e., Partial Least Squares-Discriminant Analysis (PLS-DA). In this last step, the principles of sensomics, which assigns higher information potential to analytes with lower odor threshold proved to be successful, and a much more powerful discrimination of samples was obtained in view of a sensory quality assessment.

Keywords: olive oil; sensomics; volatile compounds; comprehensive multidimensional gas chromatography (GCxGC); fingerprinting; pattern recognition; aroma defects

1. Introduction

Modern -omics disciplines investigate sample constituents considered collectively (primary and secondary metabolites, compounds generated by thermal treatments and/or enzymatic activity) [1,2] and represent the route of choice for a comprehensive evaluation of food attributes, in particular if sensory properties (sensomics and flavoromics) are unequivocally correlated to a specific distribution of chemicals [3-5].

Sensomics focuses analytical efforts on revealing sensorially-active compounds (odorants and tastants) in order to link the chemical composition of a food with its sensory quality (flavor) [6]. From this perspective, the analytical process through which the peculiar product profile (i.e. the blueprint) is revealed closely fits the aims of the above mentioned - omics disciplines, in addition to comprehensive information on chemical composition, quality and authenticity,.

However, the number of volatiles effectively contributing to food aroma is relatively small, and the detection, identification, and quantitation of odor-active components, sometimes occurring at trace (ng/kg) levels, require complex procedures that combine effective extraction, powerful separation and highly informative (identity confirmation and mass quantitation) detection.

In this perspective, two-dimensional comprehensive gas chromatography (GC×GC) coupled to mass spectrometry (MS) can be successful in revealing the sensory-related chemical blueprint [7] of high quality foods since it fully exploits the informative potential of each analytical dimension (separation and detection) by (a) increasing the separation power, (b) ordering analytes through logic structures over the chromatographic plane and (c) by enhancing detection sensitivity,. In particular, the possibility of extending the analytical investigation to a broader range of chemicals, if compared to mono-dimensional approaches, improves fingerprinting sensitivity and effectiveness of sample classification [8,9].

Methods for reliable selections of 2D peak information and alignment across set of chromatograms were recently reviewed by Reichenbach *et al.* [10]. In particular, the most performing approaches are those referred to as *peak feature* methods, i.e. methods that enable not only, untargeted analytes alignment (based on peak features, i.e., retention times and MS fragmentation pattern similarity) but also pattern recognition and pair-wise comparisons based on 2D peaks distribution over the chromatographic plane.

In this study, the informative potentials of GC×GC-MS in revealing the chemical blueprint of (Extra) Virgin Olive Oil (EVO and VO) have been investigated by adopting advanced fingerprinting methods for reliable peak selection and focusing on peculiar distributions of sensory-active volatiles (aroma compounds) related to quality attributes currently regulated by international organisms (European Union (EU), International Olive Council (IOC), and the Codex Alimentarius) [11-14].

Interestingly, olive oil is presently the only food product whose sensory attributes are officially regulated, and sensory assessment carried out through well-standardized protocols [11-14] by highly and continuously (permanently) trained panelists. On the basis of the presence (or absence) and the intensities of specific

84 defects and the intensity of the “fruity” perception, evaluated by smelling and tasting, Virgin Olive Oil can
85 be classified into three categories (determining its economic value) namely EVO, VO, and *lampante* oil [15].
86 Undoubtedly, volatile compounds play a crucial role in defining olive oil sensory quality and strong efforts
87 have been dedicated to unravel the composition of this informative fraction, to understand correlations
88 with quality attributes [16-18]. The relative distribution of volatiles depends on several parameters (i.e.
89 cultivar, geographical origin, fruit ripeness, processing practices, and storage [18-27]) and the identification
90 of an unequivocal fingerprint correlated to quality may be a difficult task. Most of these variables concur to
91 determine the intensity and quality of the green and fruity perception, while the presence of defects is
92 mainly due to inappropriate manufacturing practices. Olive storage has to last no longer than 24-48h under
93 appropriate conditions to avoid the formation of three main defects, namely mustiness or mold, winey or
94 vinegary, and fusty. Such defects are mainly due to yeast and fungi metabolic activity and aerobic or
95 anaerobic fermentations, respectively [13]. While inappropriate oil storage conditions after production may
96 lead to muddy and rancid defects, the former may be originated by an excessive contact of the oil with
97 sediment in tanks, and the latter by an auto-oxidation process. The descriptor of the fusty and muddy
98 defects was recently unified due to several misinterpretations [12], although the origin of such defects is
99 very different both from the technological and chemical viewpoint.

100 In the panorama of existing studies, just few of them focused on the identification of markers of the
101 specific defects by adopting reference standard samples provided by the IOC [27, 28]. Morales and co-
102 workers [26], developed a dynamic headspace (D-HS) sampling procedure to extract and concentrate
103 informative volatiles from oil samples, to establish correspondences between specific odorants and related
104 sensory attributes, with known metabolic pathways from endogenous or exogenous (from microorganism)
105 enzymes. The results, normalized and interpreted on the basis of target analytes odor thresholds (OT),
106 placed in a more realistic perspective the sensory significance of specific volatiles suggesting an effective
107 strategy of investigation on this subject.

108 However, despite an extensive knowledge on the volatile composition of olive oil, sensory evaluation has
109 not yet been replaced nor supported by any recognized analytical procedure. This is probably due to
110 several aspects, recently discussed by Aparicio *et al.* [29], among which the complexity of sensory
111 perception, the fuzziness of the semantics of some descriptors, and at the same time, the necessity to
112 improve analytical selectivity and resolution along with the pre-concentration procedure, represent
113 objective limits of current methodologies adopted for volatile profiling of olive oil.

114 In this perspective, the present work aims at evaluating whether GC×GC-MS, in combination with both
115 headspace sampling carried out with Solid Phase Microextraction (HS-SPME) and implemented with
116 advanced data elaboration tools, can be a valuable and more informative analytical approach for a high
117 throughput detailed profiling and an effective fingerprinting of the volatile fraction of olive oil in view of its
118 quality assessment.

119 The chemical blueprint of samples corresponding to regulated defects and authentic EVO samples,
120 provided by IOC, has been adopted as a probe to classify commercial products and to guide an iterative
121 process including sensory evaluation (performed by both an official panel and an internal one), volatiles
122 distribution and their sensory relevance (OT).

123

124 **2. Materials and Methods**

125 **2.1 Reference compounds and solvents**

126 Pure reference compounds for analyte identity confirmation and *n*-alkanes (*n*-C7 to *n*-C30) for Linear
127 Retention Index (I^T_s) determination were from Sigma-Aldrich and Supelco (Milan, Italy). Solvents for *n*-
128 alkanes dilution (toluene and cyclohexane) were HPLC-grade from Riedel-de Haen (Seelze, Germany).

129

130 **2.2 Samples**

131 The IOC supplied five samples of standard defected oils (*muddy, vinegary, mold, fusty, and rancid*); five
132 samples (EVO-1 – EVO-5) were provided by an Official Panel (Stazione Sperimentale per le Industrie degli
133 Oli e dei Grassi, Milano, Italy) after sensory evaluation; 9 samples (EVO-6 – EVO14) were purchased from
134 the local market and evaluated by a trained internal panel (5 people), according to the standardized official
135 method approved by the IOC [14] and reported in Annex XII of ref. [12]. The list of samples analyzed by
136 both research units, along with label descriptions and sensory evaluation is reported in **Table 1**.

137

138 **2.3 SPME devices and validation**

139 The SPME devices and fibers were from Supelco (Bellefonte, PA, USA). A
140 Divinylbenzene/Carboxen/Polydimethylsiloxane (DVB/CAR/PDMS) d_f 50/30 μ m, 2 cm length fiber was
141 chosen and conditioned before use as recommended by the manufacturer.

142 The internal standard (ISTD) loading procedure onto the SPME fiber [30,31] preceded each sampling and
143 was as follows: fiber was inserted into a 20 mL sealed vial containing 500 μ L of ISTD (menthol standard
144 solution at 50 ng/mL in dibutyl phthalate) and exposed to the headspace at 50°C for 40 min. After ISTD
145 loading, the fiber was exposed to the oil sample headspace(2.000 g) at 50°C for 40 min and subsequently
146 introduced into the GC injector for thermal desorption for 10 min at 250 °C. Samples were analyzed in
147 duplicate from the two research units. Precision data referred to intra-laboratory relative standard
148 deviations (RSD %) for retention times of 5% maximum and of 18% on Normalized 2D Peak Volumes (i.e.,
149 Cumulative 2D Peak Area normalized over ISTD) for a selection of target analytes (**Table 2**). Inter-laboratory
150 precision was evaluated on Normalized 2D Peak Volumes and was always within a 25% (RSD) range.

151

152 **2.4. HS-SPME-GC×GC-MS Apolar × Polar set-up**

153 **2.4.1 GC×GC-MS cryogenic modulation optimization**

154 A loop-type single-jet prototype developed at the University of Udine, was adopted as a cryogenic
155 modulator for the apolar \times polar set-up. The heart of the system is a liquid CO₂ valve assembly, equipped
156 with a fixed support to locate easily the loop just at the exit of the CO₂ jet (**Supplementary Figure 1 SF1**),
157 the latter provided by SGE (Australia). The CO₂ jet was periodically activated and deactivated (according to
158 the modulation time set) through the use of compressed air controlled by a 5-way 2 exit solenoid valve
159 (cod. SY3120-5LOU-M5, Clippard, Louvain-La-Neuve, Belgium). The modulation parameters were regulated
160 by a timing controller (H5CX, Omron, Hoofddorp, Netherland).

161 It is well known that some problems related to the focusing of highly volatile components may occur using
162 CO₂ (instead of liquid nitrogen) as cooling fluid. To circumvent such a limitation, a thick film first dimension
163 column (0.50 $\mu\text{m } d_f$) was employed to increase the elution temperature of highly volatile compounds.
164 Furthermore, a 1-m capillary segment, with a thin coating (0.10 $\mu\text{m } d_f$) was used as the modulator loop, to
165 increase the retention capacity during the focusing step. The flow exiting the loop was split between the
166 analytical column and a 1 m \times 0.1 mm i.d. of an uncoated capillary segment, thus diverting a part of the
167 first-column effluent to waste (about 50%). Using such a configuration, two main goals were achieved: (a)
168 closer to optimal linear velocities generated in both dimensions (about 20 and 115 cm/sec in the ¹D and ²D
169 column, respectively) [32], and (b) a better reinjection efficiency onto the ²D column, as recently discussed
170 by Tranchida *et al.* [33]. Furthermore, the first modulation stage was performed at the head of the loop
171 (0.25 mm i.d.), while the second was achieved on the 0.1 mm i.d. analytical column (after the split
172 connection), to further reduce the injected bandwidth.

173 A common orthogonal set, namely apolar \times polar, was employed to simplify the peak recognition process,
174 based on the use of Linear Retention Index (I_s^T) values reported in literature and in some commercially
175 available MS databases. A ± 30 range for I_s^T match was applied to confirm or discard less probable mass
176 spectrum similarity match, except for few compounds (10) for which a ± 50 range was considered
177 acceptable, as reported elsewhere for GC \times GC application [34]. In fact, the presence of the polar ²D column,
178 in terms of retention, is negligible for the less polar analytes, on the contrary, the more polar components,
179 which have intense interactions, can be strongly retained, causing a significantly shifting in the retention
180 time.

181

182 **2.4.2. Instrument set up**

183 At the University of Udine, sampling was performed manually following the SPME procedure described in
184 section 2.3 .

185 The GC \times GC system consisted of two GC2010 gas chromatographs on-line coupled to a QP2010 Ultra
186 quadrupole MS (Shimadzu, Kyoto, Japan) operating in EI mode at 70 eV. The transfer line was set to 250°C.
187 The scan range was set to m/z 40-350 with a scanning rate of 20,000 amu/s, obtaining a spectra generation
188 frequency of 33 Hz.

The system was equipped with a prototype of loop-type dual-stage thermal modulator, laboratory-made as described in section 2.4.1., cooled with liquid CO₂ with a modulation time of 5 s (4.1 s CO₂ flow on and 0.9 s off) adopted for all experiments. The column set was configured as follows: ¹D Rxi-5ms column (similar to 5% phenyl/95% dimethyl polysiloxane)(30 m × 0.25 mm i.d., 0.50 μm *d_f*) from Restek (Bellefonte, USA) (located in GC oven 1) connected to an Rxi-5ms capillary segment (1.0 m × 0.25 mm i.d. × 0.1 μm *d_f*) and then to a Supelcowax-10 (polyethylene glycol) segment (1.2 m × 0.1 mm i.d., 0.10 μm *d_f*) and to a 1 m × 0.1 mm i.d. uncoated capillary (both columns, provided by Supelco, were located in GC2), by using a fixed outlet capillary column splitter (SGE, Ringwood, Australia). A split ratio of about 1:1 occurred at this point, thus a splitless injection was performed to maintain a satisfactory sensitivity of the entire method [35].

Fibers thermo desorption into the GC injector port occurred using the following set-up: split/splitless in splitless mode held for 1 min, injector temperature 250°C. The carrier gas was helium at a constant flow of 0.7 mL/min (initial head pressure 243 KPa). The temperature program was 40°C (2 min) to 200°C at 3°C/min and to 320°C at 15°C/min. An offset of +10°C was applied in the second oven to reduce wrap-around.

Data were acquired by GCMS Solution (ver. 2.6, Shimadzu, Japan) and processed using GC Chromsquare ver. 1.6 (Shimadzu Europe, Duisburg, Germany). FFNSC 1.6 (Chromaleont, Messina, Italy) and NIST08s MS commercial libraries were used for identification. Statistical analysis was performed with SPSS 14.0 (SPSS Inc. Chicago, Illinois, USA).

For *I_T⁷* determination, the *n*-alkane liquid sample solution was injected using an AOC-20i autosampler (Shimadzu, Japan) under the following conditions: split/splitless injector, split mode, split ratio 1/50, injector temperature 280°C, injection volume 1 μL.

209

210 **2.5 HS-SPME-GC×GC-MS Polar × Apolar set-up**

At the University of Torino, volatiles from 2.000 g of oil placed in a 20 mL headspace vial were processed using a HT280T multipurpose sampler (HTA, Brescia, Italy) integrated with an Agilent 6890 GC unit coupled to an Agilent 5975C MS detector (Agilent, Little Falls, DE, USA) operating in EI mode at 70 eV. The GC transfer line was set at 270°C. A *Standard Tune* was used and the scan range was set to *m/z* 40-240 with a scanning rate of 10,000 amu/s, obtaining a spectra generation frequency of 24 Hz.

The system was equipped with a two-stage KT 2004 loop thermal modulator (Zoex Corporation, Houston, TX) cooled with liquid nitrogen and with the hot jet pulse time set at 250 ms with a modulation time of 5 s adopted for all experiments. Fused silica capillary loop dimensions were 1.0 m length and 0.1 m inner diameter. The column set was configured as follows: ¹D SolGel-Wax column (100% polyethylene glycol)(30 m × 0.25 mm i.d., 0.25 μm *d_f*) from SGE Analytical Science (Ringwood, Australia) coupled with a ²D OV1701 column (86% polydimethylsiloxane, 7% phenyl, 7% cyanopropyl) (1 m × 0.1 mm i.d., 0.10 μm *d_f*) from Mega (Legnano, Milan, Italy).

Fibers thermo desorption into the GC injector port occurred using the following set-up: split/splitless in the split mode, split ratio 1:20, injector temperature 250°C. The carrier gas was helium at a constant flow of 0.7 mL/min (initial head pressure 260 KPa). The temperature program was 40°C (1 min) to 180°C at 3°C/min and to 260°C at 20°C/min (5 min).

Data were acquired by Agilent MSD ChemStation ver D.02.00.275 and processed using GC Image GC×GC Software version 2.1b1 (GC Image, LLC Lincoln NE, USA). Statistical analysis was performed with SPSS 14.0 (SPSS Inc. Chicago, Illinois, USA).

For I_s^T determination, the *n*-alkane liquid sample solution was injected using the HT280T sampler (HTA, Brescia, Italy) under the following conditions: split/splitless injector, split mode, split ratio 1/50, injector temperature 280°C, injection volume 1μL.

2.6 Statistical Elaborations and visualization

Principal Component Analysis (PCA), Partial Least Square Discriminant Analysis (PLS-DA) and Orthogonal Partial Least Square Discriminant Analysis (OPLS-DA) were performed with Pirouette software ver. 4.0 (Infometrix, Inc. Bothell, WA, USA). Data were pre-treated by baseline correction through noise subtraction and subsequently pre-processed [36]. Data visualization as heat map was obtained by GENE-E v 3.0.77 (Broad Institute, Inc. Cambridge, MA, USA).

3. Results and Discussion

The analysis strategy started with the untargeted processing of chromatographic data (on raw and pre-processed variables) for sample pre-classification; thus data reduction followed a targeted basis, selecting the most informative analytes whose correlation with sensory defects or technological processing is documented in literature. 2D Patterns obtained with the two-different column combinations adopted for the experiments were successively compared pair-wise to define peculiar distributions of analytes correlated with sensory classification and underline chemical blueprints to be used for a more effective discrimination.

Emphasis will be put on the advantages of the data reduction strategy based on both: (a) structured and detailed 2D patterns provided by GC×GC-MS and (b) higher informative potential of potent odorants within the group of volatiles.

3.2 Untargeted fingerprinting

To fully exploit the informative content of the three dimensional GC×GC-MS data set (1D and 2D retention times and MS fragmentation patterns) the so-called *peak features* methods, have demonstrated

to be effective and, although to a different extent, sensitive enough in monitoring variations in the chemical pattern across complex samples [10].

Among the existing methods, *Comprehensive Template Matching* fingerprinting (CTMF) [34] offers the possibility to fully exploit the data matrix generated by GC×GC-MS, by extending correspondences not only to ¹D and ²D retention, but also to the similarities with the MS fragmentation pattern adopting the NIST MS Search algorithm [37].

In the present study, 2D patterns from oil samples analyzed with the two column combinations (apolar × polar and polar × apolar) in different laboratories, for a total of 76 2D plots (19 sample × 2 replicates × 2 laboratories), were processed by CTMF and reliable *peak features* collected in a cumulative template used to align 2D peaks across sample chromatograms and extract mass quantitative descriptors (Normalized 2D Peak Volumes) for the comparative analysis.

The polar × apolar column set provided 452 *peak features*, while the apolar × polar set gave 395. The list of *peak features* and related information (¹D and ²D retention times, I_s^T in the ¹D, Normalized 2D Peak Volumes and MS fragmentation pattern) are provided as supplementary data (**Supplementary Table ST1**). The resulting Normalized 2D Peak Volumes, referred to the 452 and 395 untargeted analytes, are visualized as a heat map in **Supplementary Figure 2** (**SF2a** apolar × polar and **SF2b** polar × apolar). Variables for the heat map construction (i.e., 2D Peak Volumes) were pre-processed by subtracting row mean and dividing by standard deviation.

The two raw dataset (corresponding to the apolar × polar and polar × apolar) were submitted to Principal Component Analysis (PCA) to see whether samples would be successfully clustered according to their sensory quality (*lampante*, VO and EVO).

As shown in **Figure 1**, considering raw data (without scaling), for both column configurations the variance described by the first three components is respectively of 63% (**Figure 1a** apolar × polar) and 86% (**Figure 1b** polar × apolar). However, the model based on the untargeted approach does not discriminate samples according to their commercial classification: i.e., *lampante* oils (*muddy*, *mold*, *vinegary* and *fusty*) are not well separated from commercial EVOs (EVO-1 to 14). When the PCA is performed on pre-treated variables (baseline correction and auto-scaling) the total variance approaches the 72 % (**Figure 1c** and **1d**) but, again, the related models do not improve discrimination attitudes. This is probably due to a high noise level from variables that are not representative in the description of the sensory quality of samples under study.

The successive steps followed a different strategy: untargeted *peak features* were identified on the basis of MS fragmentation patterns, I_s^T in the ¹D and, where possible, authentic standard confirmation. The list of 261 target analytes is provided as **Supplementary Table 2**, along with their odor descriptors, known correlation with sensory defects or quality (Mo for *mold*; M for *muddy*; F for *fusty*; V for *vinegary*; R for *rancid*; Fr for *fruity*) [27] and the MS matching %.

The resulting 2D patterns were submitted to sequential pair-wise comparisons to reveal peculiar qualitative distributions of informative chemicals to be correlated with the sensory quality of samples.

3.2 Potentials of 2D pattern recognition in revealing the chemical blueprint of olive oil

Volatiles that contribute to the whole aroma of virgin olive oil belong to several classes with a prevalence of polar functionalities: aldehydes, ketones, alcohols, esters, hydrocarbons, furans, lactones, and their relative distribution depends on several, concurrent, parameters. The resulting chemical pattern is quite complex and the identification of an unequivocal sample fingerprint for sensory quality assessment is challenging.

Cultivar, geographical origin, fruit ripeness, processing practices, and storage significantly modify the distribution of volatiles [39], as well as enzymatic activity, such as that of the lipoxygenase pathway (LOX) [40]. In addition, the technological process, and in particular the crushing and malaxation steps, may influence volatiles biosynthesis by promoting/inhibiting different pathways [24,25,41,42].

Sensory-relevant and diagnostic analytes may be hidden, or their reliable identification limited by method LOAV (Limit of Odour Activity Value) [4], a useful parameter recently introduced in sensomics, given by the OT divided by analyte LOQ. By definition, LOAV values above 1 refer to sensitive methods capable of efficient and quantitative assessment of odorants, above their odor threshold, while LOAV below 1 indicates the concentration limit down to which an odorant can be identified but not quantified.

In this context, GC×GC undoubtedly is a strategic choice, able to provide highly detailed, sensitive and structured 2D pattern that acts as a peculiar sample signature for chemical fingerprinting [9,43].

In this perspective, the visual comparison tool of both elaboration software adopted in this study, was used as a support for the 2D pattern recognition and peak finding; each reference defected standard oil from IOC (Mo, M, F, V and R) was pair-wise compared to authentic EVO samples certified by the panel (EVO-1, EVO-5, EVO-8, EVO-10, EVO-11 and EVO-14) to localize relevant discriminating *peak features* across patterns. **Figure 2** reports a schematic summary of the sequential process adopted together with two GC×GC patterns of volatiles. Analytes correlated with specific defects were thus collected in a fingerprint template that was successively applied to the 2D pattern of slightly defected oils (EVO-2, EVO-3, EVO-4, EVO-6, EVO-7, EVO-9, EVO-12 and EVO-13). This last step was implemented to validate the effectiveness of selected analytes/probes, and to see whether or not the prediction of a certain combination of defects was possible on the basis of the chemical fingerprint.

The results show good consistency within the column set in terms of quantitative descriptors (see also section 2.3). Target analytes relative abundance (normalization was done on menthol as ISTD) did not exceed the $\pm 25\%$ of RSD, allowing the cross-validation of the comparative analysis, although, several diagnostic analytes were not detected by the polar \times apolar set; as in the case of ethyl propanoate (25),

325 methyl butanoate (28), (3Z)-hex-3-enal (45), 6-methylhept-5-en-2-ol (137), o-guaiacol (178) and (2Z)-non-2-
 326 enal (217).

327 Looking at fingerprint analytes related to defected and/or to EVO samples (**Supplementary Table 2** -
 328 Sensory defect column header), the pattern recognition strategy proved to be highly reliable: known
 329 marker compounds already reported in literature [21-27,29] were all positively matched and several
 330 others, with a discriminating potential and sensory contribution (very low OTs) could be added.

331 Some alcohols, like butan-1-ol (12), pentan-3-ol (23) and (3Z)-hex-3-en-1-ol (66) were found to be also
 332 present in defected samples (e.g. *vinegary*); carbonyls like pentan-2-one (16), 2-methylbutanal (20), (3E)-
 333 octa-3,5-dien-2-one (185), (3E)-6-methylhepta-3,5-dien-2-one (191) were included, with the last two
 334 targets with very low OTs (0.0005 and 0.38 ng/g respectively), and two esters: ethyl octanoate (230), ethyl
 335 decanoate (256). To the *fusty* fingerprint, (3Z)-hex-3-en-1-ol (66) and propyl acetate (26) were also added
 336 while for *mold*, (2E)-2-methylbut-2-enal (40), (3E)-octa-3,5-dien-2-one (185) and (3E)-6-methylhepta-3,5-
 337 dien-2-one (191) were discriminant.

338 Positive matches and relative Normalized 2D Peak Volumes collected for all the data set and reported in
 339 **Table 2** were submitted to a Partial Least Square Discriminant Analysis (PLS-DA) to prove the effectiveness
 340 of the blueprint analytes selected (for whom experimental OT values in lipophilic media were available in
 341 literature [16, 27]) in predicting the presence of certain quality attributes in a sample,.

342 This algorithm is a PLS regression where the response variable is categorical and indicates samples
 343 classes/categories. Supervised pattern recognition techniques are effective to build a model to be used for
 344 future classifications when classification is known "*a priori*". PLS-DA enables both classification and
 345 discrimination, and compared to other algorithms (SIMCA or KNN) its power rely on the ability to explain
 346 differences between overall class properties and related most informative class analytes.

347 The results are visualized in **Figure 3**, and refer a good separation between defected (arbitrarily categorized
 348 as *class 1*) grouped on the right side of the graph and authentic EVO samples (*class 2*) on the left side. In
 349 fact, according to the IOC classification codes of regulated sensory defects (M, V, Mo and F), PLS-DA results
 350 show that the first three PCs (**Figure 3a**) afford to effectively classify the two classes with the 56.17% of
 351 explained variance, although the chemical pattern of the IOC defected samples differs significantly from
 352 defected EVOs. **Figure 3b** (loading plot) shows the variables that are more correlated to the defects and
 353 *vice versa*.

354 Interestingly, analytes that prevail are those strongly correlated with oxidation and enzymatic activity
 355 (mono and poly-unsaturated aldehydes, esters formed as secondary products of fermentation and some 2
 356 methyl-ketones).

357 The Orthogonal Signal Correction (OPLS) on PLS-DA results enables to isolate within the data matrix, the
 358 most relevant information according to intra and inter-class samples classification.

OPLS-DA results (not shown) indicated that IOC defected samples were highly correlated with hexanal (46) (Mo and V), octane (53) (V), 6-methylhept-5-en-2-one (126) (V, Mo and F), nonanal (196) (M), while for the fusty (F) reference sample and EVOs (3-5-8-10-11-14) other variables referred a better correlation (i.e., 2E,4Z-hept-2,4-dienal (125), 2E-hept-2-enal (110), octan-2-one (134), pentan-3-ol (23), 3z-hex-3-1-ol (66) and octanal (138)).

The model was internally full cross-validated (leave-one-out method) and calculated sensitivity and specificity were of 100%. The model sensitivity refers to the percentage of acceptance related to its own objects while specificity refers to the percentage of the objects of another class rejected by the model. Such an over fitting was therefore expected because of the statistically limited number of samples considered; this limitation would be overcome by appropriate sampling extending the number of informative samples per each class.

However, although the combination of HS-SPME and GC×GC-MS allowed to improve method LOD and LOQs, compared to conventional 1DGC-MS approaches, resulting in an extended pattern of informative volatiles, a further step was necessary to develop a predictive method for quality assessment. In this perspective, the principles of sensomics were implemented to see whether or not it may be successful in this context.

375

3.3 Improving the informative potential of the chemical blueprint: toward sensomics

The sensory quality of each sample not only depends on the presence or absence of specific compounds, but it is related to their odour potency and the resulting concentration in the final product. In this perspective, OT can be used to rank analytes as a function of the sensory potency, although with some limits due to the non-linear dose-response relationship between the perception and the odorant concentration in the sample as demonstrated, and critically discussed, by several authors [45-47].

Therefore, keeping in mind that the pattern of volatiles obtained by HS-SPME sampling reflects, not only analytes concentration, but also their volatility and/or affinity for fiber's coating; 2D Peak Volumes were normalized by OT values from literature. This approach that impacts on the magnitude of variables according with their odour potency, emphasizes the role played in the classification by components with very low OT values or, on the other hand, those with a higher response, reasonably more abundant in the sample.

PLS-DA performed on the resulting data matrix (**Figure 4a** and **4b**) shows a better classification of samples; the first three PCs explain 71.28% of the total and the most informative analytes (**Figure 4b**) for authentic EVOs resulted in: pent-1-en-3-one (14), pent-1-en-3-ol (15), (2E)-hex-2-enal (63) while for defected samples were: butan-1-ol (12), 3-methylbutan-1-ol (31), hexanal (46), octane (53), heptanal (85), (2E)-hept-2-enal (110), oct-1-en-3-ol (131), nonanal (196), (2E)-non-2-enal (218), (2E)-dec-2-enal (245).

393 It is worthy of note that this classification model is mostly driven by markers of defects rather than quality
394 indicators, due to the contribution of IOC reference samples to the defected class. These samples are in
395 fact characterized by a peculiar chemical pattern. Therefore, it can be stated that, although the study of
396 reference defected samples is important to understand any phenomenon, a more “comprehensive”
397 sampling (including slightly defected oils and real-world samples) gives more realistic pictures facilitating
398 the identification of relevant analytes responsible for a particular fingerprint.

399

400 **4. Conclusions and future perspectives**

401 The present study, although not conclusive, emphasizes the analytical advantages of GC×GC-MS in terms of
402 sensitivity, reproducibility (evaluated by cross-validation between laboratories and instrumental set-up)
403 and information potential of the 2D patterns produced, proposing a productive investigation strategy for a
404 reliable quality assessment of EVO and VO oils.

405 Reliable results were obtained on a larger pattern of chemicals in the perspective of samples classification
406 based on sensory quality by combining the effectiveness of the analytical procedure (including volatile
407 extraction by HS-SPME) to provide a large and reliable data set of potentially discriminating analytes and
408 adopting sensomics data interpretation principles (mostly based on sensory active compounds and their
409 related odour potency).

410 The results were validated from two different points of view: (a) by cross-comparison between laboratories
411 working with different instrumentation and software packages and (b) by confirming analytes informative
412 potential with existing studies and available literature. The role played by those chemicals specifically
413 revealed by the presented approach will be the focus of future research aimed at their quantification to
414 obtain reference blueprints of defects. In this perspective the adoption of alternative HCC-HS sampling
415 techniques, connoted by a higher concentration capacity (such as Headspace Sorptive Extraction – HSSE
416 with different polar and/or apolar extraction polymers), which improve fingerprinting sensitivity and enable
417 true quantitation (as for example by Multiple Headspace Extraction – MHE) of key-odorants [43] would be
418 of interest in the perspective of an objective evaluation of EVO oils quality.

419 This approach, applied to a wider and strictly structured experimental design (considering more variables,
420 such as cultivar, geographical origin, different defect intensity etc), supported by more than one official
421 Panel, can allow to robustly and reliably characterize specific markers and related characteristics
422 concentration windows, and thus be used as a valuable tool to support, or even replace, sensory
423 evaluation.

424

425

426

427 **Acknowledgements**

428 University of Udine Unit was financially supported by the Italian Ministry for University and Research
429 (MIUR) with a FIRB “Futuro in Ricerca” project n. RBFR10GSJK while University of Turin Unit was supported
430 by Ricerca finanziata da Università - Fondo per la Ricerca Locale (Ex 60%) Anno 2013.
431
432 *The authors declare no conflict of interest*

433 **References**

- 434 1. S. G. Oliver, M. K. Winson, D. B. Kell, F. Baganz, *Trends Biotechnol.* 16 (1998) 373.
- 435 2. C. H. R. de Vos, Y. Tikunov, A. G. Bovy, R. D. Hall. In I. Blank, M. Wust, C. Yeretian (eds.), *Expression*
- 436 *of Multidisciplinary Flavour Science*, Proceedings of the 12th Weurman Symposium, Zürcher Hochschule für
- 437 *Angewandte and Institut Für Chemie und Biologische Chemie*, Interlaken, Switzerland, 2008, pp. 573–580
- 438 3. M. Herrero, C. Simó, V. García-Cañas, E. Ibáñez, A. Cifuentes, *Mass Spectrom. Rev.* 31(1) (2012) 49.
- 439 4. J. Kiefl, G. Pollner, P. Schieberle, *J. Agric. Food Chem.* 61(22) (2013) 5226
- 440 5. J. Charve, C. Chen, A.D. Hegeman, G.A. Reineccius, *Flavour Frag. J.* 26(6) 2011 429.
- 441 6. W. Grosch, *Chem. Senses* 26 (2001) 533.
- 442 7. M. Christlbauer, P. Schieberle, *J. Agric. Food Chem.* 57 (19) (2009) 9114.
- 443 8. J. Kiefl, C. Cordero, L. Nicolotti, P. Schieberle, S.E. Reichenbach, C. Bicchi. *J. Chromatogr. A.* 1243
- 444 (2012) 81.
- 445 9. C. Cordero, E. Liberto, C. Bicchi, P. Rubiolo, P. Schieberle, S.E. Reichenbach, X. Tian, Q. Tao, *J*
- 446 *Chromatogr. A.* 1217 (2010) 5848.
- 447 10. S.E. Reichenbach, X. Tian, C. Cordero, Q. Tao, *J. Chromatogr. A.* 1226 (2012) 140.
- 448 11. Reg. CE n. 2568/91 Off. J. Eur. Comm. L347 (1991) 69
- 449 12. Reg CE 796/02 Official Journal of European Community L128 2002, 8-28
- 450 13. International Olive Council (IOC). Document COI/T.20/Doc. No 15/Rev.3. Madrid, Spain, 2010
- 451 14. International Olive Council (IOC). Document COI/T.20/Doc. No 15/Rev. 4 November 2011
- 452 15. Reg. CE n 1989/03 Off. J Eur. Union L292 (2003) 57
- 453 16. C.M. Kalua, M.S. Allen, D.R. Jr Bedgood, A.G. Bishop, P.D. Prenzler, K. Robards, *Food Chem.* 100
- 454 (2007) 273.
- 455 17. F. Angerosa, *Eur. J. Lipid Sci. Technol.* 104 (2002) 639.
- 456 18. G. Luna, M.T. Morales, R. Aparicio, *Food Chem.* 98 (2006) 243.
- 457 19. F. Angerosa, M. Servili, R. Selvaggini, A. Taticchi, S. Esposito, G. Montedoro, *J. Chromatogr. A* 1054
- 458 (2004) 17.
- 459 20. M.T. Morales, J.J. Calvente, R. Aparicio, *Flav. Fragr. J.* 11 (1996) 171.
- 460 21. R. Aparicio, M.T. Morales, *J. Agric. Food Chem.* 46 (1998) 1116.
- 461 22. F. Angerosa, R. Mostallino, C. Basti, R. Vito, *Food Chem.* 72 (2001) 19.
- 462 23. M.T. Morales, F. Angerosa, R. Aparicio, *Grasas Aceites.* 50 (1999) 114.
- 463 24. M. Servili, R. Selvaggini, A. Taticchi, S. Esposto, G.F. Montedoro, *J. Am. Oil Chem. Soc.* 80 (2003)
- 464 685.
- 465 25. M. Servili, R. Selvaggini, A. Taticchi, S. Esposto, G.F. Montedoro, *J. Am. Oil Chem. Soc.* 51 (2003)
- 466 7980.

467 26. Servili, M., Taticchi, A., Esposito, S., Urbani, S., Selvaggini, R., & Montedoro, G. J. Agric. Food Chem.
468 56 (2008) 10048.

469 27. M.T. Morales, G. Luna, R. Aparicio, Food. Chem. 91 (2005) 293

470 28. G. Procida, A. Giomo, A. Cichelli, L.S. Conte, J. Sci. Food Agric. 85 (2005) 2175

471 29. R. Aparicio, M.T. Morales, D.L. García-González, Eur. J. Lipid Sci. Technol. 114 (2012) 1114

472 30. Y. Wang, J. O'Reilly, Y. Chen, J. Pawliszyn, J. Chromatogr. A, 1072 (2005) 13.

473 31. L. Setkova, S. Risticvic, C.M. Linton, G. Ouyang, L.M. Bragg, J. Pawliszyn, Anal. Chim. Acta 581
474 (2007) 221.

475 32. G. Purcaro, P.Q. Tranchida, P. Dugo, E. La Camera, G. Bisignano, L. Conte, L. Mondello, J. Sep. Sci.
476 33 (2010) 2334.

477 33. P.Q. Tranchida, M. Zoccali, F.A. Franchina, P. Dugo, L. Mondello, J. Sep. Sci. 36 (2013) 212

478 34. L. Mondello, A. Casilli, P.Q. Tranchida, G. Dugo, P. Dugo, J. Chromatogr. A, 1067 (2005) 235-243.

479 35. P.Q. Tranchida, G. Purcaro, C. Fanali, P. Dugo, G. Dugo, L. Mondello, J. Chromatogr. A, 1217 (2010)
480 4160.

481 36. R. Beebe, R.J. Pell, M. B. Seasholtz, Chemometrics a practical guide (1998) John Wiley & Sons, Inc.
482 Chichester, UK

483 37. C. Cordero, E. Liberto, C. Bicchi, P. Rubiolo, S.E. Reichenbach, X. Tian, Q. Tao, J. Chromatogr. Sci. 48
484 (2010) 251.

485 38. S. Stein NIST Mass Spectral Search Program (Version 2.0f). NIST Mass Spectrometry Data Center,
486 2008.

487 39. F. Angerosa, M. Servili, R. Selvaggini, A. Taticchi, S. Esposito, G. Montedoro, J. Chromatogr. A. 1054
488 (2004) 17.

489 40. F. Angerosa, C. Basti, R. Vito, J. Agric. Food Chem. 47 (1999) 836.

490 41. F. Angerosa, R. Mostallino, C. Basti, R. Vito, Food Chem. 72 (2001) 19.

491 42. M.T. Morales, F. Angerosa, R. Aparicio, Grasas Aceites. 50 (1999) 114.

492 43. C. Cordero, C. Cagliero, E. Liberto, L. Nicolotti, P. Rubiolo, B. Sgorbini, C. Bicchi J. Chromat. A. 1318
493 (2013) 1.

494 44. M. Servili, A. Taticchi, S. Esposito, S. Urbani, R. Selvaggini, G. Montedoro, J. Agric. Food Chem. 56
495 (2008) 10048.

496 45. V. Audouin, F. Bonnet, Z.M. Vickers, G.A. Reineccius ACS Symposium Series Volume 782 (2001) 156

497 46. J.C. Stevens, W.S. Cain, R.J. Burke. Chem Senses. 13(4) (1988) 643.

498 47. J.E.R. Frijters. Chemical Senses and Flavor. 3(2) (1978) 227.

499

500

501 **Supplementary data**

502 **Supplementary Figure 1.** Loop-type single-jet prototype developed at the University of Udine. For details
503 see text.

504 **Supplementary Figure 2.** Heat map based on untargeted *peak features* detected by the apolar × polar
505 (ST1a) and the polar × apolar (ST1b) column set. Absolute 2D Peak Volumes are normalized across each
506 sample set by subtracting row mean and dividing by standard deviation. *Peak features* along columns
507 follow the retention time ordering.

508 **Supplementary Table 1** List of *peak features* and related information (¹D and ²D retention times, I_s^T in the
509 ¹D, Normalized 2D Peak Volumes and MS fragmentation pattern).

510 **Supplementary Table 2** List of identified compounds in the apolar × polar set, along with ¹D and ²D
511 retention times, experimental and library reported I_s^T , known correlation with sensory defects, odor
512 descriptors and MS similarity matching (%). Mo: *mold*; M: *muddy*; F: *fusty*; V: *vinegary*; R: *rancid*; Fr: *fruity*.
513 Defects marked with * refer to correlation highlighted using the proposed sequential process.

514

515 **Figure captions**

516 **Figure 1.** PCA score plots performed on untargeted *peak features* resulted from GC×GC-MS carried out on
517 the apolar x polar (1a and 1c) and polar x apolar column set (1b and 1d). Raw data visualization 1a)
518 explained variance 63%, 1b) explained variance 86.2%. PCA on pre-treated (baseline correction) and auto-
519 scaled data 1c) explained variance 71.95% and 1d) explained variance 72.56%. For data interpretations see
520 text.

521 **Figure 2.** Scheme of the sequential process applied in this study.

522 **Figure 3.** Internally cross-validated PLS-DA model and resulting score (3a) and loadings (3b) plots
523 performed on the Normalized 2D peak volumes of the compounds listed in Table 2. Explained variance on
524 the first three PCs is 56.17%; baseline correction and pareto pre-processing were performed.

525 **Figure 4.** Internally cross-validated PLS-DA model and resulting score (4a) and loadings (4b) plots
526 performed on the Normalized 2D peak volumes divided by the OTs on compounds listed in Table 2.
527 Explained variance on the first three PCs is 71.28%, baseline correction and pareto pre-processing were
528 performed.

529

530 **Table captions**

531 **Table 1.** Description of the olive oil samples.

532 **Table 2.** List of compounds adopted in the chemical blueprint of samples. ¹D and ²D retention times,
533 experimental I_s^T , known correlation with sensory defect, odor threshold values (OT in mg/kg) and
534 normalized peak volume (multiplied per 10⁵) are reported for both apolar × polar and polar × apolar (in
535 *italic*) column set.

Table 1. Description of the olive oil samples.

| Sample Name | Type | Sensory Analysis | | Additional label information |
|-----------------|--|------------------|------------------------|---|
| | | Panel | Descriptors | |
| <i>Musty</i> | Lampante oil | Official | Intense Musty | IOC standard |
| <i>Vinegary</i> | Lampante oil | Official | Intense Vinegary | IOC standard |
| <i>Fusty</i> | Lampante oil | Official | Intense Fusty | IOC standard |
| <i>Mould</i> | Lampante oil | Official | Intense Mould | IOC standard |
| <i>Rancid</i> | Lampante oil | Official | Moderate Rancid | IOC standard |
| EVO-1 | Commercial EVO (small production) | Official | EVO | Declared to be produced with Italian olives |
| EVO-2 | Commercial EVO (small production) | Official | Mould and rancid | Declared to be produced with Italian olives |
| EVO-3 | Commercial EVO (small production) | Official | Fusty and muddy | Declared to be produced with Spanish olives |
| EVO-4 | Commercial EVO (small production) | Official | Fusty, muddy, vinegary | Declared to be produced with Spanish olives |
| EVO-5 | Commercial EVO (small production) | Official | EVO | Declared to be produced with Greece olives |
| EVO-6 | Commercial EVO (industrial production) | Internal | Fusty/muddy | Declared to be produced with EU olives |
| EVO-7 | Commercial EVO (industrial production) | Internal | Fusty/muddy | Declared to be produced with EU olives |
| EVO-8 | Commercial EVO (industrial production) | Internal | EVO | Declared to be produced with EU olives |
| EVO-9 | Commercial EVO (industrial production) | Internal | Slightly fusty | Declared to be produced with Italian olives |
| EVO-10 | Commercial EVO (industrial production) | Internal | EVO | Declared to be produced with Italian olives |
| EVO-11 | Commercial EVO (industrial production) | Internal | EVO | Declared to be produced with EU olives |
| EVO-12 | Commercial EVO (industrial production) | Internal | Intense rancid | Declared to be produced with Italian olives |
| EVO-13 | Commercial EVO (industrial production) | Internal | EVO/slightly rancid | Declared to be produced with Italian olives |
| EVO-14 | Commercial EVO (small production) | Internal | EVO | Declared to be produced with Slovenian olives |

Table 2. List of compounds adopted in the chemical blueprint of samples. ¹D and ²D retention times, experimental I_s^T , known correlation with sensory defect, odor threshold values (OT in mg/kg) and normalized peak volume (multiplied per 10⁵) are reported for both apolar × polar and polar × apolar (in *italic*) column set.

| ID# | | ¹ D (min) | ² D (s) | I ¹ _s | Sensory defect | OT (mg/kg) | Muddy | Vinegary | Mould | Fusty | Rancid | EVO-01 | EVO-02 | EVO-03 | EVO-04 | EVO-05 | EVO-06 | EVO-07 | EVO-08 | EVO-09 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|--------------------|-------------------------|-----------------------|-----------------------------|-------------------|---------------|-------|----------|-------|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 1 | ethanol | 2.78 | 1.72 | 552 | V | 30 | 3.7 | 1.0 | 2.8 | 4.3 | 6.3 | 2.1 | 6.2 | 0.0 | 114.9 | 2.6 | 79.0 | 3.4 | 29.2 | 24.1 | 18.5 | 12.3 | 29.7 | 8.5 | 5.5 |
| | | 7.67 | 1.14 | 883 | | | 4.7 | 1.3 | 3.3 | 4.0 | 6.3 | 2.8 | 5.6 | 0.0 | 65.0 | 4.0 | 75.7 | 3.1 | 37.0 | 14.9 | 18.0 | 10.0 | 20.0 | 8.8 | 6.2 |
| 6 | butanal | 6.85 | 1.05 | 607 | F/M | 0.018 | 11.3 | 0.0 | 0.3 | 2.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 7.00 | 2.48 | 857 | | | 14.2 | 0.0 | 0.4 | 3.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 7 | ethyl acetate | 7.52 | 0.96 | 613 | F/V | 0.94 | 18.9 | 17.2 | 5.3 | 4.7 | 0.0 | 0.0 | 0.0 | 9.4 | 16.9 | 0.0 | 85.7 | 40.2 | 12.4 | 10.7 | 8.5 | 2.8 | 0.0 | 3.2 | 0.0 |
| | | 6.75 | 1.35 | 850 | | | 23.8 | 18.9 | 6.2 | 4.3 | 0.0 | 0.0 | 0.0 | 9.2 | 9.6 | 0.0 | 82.1 | 37.1 | 15.8 | 6.6 | 8.3 | 2.3 | 0.0 | 3.2 | 0.0 |
| 9 | acetic acid | 8.11 | 1.66 | 622 | F/V/R | 0.5 | 12.5 | 23.7 | 18.4 | 32.9 | 12.6 | 4.6 | 12.6 | 21.0 | 194.2 | 5.6 | 3.6 | 37.0 | 13.5 | 5.0 | 2.9 | 6.7 | 0.7 | 1.1 | 8.1 |
| | | 28.50 | 0.97 | 1457 | | | 15.7 | 26.0 | 21.5 | 30.7 | 12.6 | 6.2 | 11.2 | 20.8 | 109.9 | 8.6 | 3.4 | 34.1 | 17.1 | 5.0 | 2.9 | 7.3 | 0.7 | 1.0 | 6.8 |
| 12 | butan-1-ol | 9.47 | 2.92 | 648 | V/M * | 0.15 | 1.6 | 4.0 | 0.0 | 17.1 | 0.0 | 0.3 | 1.3 | 1.8 | 0.6 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.2 | 0.0 |
| | | 14.42 | 1.26 | 1120 | | | 1.5 | 3.6 | 0.0 | 16.0 | 0.0 | 0.4 | 1.1 | 1.8 | 0.4 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.2 | 0.0 |
| 13 | 3-methylbutanal | 9.52 | 0.96 | 649 | F/V | 0.013 | 0.0 | 1.2 | 0.0 | 0.1 | 0.0 | 0.0 | 0.0 | 0.6 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 8.75 | 2.61 | 927 | | | 0.0 | 1.5 | 0.0 | 0.1 | 0.0 | 0.0 | 0.0 | 0.6 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 14 | pent-1-en-3-one | 10.61 | 1.39 | 668 | M | 0.00073 | 0.7 | 1.6 | 0.6 | 2.4 | 1.4 | 29.5 | 2.0 | 6.1 | 0.4 | 4.6 | 0.9 | 24.3 | 1.5 | 0.0 | 3.4 | 2.8 | 0.3 | 0.5 | 8.9 |
| | | 10.17 | 1.47 | 983 | | | 0.9 | 1.8 | 0.8 | 2.2 | 1.4 | 39.6 | 1.8 | 6.0 | 0.2 | 7.1 | 0.9 | 22.4 | 1.9 | 0.0 | 3.3 | 2.3 | 0.2 | 0.6 | 10.1 |
| 15 | pent-1-en-3-ol | 10.71 | 2.77 | 670 | F/Mo | 0.4 | 1.6 | 6.8 | 1.2 | 11.3 | 11.6 | 18.0 | 8.4 | 7.5 | 9.6 | 15.1 | 6.0 | 8.3 | 9.4 | 11.7 | 11.8 | 10.8 | 7.4 | 5.3 | 4.5 |
| | | 15.08 | 1.30 | 1135 | | | 2.0 | 7.5 | 1.4 | 10.6 | 11.6 | 24.2 | 7.5 | 7.4 | 8.2 | 20.2 | 5.8 | 7.6 | 11.9 | 7.3 | 11.5 | 8.8 | 6.1 | 5.4 | 5.1 |
| 16 | pentan-2-one | 10.77 | 1.18 | 671 | V | 70 | 0.5 | 1.8 | 0.5 | 0.9 | 0.0 | 0.0 | 0.5 | 1.0 | 0.0 | 1.8 | 0.5 | 0.0 | 12.3 | 23.8 | 0.0 | 0.0 | 0.4 | 1.9 | 0.8 |
| | | 8.84 | 1.47 | 930 | | | 0.7 | 1.9 | 0.6 | 0.8 | 0.0 | 0.0 | 0.4 | 0.9 | 0.0 | 2.4 | 0.5 | 0.0 | 15.6 | 14.8 | 0.0 | 0.0 | 0.3 | 1.9 | 0.9 |
| 20 | 2-methylbutanal | 11.02 | 1.15 | 676 | V/M * | 0.01 | 3.8 | 16.0 | 7.5 | 9.3 | 0.7 | 16.7 | 20.5 | 18.2 | 5.1 | 0.0 | 16.3 | 0.0 | 0.0 | 0.0 | 21.1 | 18.9 | 8.4 | 19.2 | 0.0 |
| | | 4.34 | 1.09 | 750 | | | 4.7 | 17.6 | 8.7 | 8.7 | 0.7 | 22.4 | 18.3 | 18.0 | 4.4 | 0.0 | 15.6 | 0.0 | 0.0 | 0.0 | 20.5 | 15.3 | 6.9 | 19.8 | 0.0 |
| 23 | pentan-3-ol | 11.75 | 4.97 | 689 | V | 9 | 0.5 | 14.0 | 0.0 | 0.0 | 0.0 | 0.8 | 0.0 | 0.0 | 1.2 | 0.0 | 1.8 | 0.8 | 0.8 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 13.00 | 1.35 | 1078 | | | 0.6 | 11.0 | 0.0 | 0.0 | 0.0 | 1.1 | 0.0 | 0.0 | 1.0 | 0.0 | 1.7 | 0.7 | 1.0 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 25 | ethyl propanoate | 12.10 | 0.93 | 696 | F | 0.1 | 4.8 | 0.4 | 0.4 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 28 | methyl butanoate | 12.43 | 1.05 | 702 | F | 0.06 | 0.8 | 0.0 | 0.0 | 4.4 | 0.0 | 0.0 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 31 | 3-methylbutan-1-ol | 13.46 | 2.62 | 720 | F/M/Mo | 0.1 | 2.4 | 8.7 | 1.0 | 2.9 | 0.0 | 1.0 | 5.3 | 18.1 | 23.1 | 3.1 | 15.5 | 11.5 | 4.9 | 8.0 | 3.3 | 4.1 | 2.8 | 3.6 | 1.2 |
| | | 18.35 | 0.94 | 1215 | | | 3.0 | 9.6 | 1.1 | 2.7 | 0.0 | 1.4 | 4.7 | 17.9 | 19.8 | 4.1 | 14.9 | 10.6 | 6.2 | 4.9 | 3.2 | 3.3 | 2.3 | 3.7 | 1.4 |
| 32 | (2E)-pent-2-enal | 13.53 | 1.93 | 721 | V | 0.3 | 5.9 | 12.9 | 6.9 | 21.4 | 10.5 | 4.4 | 3.4 | 5.7 | 1.2 | 4.5 | 1.6 | 0.0 | 1.2 | 0.5 | 1.4 | 1.1 | 2.9 | 1.3 | 1.0 |
| | | 14.08 | 1.52 | 1110 | | | 7.4 | 14.2 | 8.1 | 20.0 | 10.5 | 6.0 | 3.0 | 5.6 | 1.0 | 6.0 | 1.5 | 0.0 | 1.5 | 0.3 | 1.4 | 0.9 | 2.4 | 1.3 | 1.1 |
| 39 | pentan-1-ol | 14.88 | 2.92 | 746 | F/M/V | 3 | 1.0 | 2.1 | 0.9 | 4.5 | 1.9 | 1.0 | 3.5 | 2.1 | 2.0 | 3.7 | 4.0 | 5.8 | 1.0 | 6.0 | 1.3 | 1.4 | 3.4 | 7.4 | 0.0 |
| | | 19.00 | 1.35 | 1231 | | | 1.0 | 2.1 | 0.9 | 4.5 | 1.9 | 1.0 | 3.5 | 2.1 | 2.0 | 3.7 | 4.0 | 5.8 | 1.0 | 6.0 | 1.3 | 1.4 | 3.4 | 7.4 | 0.0 |
| 45 | (3Z)-hex-3-enal | 15.86 | 1.66 | 764 | Fr | 0.003 | 0.0 | 3.2 | 1.3 | 4.7 | 11.8 | 14.2 | 7.6 | 2.7 | 0.0 | 5.3 | 1.5 | 1.7 | 3.4 | 0.7 | 1.8 | 3.9 | 28.9 | 15.8 | 9.4 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 46 | Hexanal | 15.94 | 1.33 | 773 | F/Mo/V/R | 0.08 | 92.8 | 170.7 | 209.1 | 7.6 | 59.8 | 20.6 | 37.1 | 44.7 | 14.7 | 22.2 | 15.8 | 16.3 | 25.1 | 32.9 | 29.4 | 33.4 | 212.8 | 37.0 | 13.0 |
| | | 12.25 | 1.77 | 1054 | | | 116.7 | 187.7 | 244.7 | 7.1 | 59.8 | 27.7 | 33.1 | 44.2 | 12.6 | 29.6 | 15.2 | 15.0 | 31.9 | 20.4 | 28.6 | 27.1 | 175.8 | 38.0 | 14.8 |
| 48 | ethyl butanoate | 17.10 | 0.93 | 786 | F | 0.03 | 5.2 | 0.7 | 0.0 | 31.3 | 0.0 | 0.0 | 0.0 | 1.6 | 1.6 | 0.0 | 0.0 | 0.0 | 0.3 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 10.59 | 1.77 | 1000 | | | 6.6 | 0.8 | 0.0 | 29.2 | 0.0 | 0.0 | 0.0 | 1.6 | 1.4 | 0.0 | 0.0 | 0.0 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |

Table 3. continued.

| ID# | | ¹ D (min) | ² D (s) | ¹ T _s | Sensory defect | OT (mg/kg) | Muddy | Vinegary | Mould | Fusty | Rancid | EVO-01 | EVO-02 | EVO-03 | EVO-04 | EVO-05 | EVO-06 | EVO-07 | EVO-08 | EVO-09 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|--------------------------|-------------------------|-----------------------|-----------------------------|-------------------|---------------|-------|----------|-------|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 49 | oct-1-ene | 17.18 | 0.54 | 787 | M | 0.08 | 3.9 | 1.6 | 2.7 | 1.5 | 0.0 | 0.0 | 3.9 | 4.8 | 7.9 | 0.7 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 8.4 | 0.0 | 0.0 |
| | | 6.09 | 1.68 | 820 | | | 4.9 | 1.8 | 3.2 | 1.4 | 0.0 | 0.0 | 3.5 | 4.7 | 6.7 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 6.9 | 0.0 | 0.0 | |
| 52 | butyl acetate | 17.60 | 1.05 | 795 | F | 0.3 | 0.0 | 0.0 | 0.0 | 9.3 | 0.0 | 0.0 | 0.3 | 1.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 12.00 | 1.73 | 1046 | | | 0.0 | 0.0 | 0.0 | 8.7 | 0.0 | 0.0 | 0.2 | 1.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 53 | octane | 17.68 | 0.51 | 800 | F/V/R | 0.94 | 51.9 | 115.2 | 65.3 | 75.4 | 0.4 | 0.9 | 21.9 | 80.2 | 51.1 | 9.3 | 24.1 | 31.8 | 8.6 | 7.2 | 7.5 | 12.2 | 20.4 | 42.7 | 0.0 |
| | | 5.59 | 1.89 | 800 | | | 65.3 | 126.7 | 76.4 | 70.3 | 0.4 | 1.2 | 19.5 | 79.2 | 43.8 | 12.4 | 23.1 | 29.3 | 10.9 | 4.5 | 7.3 | 9.9 | 16.9 | 43.9 | 0.0 |
| 63 | (2E)-hex-2-enal | 19.11 | 1.84 | 822 | Mo/V/F/R | 0.42 | 1.0 | 1.7 | 1.9 | 8.2 | 565.4 | 582.0 | 396.0 | 37.0 | 10.4 | 209.1 | 218.4 | 209.0 | 545.4 | 444.5 | 565.9 | 425.7 | 406.8 | 391.3 | 782.0 |
| | | 18.00 | 1.64 | 1208 | | | 1.2 | 1.8 | 2.3 | 7.7 | 565.4 | 781.4 | 353.0 | 36.5 | 8.9 | 279.5 | 209.2 | 192.6 | 692.1 | 275.2 | 551.0 | 345.5 | 336.1 | 402.5 | 886.9 |
| 64 | (3E)-hex-3-en-1-ol | 19.31 | 3.69 | 826 | F/R/V | 6 | 0.7 | 0.4 | 0.6 | 0.8 | 1.6 | 0.6 | 0.7 | 1.4 | 2.6 | 2.6 | 3.7 | 2.6 | 1.5 | 2.1 | 1.3 | 3.0 | 1.9 | 5.9 | 2.2 |
| | | 24.08 | 1.30 | 1350 | | | 0.9 | 0.5 | 0.6 | 0.7 | 1.6 | 0.8 | 0.6 | 1.3 | 2.2 | 3.5 | 3.6 | 2.4 | 2.0 | 1.3 | 1.3 | 2.4 | 1.6 | 6.1 | 2.5 |
| 66 | (3Z)-hex-3-en-1-ol | 19.55 | 3.25 | 830 | V/F | 1.5 | 1.4 | 0.6 | 0.0 | 15.8 | 20.0 | 9.3 | 1.2 | 73.1 | 75.8 | 58.2 | 71.8 | 83.2 | 34.5 | 52.1 | 49.8 | 42.7 | 5.1 | 6.8 | 18.4 |
| | | 24.92 | 1.35 | 1369 | | | 1.8 | 0.6 | 0.0 | 14.7 | 20.0 | 12.5 | 1.1 | 72.1 | 64.9 | 77.8 | 68.8 | 76.7 | 43.8 | 32.2 | 48.5 | 34.6 | 4.2 | 7.0 | 20.8 |
| 68 | (2E)-hex-2-en-1-ol | 19.89 | 3.58 | 837 | V | 5 | 0.9 | 3.9 | 0.0 | 0.0 | 49.0 | 55.8 | 115.9 | 9.8 | 26.1 | 80.3 | 61.0 | 87.9 | 54.9 | 81.7 | 33.4 | 65.3 | 14.8 | 194.5 | 19.1 |
| | | 25.92 | 1.26 | 1392 | | | 1.1 | 4.3 | 0.0 | 0.0 | 49.0 | 75.0 | 103.3 | 9.7 | 22.4 | 107.3 | 58.5 | 81.0 | 69.6 | 50.6 | 32.6 | 53.0 | 12.2 | 200.0 | 21.6 |
| 74 | hexan-1-ol | 20.63 | 2.71 | 850 | Fr | 0.4 | 1.3 | 4.6 | 1.3 | 19.2 | 0.0 | 25.1 | 33.6 | 23.8 | 40.3 | 34.9 | 59.4 | 98.5 | 54.0 | 82.0 | 39.3 | 54.1 | 11.0 | 79.5 | 17.5 |
| | | 18.17 | 1.81 | 1212 | | | 1.7 | 5.0 | 1.5 | 17.9 | 0.0 | 33.7 | 30.0 | 23.5 | 34.6 | 46.6 | 56.9 | 90.8 | 68.6 | 50.8 | 38.2 | 43.9 | 9.0 | 81.8 | 19.8 |
| 80 | heptan-2-one | 21.85 | 1.27 | 872 | V | 0.3 | 3.8 | 10.0 | 7.6 | 7.5 | 0.5 | 0.3 | 0.7 | 2.5 | 0.9 | 3.4 | 2.7 | 1.6 | 0.6 | 0.6 | 0.0 | 1.4 | 2.1 | 0.0 | 0.0 |
| | | 16.42 | 1.89 | 1169 | | | 4.8 | 10.9 | 8.9 | 7.0 | 0.5 | 0.4 | 0.6 | 2.5 | 0.8 | 4.5 | 2.6 | 1.5 | 0.8 | 0.4 | 0.0 | 1.1 | 1.7 | 0.0 | 0.0 |
| 85 | heptanal | 22.35 | 1.24 | 881 | R | 0.5 | 10.8 | 6.9 | 16.7 | 7.3 | 1.6 | 0.3 | 3.5 | 5.6 | 1.6 | 1.1 | 2.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 8.0 | 2.0 | 0.0 |
| | | 16.50 | 1.89 | 1161 | | | 13.6 | 7.6 | 19.5 | 6.8 | 1.6 | 0.4 | 3.1 | 5.5 | 1.4 | 1.5 | 2.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 6.6 | 2.0 | 0.0 |
| 86 | ethyl pentanoate | 22.68 | 0.96 | 887 | M | 0.0015 | 1.3 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.8 | 1.3 | 0.0 |
| | | 23.08 | 2.27 | 1327 | | | 1.7 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 1.3 | 0.0 |
| 103 | (2Z)-hept-2-enal | 24.28 | 1.72 | 917 | R | 0.042 | 0.7 | 0.8 | 1.2 | 2.4 | 0.5 | 0.0 | 0.2 | 0.0 | 0.0 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.1 | 1.3 | 3.1 | 0.0 |
| | | 22.58 | 1.77 | 1315 | | | 0.8 | 0.9 | 1.4 | 2.2 | 0.5 | 0.0 | 0.2 | 0.0 | 0.0 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.7 | 1.0 | 3.1 | 0.0 |
| 110 | (2E)-hept-2-enal | 24.78 | 1.87 | 926 | Mo/R | 0.005 | 14.8 | 14.0 | 39.6 | 2.8 | 12.0 | 2.3 | 9.9 | 12.3 | 4.6 | 4.8 | 3.6 | 3.8 | 0.9 | 3.5 | 2.4 | 0.0 | 21.1 | 0.0 | 0.4 |
| | | 22.67 | 1.81 | 1317 | | | 18.6 | 15.5 | 46.3 | 2.6 | 12.0 | 3.1 | 8.8 | 12.2 | 4.0 | 6.4 | 3.5 | 3.5 | 1.2 | 2.2 | 2.3 | 0.0 | 17.4 | 0.0 | 0.5 |
| 122 | (2Z,4E)-hepta-2,4-dienal | 26.30 | 2.98 | 955 | R/Mo/F | 0.36 | 0.7 | 1.7 | 0.6 | 0.7 | 1.2 | 0.0 | 2.7 | 11.6 | 0.0 | 0.0 | 4.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.7 | 0.0 |
| | | 28.58 | 1.60 | 1455 | | | 0.9 | 1.9 | 0.8 | 0.7 | 1.5 | 0.0 | 2.4 | 11.5 | 0.0 | 0.0 | 4.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.8 | 0.0 |
| 123 | oct-1-en-3-one | 26.44 | 1.42 | 958 | Mo | 0.01 | 13.5 | 11.0 | 37.8 | 27.1 | 9.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 14.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 21.58 | 1.89 | 1292 | | | 16.9 | 12.1 | 44.3 | 25.3 | 9.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 8.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 125 | (2E,4Z)-hepta-2,4-dienal | 26.63 | 2.77 | 961 | R/Mo/F | 10 | 6.6 | 15.7 | 8.3 | 24.4 | 25.7 | 4.2 | 0.3 | 0.6 | 2.5 | 4.0 | 2.6 | 4.3 | 1.0 | 4.3 | 4.3 | 2.9 | 5.5 | 0.0 | 0.0 |
| | | 28.66 | 3.24 | 1457 | | | 8.3 | 17.2 | 9.7 | 22.8 | 25.7 | 5.7 | 0.3 | 0.6 | 2.2 | 5.3 | 2.5 | 4.0 | 1.3 | 2.6 | 4.2 | 2.4 | 4.6 | 0.0 | 0.0 |
| 126 | 6-methylhept-5-en-2-one | 26.69 | 1.63 | 962 | Mo/F/R | 1 | 148.2 | 88.6 | 88.6 | 70.7 | 2.1 | 0.8 | 4.7 | 11.0 | 4.4 | 4.0 | 2.9 | 3.0 | 0.0 | 0.0 | 0.0 | 1.9 | 4.5 | 2.0 | 0.0 |
| | | 23.17 | 1.85 | 1329 | | | 186.3 | 97.5 | 103.7 | 66.0 | 2.1 | 1.1 | 4.2 | 10.8 | 3.7 | 5.3 | 2.8 | 2.8 | 0.0 | 0.0 | 0.0 | 1.5 | 3.7 | 2.1 | 0.0 |
| 131 | oct-1-en-3-ol | 27.04 | 2.26 | 969 | Mo | 0.05 | 4.0 | 3.1 | 11.2 | 26.2 | 2.6 | 0.5 | 2.1 | 2.0 | 1.2 | 1.3 | 1.6 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 1.2 | 2.0 | 0.0 |
| | | 27.83 | 1.47 | 1437 | | | 0.6 | 0.3 | 1.0 | 1.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 134 | octan-2-one | 27.35 | 1.27 | 975 | V | 0.51 | 4.4 | 4.9 | 7.3 | 5.9 | 0.5 | 0.4 | 0.6 | 1.0 | 1.6 | 3.6 | 2.6 | 0.9 | 0.0 | 0.9 | 1.7 | 0.8 | 1.3 | 0.5 | 0.0 |
| | | 20.83 | 2.06 | 1274 | | | 0.0 | 0.0 | 0.1 | 0.2 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |

Table 3. continued.

| ID# | | ¹ D (min) | ² D (s) | I ^T _s | Sensory defect | OT (mg/kg) | Muddy | Vinegary | Mould | Fusty | Rancid | EVO-01 | EVO-02 | EVO-03 | EVO-04 | EVO-05 | EVO-06 | EVO-07 | EVO-08 | EVO-09 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|-----------------------------------|-------------------------|-----------------------|-----------------------------|-------------------|---------------|-------|----------|-------|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 135 | (2E,4E)-hepta-2,4-dienal | 27.38 | 2.98 | 976 | R/Mo/F | 0.36 | 3.3 | 6.7 | 3.1 | 12.8 | 16.3 | 1.0 | 1.0 | 3.6 | 0.8 | 0.0 | 0.0 | 2.1 | 0.0 | 1.4 | 2.6 | 1.6 | 3.4 | 1.6 | 0.0 |
| | | 28.75 | 1.73 | 1459 | | | 4.1 | 7.3 | 3.6 | 11.9 | 16.3 | 1.3 | 0.9 | 3.6 | 0.7 | 0.0 | 0.0 | 1.9 | 0.0 | 0.8 | 2.6 | 1.3 | 2.8 | 1.6 | 0.0 |
| 137 | 6-methylhept-5-en-2-ol | 27.70 | 2.23 | 982 | Mo | 2 | 1.0 | 0.6 | 1.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 138 | octanal | 27.85 | 1.27 | 984 | Mo/R | 0.32 | 5.2 | 8.0 | 13.1 | 12.9 | 1.2 | 0.3 | 0.0 | 79.6 | 1.1 | 58.8 | 54.1 | 60.6 | 48.1 | 80.9 | 10.7 | 63.9 | 4.8 | 2.7 | 0.0 |
| | | 21.08 | 2.02 | 1280 | | | 6.5 | 8.8 | 15.3 | 12.0 | 1.2 | 0.3 | 0.0 | 78.6 | 0.9 | 78.6 | 51.8 | 55.9 | 61.1 | 50.1 | 10.4 | 51.8 | 3.9 | 2.8 | 0.0 |
| 140 | octan-2-ol | 27.94 | 1.63 | 986 | Mo | 0.1 | 0.7 | 0.0 | 3.5 | 1.3 | 0.0 | 0.0 | 0.0 | 0.0 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 25.50 | 1.60 | 1383 | | | 0.8 | 0.0 | 4.1 | 1.2 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 149 | hexyl acetate | 28.52 | 1.08 | 997 | Fr | 1.04 | 5.4 | 5.2 | 4.4 | 1.1 | 2.0 | 0.8 | 5.5 | 20.6 | 8.5 | 32.7 | 20.2 | 26.5 | 22.8 | 31.7 | 17.8 | 34.7 | 5.9 | 15.0 | 4.4 |
| | | 20.33 | 2.02 | 1263 | | | 6.8 | 5.7 | 5.1 | 1.1 | 2.0 | 1.1 | 4.9 | 20.3 | 7.3 | 43.7 | 19.4 | 24.4 | 28.9 | 19.6 | 17.3 | 28.2 | 4.9 | 15.5 | 5.0 |
| 169 | (2E)-oct-2-enal | 30.36 | 1.78 | 1034 | R | 0.004 | 2.4 | 2.7 | 5.2 | 22.4 | 1.0 | 0.0 | 1.9 | 2.6 | 1.0 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 0.0 |
| | | 27.25 | 1.89 | 1424 | | | 3.0 | 2.9 | 6.1 | 20.9 | 1.0 | 0.0 | 1.7 | 2.6 | 0.8 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 0.0 |
| 178 | o-guaiacol | 31.07 | 4.19 | 1049 | Mo | 0.02 | 5.3 | 0.0 | 28.7 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 185 | (3E)-octa-3,5-dien-2-one | 31.88 | 2.65 | 1065 | V/Mo | 0.0005 | 1.8 | 4.2 | 1.6 | 4.1 | 0.7 | 0.3 | 0.0 | 0.6 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.6 | 0.0 | 0.0 |
| | | 30.91 | 1.73 | 1561 | | | 2.3 | 4.6 | 1.9 | 3.8 | 0.7 | 0.4 | 0.0 | 0.6 | 0.0 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.3 | 0.0 | 0.0 |
| 191 | (3E)-6-methylhepta-3,5-dien-2-one | 32.38 | 2.77 | 1076 | V/Mo | 0.38 | 4.6 | 4.4 | 1.9 | 1.9 | 0.0 | 0.0 | 0.6 | 0.6 | 0.5 | 0.3 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 33.83 | 1.64 | 1582 | | | 5.8 | 4.8 | 2.2 | 1.8 | 0.0 | 0.0 | 0.5 | 0.6 | 0.5 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 196 | nonanal | 32.77 | 1.33 | 1084 | R | 0.15 | 110.8 | 64.4 | 137.2 | 56.4 | 11.7 | 3.4 | 31.5 | 34.6 | 15.9 | 13.9 | 8.2 | 10.4 | 2.8 | 5.8 | 7.0 | 9.4 | 33.2 | 7.4 | 2.3 |
| | | 25.75 | 2.19 | 1388 | | | 139.3 | 70.8 | 160.5 | 52.6 | 11.7 | 4.6 | 28.1 | 34.1 | 13.6 | 18.6 | 7.8 | 9.6 | 3.5 | 3.6 | 6.8 | 7.7 | 27.4 | 7.6 | 2.6 |
| 213 | ethyl cyclohexanecarboxylate | 34.44 | 1.33 | 1119 | Fr | 0.00016 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 217 | (2Z)-non-2-enal | 35.44 | 1.57 | 1141 | R | 0.0045 | 0.0 | 0.0 | 0.0 | 0.0 | 2.3 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 218 | (2E)-non-2-enal | 35.45 | 1.75 | 1141 | R | 0.9 | 1.3 | 1.0 | 2.2 | 4.3 | 0.5 | 0.0 | 0.8 | 0.6 | 0.4 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 0.4 | 0.0 |
| | | 31.75 | 1.98 | 1531 | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 230 | ethyl octanoate | 37.68 | 1.02 | 1189 | V | 10 | 6.7 | 2.8 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.5 | 0.6 | 0.0 | 2.9 | 0.0 | 0.8 | 0.0 | 0.0 | 0.4 | 0.0 | 0.0 | 0.0 |
| | | 27.50 | 2.36 | 1429 | | | 8.4 | 3.1 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.5 | 0.5 | 0.0 | 2.8 | 0.0 | 1.0 | 0.0 | 0.0 | 0.3 | 0.0 | 0.0 | 0.0 |
| 235 | Decanal | 37.94 | 1.27 | 1195 | R | 0.65 | 1.6 | 1.3 | 2.8 | 1.7 | 0.4 | 0.0 | 0.3 | 1.0 | 0.5 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.0 | 1.2 | 0.0 |
| | | 30.25 | 2.23 | 1494 | | | 2.0 | 1.5 | 3.3 | 1.6 | 0.4 | 0.0 | 0.3 | 1.0 | 0.4 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.6 | 1.2 | 0.0 |
| 245 | (2E)-dec-2-enal | 40.03 | 1.75 | 1243 | R | 0.01 | 7.6 | 3.6 | 7.8 | 7.6 | 1.0 | 0.4 | 2.2 | 6.6 | 2.8 | 1.3 | 1.4 | 3.2 | 1.7 | 0.8 | 0.6 | 3.1 | 2.4 | 2.3 | 0.9 |
| | | 36.08 | 2.02 | 1638 | | | 9.6 | 4.0 | 9.2 | 7.1 | 1.0 | 0.5 | 1.9 | 6.5 | 2.4 | 1.8 | 1.3 | 3.0 | 2.2 | 0.5 | 0.5 | 2.5 | 2.0 | 2.4 | 1.0 |
| 249 | (2E, 4Z)-deca-2,4-dienal | 41.29 | 2.32 | 1273 | R | 0.01 | 0.3 | 0.0 | 0.5 | 1.2 | 2.3 | 0.0 | 0.0 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.4 | 0.0 | 0.0 |
| | | 40.58 | 1.73 | 1756 | | | 0.4 | 0.0 | 0.6 | 1.2 | 2.3 | 0.0 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.3 | 0.0 | 0.0 |
| 251 | (2E,4E)-deca-2,4-dienal | 42.29 | 2.53 | 1296 | R | 0.18 | 0.3 | 0.0 | 0.4 | 1.7 | 3.6 | 0.0 | 0.4 | 0.6 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.3 | 0.0 | 0.0 |
| | | 42.25 | 1.64 | 1800 | | | 0.3 | 0.0 | 0.4 | 1.7 | 4.8 | 0.0 | 0.4 | 0.6 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.4 | 0.0 | 0.0 |
| 256 | ethyl decanoate | 46.19 | 1.08 | 1393 | V | 10 | 1.0 | 1.1 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.0 | 1.6 | 0.5 | 1.6 | 2.4 | 3.7 | 1.1 | 0.7 | 0.5 |
| | | 35.91 | 2.48 | 1634 | | | 1.2 | 1.2 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.0 | 1.4 | 0.6 | 1.0 | 2.3 | 3.0 | 0.9 | 0.7 | 0.6 |

Table 2. List of compounds adopted in the chemical blueprint of samples. ¹D and ²D retention times, experimental I_s^T , known correlation with sensory defect, odor threshold values (OT in mg/kg) and normalized peak volume (multiplied per 10⁵) are reported for both apolar × polar and polar × apolar (in *italic*) column set.





| ID# | | ¹ D (min) | ² D (s) | I_s^T | Sensory defect | OT (mg/kg) | Muddy | Vinegary | Mould | Fusty | Rancid | EVO-01 | EVO-02 | EVO-03 | EVO-04 | EVO-05 | EVO-06 | EVO-07 | EVO-08 | EVO-09 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|--------------------|-------------------------|-----------------------|-------------|-------------------|---------------|--------------|--------------|--------------|-------------|-------------|-------------|-------------|-------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|--------------|-------------|-------------|
| 1 | ethanol | 2.78 | 1.72 | 552 | V | 30 | 3.7 | 1.0 | 2.8 | 4.3 | 6.3 | 2.1 | 6.2 | 0.0 | 114.9 | 2.6 | 79.0 | 3.4 | 29.2 | 24.1 | 18.5 | 12.3 | 29.7 | 8.5 | 5.5 |
| | | <i>7.67</i> | <i>1.14</i> | <i>883</i> | | | <i>4.7</i> | <i>1.3</i> | <i>3.3</i> | <i>4.0</i> | <i>6.3</i> | <i>2.8</i> | <i>5.6</i> | <i>0.0</i> | <i>65.0</i> | <i>4.0</i> | <i>75.7</i> | <i>3.1</i> | <i>37.0</i> | <i>14.9</i> | <i>18.0</i> | <i>10.0</i> | <i>20.0</i> | <i>8.8</i> | <i>6.2</i> |
| 6 | butanal | 6.85 | 1.05 | 607 | F/M | 0.018 | 11.3 | 0.0 | 0.3 | 2.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | <i>7.00</i> | <i>2.48</i> | <i>857</i> | | | <i>14.2</i> | <i>0.0</i> | <i>0.4</i> | <i>3.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> |
| 7 | ethyl acetate | 7.52 | 0.96 | 613 | F/V | 0.94 | 18.9 | 17.2 | 5.3 | 4.7 | 0.0 | 0.0 | 0.0 | 9.4 | 16.9 | 0.0 | 85.7 | 40.2 | 12.4 | 10.7 | 8.5 | 2.8 | 0.0 | 3.2 | 0.0 |
| | | <i>6.75</i> | <i>1.35</i> | <i>850</i> | | | <i>23.8</i> | <i>18.9</i> | <i>6.2</i> | <i>4.3</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>9.2</i> | <i>9.6</i> | <i>0.0</i> | <i>82.1</i> | <i>37.1</i> | <i>15.8</i> | <i>6.6</i> | <i>8.3</i> | <i>2.3</i> | <i>0.0</i> | <i>3.2</i> | <i>0.0</i> |
| 9 | acetic acid | 8.11 | 1.66 | 622 | F/V/R | 0.5 | 12.5 | 23.7 | 18.4 | 32.9 | 12.6 | 4.6 | 12.6 | 21.0 | 194.2 | 5.6 | 3.6 | 37.0 | 13.5 | 5.0 | 2.9 | 6.7 | 0.7 | 1.1 | 8.1 |
| | | <i>28.50</i> | <i>0.97</i> | <i>1457</i> | | | <i>15.7</i> | <i>26.0</i> | <i>21.5</i> | <i>30.7</i> | <i>12.6</i> | <i>6.2</i> | <i>11.2</i> | <i>20.8</i> | <i>109.9</i> | <i>8.6</i> | <i>3.4</i> | <i>34.1</i> | <i>17.1</i> | <i>5.0</i> | <i>2.9</i> | <i>7.3</i> | <i>0.7</i> | <i>1.0</i> | <i>6.8</i> |
| 12 | butan-1-ol | 9.47 | 2.92 | 648 | V/M * | 0.15 | 1.6 | 4.0 | 0.0 | 17.1 | 0.0 | 0.3 | 1.3 | 1.8 | 0.6 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.2 | 0.0 |
| | | <i>14.42</i> | <i>1.26</i> | <i>1120</i> | | | <i>1.5</i> | <i>3.6</i> | <i>0.0</i> | <i>16.0</i> | <i>0.0</i> | <i>0.4</i> | <i>1.1</i> | <i>1.8</i> | <i>0.4</i> | <i>0.8</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>1.2</i> | <i>0.0</i> |
| 13 | 3-methylbutanal | 9.52 | 0.96 | 649 | F/V | 0.013 | 0.0 | 1.2 | 0.0 | 0.1 | 0.0 | 0.0 | 0.0 | 0.6 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | <i>8.75</i> | <i>2.61</i> | <i>927</i> | | | <i>0.0</i> | <i>1.5</i> | <i>0.0</i> | <i>0.1</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.6</i> | <i>1.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> |
| 14 | pent-1-en-3-one | 10.61 | 1.39 | 668 | M | 0.00073 | 0.7 | 1.6 | 0.6 | 2.4 | 1.4 | 29.5 | 2.0 | 6.1 | 0.4 | 4.6 | 0.9 | 24.3 | 1.5 | 0.0 | 3.4 | 2.8 | 0.3 | 0.5 | 8.9 |
| | | <i>10.17</i> | <i>1.47</i> | <i>983</i> | | | <i>0.9</i> | <i>1.8</i> | <i>0.8</i> | <i>2.2</i> | <i>1.4</i> | <i>39.6</i> | <i>1.8</i> | <i>6.0</i> | <i>0.2</i> | <i>7.1</i> | <i>0.9</i> | <i>22.4</i> | <i>1.9</i> | <i>0.0</i> | <i>3.3</i> | <i>2.3</i> | <i>0.2</i> | <i>0.6</i> | <i>10.1</i> |
| 15 | pent-1-en-3-ol | 10.71 | 2.77 | 670 | F/Mo | 0.4 | 1.6 | 6.8 | 1.2 | 11.3 | 11.6 | 18.0 | 8.4 | 7.5 | 9.6 | 15.1 | 6.0 | 8.3 | 9.4 | 11.7 | 11.8 | 10.8 | 7.4 | 5.3 | 4.5 |
| | | <i>15.08</i> | <i>1.30</i> | <i>1135</i> | | | <i>2.0</i> | <i>7.5</i> | <i>1.4</i> | <i>10.6</i> | <i>11.6</i> | <i>24.2</i> | <i>7.5</i> | <i>7.4</i> | <i>8.2</i> | <i>20.2</i> | <i>5.8</i> | <i>7.6</i> | <i>11.9</i> | <i>7.3</i> | <i>11.5</i> | <i>8.8</i> | <i>6.1</i> | <i>5.4</i> | <i>5.1</i> |
| 16 | pentan-2-one | 10.77 | 1.18 | 671 | V | 70 | 0.5 | 1.8 | 0.5 | 0.9 | 0.0 | 0.0 | 0.5 | 1.0 | 0.0 | 1.8 | 0.5 | 0.0 | 12.3 | 23.8 | 0.0 | 0.0 | 0.4 | 1.9 | 0.8 |
| | | <i>8.84</i> | <i>1.47</i> | <i>930</i> | | | <i>0.7</i> | <i>1.9</i> | <i>0.6</i> | <i>0.8</i> | <i>0.0</i> | <i>0.0</i> | <i>0.4</i> | <i>0.9</i> | <i>0.0</i> | <i>2.4</i> | <i>0.5</i> | <i>0.0</i> | <i>15.6</i> | <i>14.8</i> | <i>0.0</i> | <i>0.0</i> | <i>0.3</i> | <i>1.9</i> | <i>0.9</i> |
| 20 | 2-methylbutanal | 11.02 | 1.15 | 676 | V/M * | 0.01 | 3.8 | 16.0 | 7.5 | 9.3 | 0.7 | 16.7 | 20.5 | 18.2 | 5.1 | 0.0 | 16.3 | 0.0 | 0.0 | 0.0 | 21.1 | 18.9 | 8.4 | 19.2 | 0.0 |
| | | <i>4.34</i> | <i>1.09</i> | <i>750</i> | | | <i>4.7</i> | <i>17.6</i> | <i>8.7</i> | <i>8.7</i> | <i>0.7</i> | <i>22.4</i> | <i>18.3</i> | <i>18.0</i> | <i>4.4</i> | <i>0.0</i> | <i>15.6</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>20.5</i> | <i>15.3</i> | <i>6.9</i> | <i>19.8</i> | <i>0.0</i> |
| 23 | pentan-3-ol | 11.75 | 4.97 | 689 | V | 9 | 0.5 | 14.0 | 0.0 | 0.0 | 0.0 | 0.8 | 0.0 | 0.0 | 1.2 | 0.0 | 1.8 | 0.8 | 0.8 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | <i>13.00</i> | <i>1.35</i> | <i>1078</i> | | | <i>0.6</i> | <i>11.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>1.1</i> | <i>0.0</i> | <i>0.0</i> | <i>1.0</i> | <i>0.0</i> | <i>1.7</i> | <i>0.7</i> | <i>1.0</i> | <i>0.6</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> |
| 25 | ethyl propanoate | 12.10 | 0.93 | 696 | F | 0.1 | 4.8 | 0.4 | 0.4 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | <i>nr</i> | | | | | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> |
| 28 | methyl butanoate | 12.43 | 1.05 | 702 | F | 0.06 | 0.8 | 0.0 | 0.0 | 4.4 | 0.0 | 0.0 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | <i>nr</i> | | | | | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> |
| 31 | 3-methylbutan-1-ol | 13.46 | 2.62 | 720 | F/M/Mo | 0.1 | 2.4 | 8.7 | 1.0 | 2.9 | 0.0 | 1.0 | 5.3 | 18.1 | 23.1 | 3.1 | 15.5 | 11.5 | 4.9 | 8.0 | 3.3 | 4.1 | 2.8 | 3.6 | 1.2 |
| | | <i>18.35</i> | <i>0.94</i> | <i>1215</i> | | | <i>3.0</i> | <i>9.6</i> | <i>1.1</i> | <i>2.7</i> | <i>0.0</i> | <i>1.4</i> | <i>4.7</i> | <i>17.9</i> | <i>19.8</i> | <i>4.1</i> | <i>14.9</i> | <i>10.6</i> | <i>6.2</i> | <i>4.9</i> | <i>3.2</i> | <i>3.3</i> | <i>2.3</i> | <i>3.7</i> | <i>1.4</i> |
| 32 | (2E)-pent-2-enal | 13.53 | 1.93 | 721 | V | 0.3 | 5.9 | 12.9 | 6.9 | 21.4 | 10.5 | 4.4 | 3.4 | 5.7 | 1.2 | 4.5 | 1.6 | 0.0 | 1.2 | 0.5 | 1.4 | 1.1 | 2.9 | 1.3 | 1.0 |
| | | <i>14.08</i> | <i>1.52</i> | <i>1110</i> | | | <i>7.4</i> | <i>14.2</i> | <i>8.1</i> | <i>20.0</i> | <i>10.5</i> | <i>6.0</i> | <i>3.0</i> | <i>5.6</i> | <i>1.0</i> | <i>6.0</i> | <i>1.5</i> | <i>0.0</i> | <i>1.5</i> | <i>0.3</i> | <i>1.4</i> | <i>0.9</i> | <i>2.4</i> | <i>1.3</i> | <i>1.1</i> |
| 39 | pentan-1-ol | 14.88 | 2.92 | 746 | F/M/V | 3 | 1.0 | 2.1 | 0.9 | 4.5 | 1.9 | 1.0 | 3.5 | 2.1 | 2.0 | 3.7 | 4.0 | 5.8 | 1.0 | 6.0 | 1.3 | 1.4 | 3.4 | 7.4 | 0.0 |
| | | <i>19.00</i> | <i>1.35</i> | <i>1231</i> | | | <i>1.0</i> | <i>2.1</i> | <i>0.9</i> | <i>4.5</i> | <i>1.9</i> | <i>1.0</i> | <i>3.5</i> | <i>2.1</i> | <i>2.0</i> | <i>3.7</i> | <i>4.0</i> | <i>5.8</i> | <i>1.0</i> | <i>6.0</i> | <i>1.3</i> | <i>1.4</i> | <i>3.4</i> | <i>7.4</i> | <i>0.0</i> |
| 45 | (3Z)-hex-3-enal | 15.86 | 1.66 | 764 | Fr | 0.003 | 0.0 | 3.2 | 1.3 | 4.7 | 11.8 | 14.2 | 7.6 | 2.7 | 0.0 | 5.3 | 1.5 | 1.7 | 3.4 | 0.7 | 1.8 | 3.9 | 28.9 | 15.8 | 9.4 |
| | | <i>nr</i> | | | | | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> |
| 46 | Hexanal | 15.94 | 1.33 | 773 | F/Mo/V/R | 0.08 | 92.8 | 170.7 | 209.1 | 7.6 | 59.8 | 20.6 | 37.1 | 44.7 | 14.7 | 22.2 | 15.8 | 16.3 | 25.1 | 32.9 | 29.4 | 33.4 | 212.8 | 37.0 | 13.0 |
| | | <i>12.25</i> | <i>1.77</i> | <i>1054</i> | | | <i>116.7</i> | <i>187.7</i> | <i>244.7</i> | <i>7.1</i> | <i>59.8</i> | <i>27.7</i> | <i>33.1</i> | <i>44.2</i> | <i>12.6</i> | <i>29.6</i> | <i>15.2</i> | <i>15.0</i> | <i>31.9</i> | <i>20.4</i> | <i>28.6</i> | <i>27.1</i> | <i>175.8</i> | <i>38.0</i> | <i>14.8</i> |
| 48 | ethyl butanoate | 17.10 | 0.93 | 786 | F | 0.03 | 5.2 | 0.7 | 0.0 | 31.3 | 0.0 | 0.0 | 0.0 | 1.6 | 1.6 | 0.0 | 0.0 | 0.0 | 0.3 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | <i>10.59</i> | <i>1.77</i> | <i>1000</i> | | | <i>6.6</i> | <i>0.8</i> | <i>0.0</i> | <i>29.2</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>1.6</i> | <i>1.4</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.4</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> | <i>0.0</i> |


Table 3. continued.

| ID# | | ¹ D (min) | ² D (s) | ¹ T _s | Sensory defect | OT (mg/kg) | Muddy | Vinegary | Mould | Fusty | Rancid | EVO-01 | EVO-02 | EVO-03 | EVO-04 | EVO-05 | EVO-06 | EVO-07 | EVO-08 | EVO-09 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|--------------------------|-------------------------|-----------------------|-----------------------------|-------------------|---------------|-------|----------|-------|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 49 | oct-1-ene | 17.18 | 0.54 | 787 | M | 0.08 | 3.9 | 1.6 | 2.7 | 1.5 | 0.0 | 0.0 | 3.9 | 4.8 | 7.9 | 0.7 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 8.4 | 0.0 | 0.0 |
| | | 6.09 | 1.68 | 820 | | | 4.9 | 1.8 | 3.2 | 1.4 | 0.0 | 0.0 | 3.5 | 4.7 | 6.7 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 6.9 | 0.0 | 0.0 | |
| 52 | butyl acetate | 17.60 | 1.05 | 795 | F | 0.3 | 0.0 | 0.0 | 0.0 | 9.3 | 0.0 | 0.0 | 0.3 | 1.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 12.00 | 1.73 | 1046 | | | 0.0 | 0.0 | 0.0 | 8.7 | 0.0 | 0.0 | 0.2 | 1.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 53 | octane | 17.68 | 0.51 | 800 | F/V/R | 0.94 | 51.9 | 115.2 | 65.3 | 75.4 | 0.4 | 0.9 | 21.9 | 80.2 | 51.1 | 9.3 | 24.1 | 31.8 | 8.6 | 7.2 | 7.5 | 12.2 | 20.4 | 42.7 | 0.0 |
| | | 5.59 | 1.89 | 800 | | | 65.3 | 126.7 | 76.4 | 70.3 | 0.4 | 1.2 | 19.5 | 79.2 | 43.8 | 12.4 | 23.1 | 29.3 | 10.9 | 4.5 | 7.3 | 9.9 | 16.9 | 43.9 | 0.0 |
| 63 | (2E)-hex-2-enal | 19.11 | 1.84 | 822 | Mo/V/F/R | 0.42 | 1.0 | 1.7 | 1.9 | 8.2 | 565.4 | 582.0 | 396.0 | 37.0 | 10.4 | 209.1 | 218.4 | 209.0 | 545.4 | 444.5 | 565.9 | 425.7 | 406.8 | 391.3 | 782.0 |
| | | 18.00 | 1.64 | 1208 | | | 1.2 | 1.8 | 2.3 | 7.7 | 565.4 | 781.4 | 353.0 | 36.5 | 8.9 | 279.5 | 209.2 | 192.6 | 692.1 | 275.2 | 551.0 | 345.5 | 336.1 | 402.5 | 886.9 |
| 64 | (3E)-hex-3-en-1-ol | 19.31 | 3.69 | 826 | F/R/V | 6 | 0.7 | 0.4 | 0.6 | 0.8 | 1.6 | 0.6 | 0.7 | 1.4 | 2.6 | 2.6 | 3.7 | 2.6 | 1.5 | 2.1 | 1.3 | 3.0 | 1.9 | 5.9 | 2.2 |
| | | 24.08 | 1.30 | 1350 | | | 0.9 | 0.5 | 0.6 | 0.7 | 1.6 | 0.8 | 0.6 | 1.3 | 2.2 | 3.5 | 3.6 | 2.4 | 2.0 | 1.3 | 1.3 | 2.4 | 1.6 | 6.1 | 2.5 |
| 66 | (3Z)-hex-3-en-1-ol | 19.55 | 3.25 | 830 | V/F | 1.5 | 1.4 | 0.6 | 0.0 | 15.8 | 20.0 | 9.3 | 1.2 | 73.1 | 75.8 | 58.2 | 71.8 | 83.2 | 34.5 | 52.1 | 49.8 | 42.7 | 5.1 | 6.8 | 18.4 |
| | | 24.92 | 1.35 | 1369 | | | 1.8 | 0.6 | 0.0 | 14.7 | 20.0 | 12.5 | 1.1 | 72.1 | 64.9 | 77.8 | 68.8 | 76.7 | 43.8 | 32.2 | 48.5 | 34.6 | 4.2 | 7.0 | 20.8 |
| 68 | (2E)-hex-2-en-1-ol | 19.89 | 3.58 | 837 | V | 5 | 0.9 | 3.9 | 0.0 | 0.0 | 49.0 | 55.8 | 115.9 | 9.8 | 26.1 | 80.3 | 61.0 | 87.9 | 54.9 | 81.7 | 33.4 | 65.3 | 14.8 | 194.5 | 19.1 |
| | | 25.92 | 1.26 | 1392 | | | 1.1 | 4.3 | 0.0 | 0.0 | 49.0 | 75.0 | 103.3 | 9.7 | 22.4 | 107.3 | 58.5 | 81.0 | 69.6 | 50.6 | 32.6 | 53.0 | 12.2 | 200.0 | 21.6 |
| 74 | hexan-1-ol | 20.63 | 2.71 | 850 | Fr | 0.4 | 1.3 | 4.6 | 1.3 | 19.2 | 0.0 | 25.1 | 33.6 | 23.8 | 40.3 | 34.9 | 59.4 | 98.5 | 54.0 | 82.0 | 39.3 | 54.1 | 11.0 | 79.5 | 17.5 |
| | | 18.17 | 1.81 | 1212 | | | 1.7 | 5.0 | 1.5 | 17.9 | 0.0 | 33.7 | 30.0 | 23.5 | 34.6 | 46.6 | 56.9 | 90.8 | 68.6 | 50.8 | 38.2 | 43.9 | 9.0 | 81.8 | 19.8 |
| 80 | heptan-2-one | 21.85 | 1.27 | 872 | V | 0.3 | 3.8 | 10.0 | 7.6 | 7.5 | 0.5 | 0.3 | 0.7 | 2.5 | 0.9 | 3.4 | 2.7 | 1.6 | 0.6 | 0.6 | 0.0 | 1.4 | 2.1 | 0.0 | 0.0 |
| | | 16.42 | 1.89 | 1169 | | | 4.8 | 10.9 | 8.9 | 7.0 | 0.5 | 0.4 | 0.6 | 2.5 | 0.8 | 4.5 | 2.6 | 1.5 | 0.8 | 0.4 | 0.0 | 1.1 | 1.7 | 0.0 | 0.0 |
| 85 | heptanal | 22.35 | 1.24 | 881 | R | 0.5 | 10.8 | 6.9 | 16.7 | 7.3 | 1.6 | 0.3 | 3.5 | 5.6 | 1.6 | 1.1 | 2.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 8.0 | 2.0 | 0.0 |
| | | 16.50 | 1.89 | 1161 | | | 13.6 | 7.6 | 19.5 | 6.8 | 1.6 | 0.4 | 3.1 | 5.5 | 1.4 | 1.5 | 2.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 6.6 | 2.0 | 0.0 |
| 86 | ethyl pentanoate | 22.68 | 0.96 | 887 | M | 0.0015 | 1.3 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.8 | 1.3 | 0.0 |
| | | 23.08 | 2.27 | 1327 | | | 1.7 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 1.3 | 0.0 |
| 103 | (2Z)-hept-2-enal | 24.28 | 1.72 | 917 | R | 0.042 | 0.7 | 0.8 | 1.2 | 2.4 | 0.5 | 0.0 | 0.2 | 0.0 | 0.0 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.1 | 1.3 | 3.1 | 0.0 |
| | | 22.58 | 1.77 | 1315 | | | 0.8 | 0.9 | 1.4 | 2.2 | 0.5 | 0.0 | 0.2 | 0.0 | 0.0 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.7 | 1.0 | 3.1 |
| 110 | (2E)-hept-2-enal | 24.78 | 1.87 | 926 | Mo/R | 0.005 | 14.8 | 14.0 | 39.6 | 2.8 | 12.0 | 2.3 | 9.9 | 12.3 | 4.6 | 4.8 | 3.6 | 3.8 | 0.9 | 3.5 | 2.4 | 0.0 | 21.1 | 0.0 | 0.4 |
| | | 22.67 | 1.81 | 1317 | | | 18.6 | 15.5 | 46.3 | 2.6 | 12.0 | 3.1 | 8.8 | 12.2 | 4.0 | 6.4 | 3.5 | 3.5 | 1.2 | 2.2 | 2.3 | 0.0 | 17.4 | 0.0 | 0.5 |
| 122 | (2Z,4E)-hepta-2,4-dienal | 26.30 | 2.98 | 955 | R/Mo/F | 0.36 | 0.7 | 1.7 | 0.6 | 0.7 | 1.2 | 0.0 | 2.7 | 11.6 | 0.0 | 0.0 | 4.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.7 | 0.0 |
| | | 28.58 | 1.60 | 1455 | | | 0.9 | 1.9 | 0.8 | 0.7 | 1.5 | 0.0 | 2.4 | 11.5 | 0.0 | 0.0 | 4.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.8 |
| 123 | oct-1-en-3-one | 26.44 | 1.42 | 958 | Mo | 0.01 | 13.5 | 11.0 | 37.8 | 27.1 | 9.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 14.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 21.58 | 1.89 | 1292 | | | 16.9 | 12.1 | 44.3 | 25.3 | 9.1 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 8.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 125 | (2E,4Z)-hepta-2,4-dienal | 26.63 | 2.77 | 961 | R/Mo/F | 10 | 6.6 | 15.7 | 8.3 | 24.4 | 25.7 | 4.2 | 0.3 | 0.6 | 2.5 | 4.0 | 2.6 | 4.3 | 1.0 | 4.3 | 4.3 | 2.9 | 5.5 | 0.0 | 0.0 |
| | | 28.66 | 3.24 | 1457 | | | 8.3 | 17.2 | 9.7 | 22.8 | 25.7 | 5.7 | 0.3 | 0.6 | 2.2 | 5.3 | 2.5 | 4.0 | 1.3 | 2.6 | 4.2 | 2.4 | 4.6 | 0.0 | 0.0 |
| 126 | 6-methylhept-5-en-2-one | 26.69 | 1.63 | 962 | Mo/F/R | 1 | 148.2 | 88.6 | 88.6 | 70.7 | 2.1 | 0.8 | 4.7 | 11.0 | 4.4 | 4.0 | 2.9 | 3.0 | 0.0 | 0.0 | 0.0 | 1.9 | 4.5 | 2.0 | 0.0 |
| | | 23.17 | 1.85 | 1329 | | | 186.3 | 97.5 | 103.7 | 66.0 | 2.1 | 1.1 | 4.2 | 10.8 | 3.7 | 5.3 | 2.8 | 2.8 | 0.0 | 0.0 | 0.0 | 1.5 | 3.7 | 2.1 | 0.0 |
| 131 | oct-1-en-3-ol | 27.04 | 2.26 | 969 | Mo | 0.05 | 4.0 | 3.1 | 11.2 | 26.2 | 2.6 | 0.5 | 2.1 | 2.0 | 1.2 | 1.3 | 1.6 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 1.2 | 2.0 | 0.0 |
| | | 27.83 | 1.47 | 1437 | | | 0.6 | 0.3 | 1.0 | 1.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 134 | octan-2-one | 27.35 | 1.27 | 975 | V | 0.51 | 4.4 | 4.9 | 7.3 | 5.9 | 0.5 | 0.4 | 0.6 | 1.0 | 1.6 | 3.6 | 2.6 | 0.9 | 0.0 | 0.9 | 1.7 | 0.8 | 1.3 | 0.5 | 0.0 |
| | | 20.83 | 2.06 | 1274 | | | 0.0 | 0.0 | 0.1 | 0.2 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |


Table 3. continued.

| ID# | | ¹ D (min) | ² D (s) | I ^T _s | Sensory defect | OT (mg/kg) | Muddy | Vinegary | Mould | Fusty | Rancid | EVO-01 | EVO-02 | EVO-03 | EVO-04 | EVO-05 | EVO-06 | EVO-07 | EVO-08 | EVO-09 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|-----------------------------------|-------------------------|-----------------------|-----------------------------|-------------------|---------------|-------|----------|-------|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| 135 | (2E,4E)-hepta-2,4-dienal | 27.38 | 2.98 | 976 | R/Mo/F | 0.36 | 3.3 | 6.7 | 3.1 | 12.8 | 16.3 | 1.0 | 1.0 | 3.6 | 0.8 | 0.0 | 0.0 | 2.1 | 0.0 | 1.4 | 2.6 | 1.6 | 3.4 | 1.6 | 0.0 |
| | | 28.75 | 1.73 | 1459 | | | 4.1 | 7.3 | 3.6 | 11.9 | 16.3 | 1.3 | 0.9 | 3.6 | 0.7 | 0.0 | 0.0 | 1.9 | 0.0 | 0.8 | 2.6 | 1.3 | 2.8 | 1.6 | 0.0 |
| 137 | 6-methylhept-5-en-2-ol | 27.70 | 2.23 | 982 | Mo | 2 | 1.0 | 0.6 | 1.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 138 | octanal | 27.85 | 1.27 | 984 | Mo/R | 0.32 | 5.2 | 8.0 | 13.1 | 12.9 | 1.2 | 0.3 | 0.0 | 79.6 | 1.1 | 58.8 | 54.1 | 60.6 | 48.1 | 80.9 | 10.7 | 63.9 | 4.8 | 2.7 | 0.0 |
| | | 21.08 | 2.02 | 1280 | | | 6.5 | 8.8 | 15.3 | 12.0 | 1.2 | 0.3 | 0.0 | 78.6 | 0.9 | 78.6 | 51.8 | 55.9 | 61.1 | 50.1 | 10.4 | 51.8 | 3.9 | 2.8 | 0.0 |
| 140 | octan-2-ol | 27.94 | 1.63 | 986 | Mo | 0.1 | 0.7 | 0.0 | 3.5 | 1.3 | 0.0 | 0.0 | 0.0 | 0.0 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 25.50 | 1.60 | 1383 | | | 0.8 | 0.0 | 4.1 | 1.2 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 149 | hexyl acetate | 28.52 | 1.08 | 997 | Fr | 1.04 | 5.4 | 5.2 | 4.4 | 1.1 | 2.0 | 0.8 | 5.5 | 20.6 | 8.5 | 32.7 | 20.2 | 26.5 | 22.8 | 31.7 | 17.8 | 34.7 | 5.9 | 15.0 | 4.4 |
| | | 20.33 | 2.02 | 1263 | | | 6.8 | 5.7 | 5.1 | 1.1 | 2.0 | 1.1 | 4.9 | 20.3 | 7.3 | 43.7 | 19.4 | 24.4 | 28.9 | 19.6 | 17.3 | 28.2 | 4.9 | 15.5 | 5.0 |
| 169 | (2E)-oct-2-enal | 30.36 | 1.78 | 1034 | R | 0.004 | 2.4 | 2.7 | 5.2 | 22.4 | 1.0 | 0.0 | 1.9 | 2.6 | 1.0 | 0.8 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 0.0 |
| | | 27.25 | 1.89 | 1424 | | | 3.0 | 2.9 | 6.1 | 20.9 | 1.0 | 0.0 | 1.7 | 2.6 | 0.8 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.7 | 0.0 |
| 178 | o-guaiacol | 31.07 | 4.19 | 1049 | Mo | 0.02 | 5.3 | 0.0 | 28.7 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 185 | (3E)-octa-3,5-dien-2-one | 31.88 | 2.65 | 1065 | V/Mo | 0.0005 | 1.8 | 4.2 | 1.6 | 4.1 | 0.7 | 0.3 | 0.0 | 0.6 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.6 | 0.0 | 0.0 |
| | | 30.91 | 1.73 | 1561 | | | 2.3 | 4.6 | 1.9 | 3.8 | 0.7 | 0.4 | 0.0 | 0.6 | 0.0 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.3 | 0.0 | 0.0 |
| 191 | (3E)-6-methylhepta-3,5-dien-2-one | 32.38 | 2.77 | 1076 | V/Mo | 0.38 | 4.6 | 4.4 | 1.9 | 1.9 | 0.0 | 0.0 | 0.6 | 0.6 | 0.5 | 0.3 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | 33.83 | 1.64 | 1582 | | | 5.8 | 4.8 | 2.2 | 1.8 | 0.0 | 0.0 | 0.5 | 0.6 | 0.5 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 196 | nonanal | 32.77 | 1.33 | 1084 | R | 0.15 | 110.8 | 64.4 | 137.2 | 56.4 | 11.7 | 3.4 | 31.5 | 34.6 | 15.9 | 13.9 | 8.2 | 10.4 | 2.8 | 5.8 | 7.0 | 9.4 | 33.2 | 7.4 | 2.3 |
| | | 25.75 | 2.19 | 1388 | | | 139.3 | 70.8 | 160.5 | 52.6 | 11.7 | 4.6 | 28.1 | 34.1 | 13.6 | 18.6 | 7.8 | 9.6 | 3.5 | 3.6 | 6.8 | 7.7 | 27.4 | 7.6 | 2.6 |
| 213 | ethyl cyclohexanecarboxylate | 34.44 | 1.33 | 1119 | Fr | 0.00016 | 1.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 217 | (2Z)-non-2-enal | 35.44 | 1.57 | 1141 | R | 0.0045 | 0.0 | 0.0 | 0.0 | 0.0 | 2.3 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| | | nr | | | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 218 | (2E)-non-2-enal | 35.45 | 1.75 | 1141 | R | 0.9 | 1.3 | 1.0 | 2.2 | 4.3 | 0.5 | 0.0 | 0.8 | 0.6 | 0.4 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.9 | 0.4 | 0.0 |
| | | 31.75 | 1.98 | 1531 | | | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 |
| 230 | ethyl octanoate | 37.68 | 1.02 | 1189 | V | 10 | 6.7 | 2.8 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.5 | 0.6 | 0.0 | 2.9 | 0.0 | 0.8 | 0.0 | 0.0 | 0.4 | 0.0 | 0.0 | 0.0 |
| | | 27.50 | 2.36 | 1429 | | | 8.4 | 3.1 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.5 | 0.5 | 0.0 | 2.8 | 0.0 | 1.0 | 0.0 | 0.0 | 0.3 | 0.0 | 0.0 | 0.0 |
| 235 | Decanal | 37.94 | 1.27 | 1195 | R | 0.65 | 1.6 | 1.3 | 2.8 | 1.7 | 0.4 | 0.0 | 0.3 | 1.0 | 0.5 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.0 | 1.2 | 0.0 |
| | | 30.25 | 2.23 | 1494 | | | 2.0 | 1.5 | 3.3 | 1.6 | 0.4 | 0.0 | 0.3 | 1.0 | 0.4 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.6 | 1.2 | 0.0 |
| 245 | (2E)-dec-2-enal | 40.03 | 1.75 | 1243 | R | 0.01 | 7.6 | 3.6 | 7.8 | 7.6 | 1.0 | 0.4 | 2.2 | 6.6 | 2.8 | 1.3 | 1.4 | 3.2 | 1.7 | 0.8 | 0.6 | 3.1 | 2.4 | 2.3 | 0.9 |
| | | 36.08 | 2.02 | 1638 | | | 9.6 | 4.0 | 9.2 | 7.1 | 1.0 | 0.5 | 1.9 | 6.5 | 2.4 | 1.8 | 1.3 | 3.0 | 2.2 | 0.5 | 0.5 | 2.5 | 2.0 | 2.4 | 1.0 |
| 249 | (2E, 4Z)-deca-2,4-dienal | 41.29 | 2.32 | 1273 | R | 0.01 | 0.3 | 0.0 | 0.5 | 1.2 | 2.3 | 0.0 | 0.0 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.4 | 0.0 | 0.0 |
| | | 40.58 | 1.73 | 1756 | | | 0.4 | 0.0 | 0.6 | 1.2 | 2.3 | 0.0 | 0.0 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.3 | 0.0 | 0.0 |
| 251 | (2E,4E)-deca-2,4-dienal | 42.29 | 2.53 | 1296 | R | 0.18 | 0.3 | 0.0 | 0.4 | 1.7 | 3.6 | 0.0 | 0.4 | 0.6 | 0.4 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.3 | 0.0 | 0.0 |
| | | 42.25 | 1.64 | 1800 | | | 0.3 | 0.0 | 0.4 | 1.7 | 4.8 | 0.0 | 0.4 | 0.6 | 0.6 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 2.4 | 0.0 | 0.0 |
| 256 | ethyl decanoate | 46.19 | 1.08 | 1393 | V | 10 | 1.0 | 1.1 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.0 | 1.6 | 0.5 | 1.6 | 2.4 | 3.7 | 1.1 | 0.7 | 0.5 |
| | | 35.91 | 2.48 | 1634 | | | 1.2 | 1.2 | 0.5 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 0.0 | 1.0 | 1.4 | 0.6 | 1.0 | 2.3 | 3.0 | 0.9 | 0.7 | 0.6 |





Apolar × Polar



Polar × Apolar

395 peak features

Untargeted analysis

452 peak features

Targeted analysis

Step 1

261 targets peaks

Identification LRI/EI-MS fragmentation

218 targets peaks

Pattern recognition

Chemica blueprint definition

Data reduction

Successive pair-wise comparisons to identify discriminant features



Targeted analysis

Step 2

Data reduction

Localization of most relevant analytes

OT-OAV







Table 2. List of identified compounds in the apolar \times polar set, along with ^1D and ^2D retention times, experimental and library reported I^T_s , known correlation with sensory defects, odor descriptors and MS similarity matching (%). Mo: *mould*; M: *muddy*; F: *fusty*; V: *vinegary*; R: *rancid*; Fr: *fruity*. Defects marked with * refer to correlation highlighted using the proposed iterative process.

| ID# | Compound | ¹ D (min) | ² D (s) | I ^T _s | I ^T _{lib} | Sensory defect | sensory descriptor | Vinegary | Fusty | Mould | Muddy | Rancids | EVO-1 | EVO-2 | EVO-3 | EVO-4 | EVO-5 | EVO-6 | EVO-7 | EVO-8 | EVO-9 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|---------------------------|-------------------------|-----------------------|-----------------------------|-------------------------------|-------------------|--|----------|-------|-------|-------|---------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------|--------|--------|--------|--------|
| 1 | ethanol | 2.78 | 1.72 | 552 | 482 | V | Alcohol | 95 | 92 | 94 | 96 | 93 | 92 | 94 | 95 | 94 | 94 | 96 | 96 | 97 | 97 | 98 | 98 | 96 | 98 | 97 |
| 2 | tri-methylamine | 4.60 | 1.15 | 560 | 577 | Mo | fish | 93 | 91 | 91 | 91 | | | | 88 | | | | | | | | | | | |
| 3 | ethyl ether | 4.85 | 0.78 | 564 | 529 | | - | | | | 94 | | | | | | | | | | | | | | | |
| 4 | Propan-1-ol | 5.96 | 2.66 | 584 | 559 | Mo | alcohol, pungent | 92 | 90 | | | 91 | | 90 | 88 | 96 | 96 | | | | | | | 95 | | |
| 5 | 2-methylprop-2-enal | 6.11 | 1.30 | 587 | 574 | F | cooked, caramel | 90 | 95 | | | | | 90 | 90 | | | | | | | | | | | |
| 6 | butanal | 6.85 | 1.05 | 607 | 607 | F/M | pungent, green | | 93 | 96 | 94 | | | | | | | | | | | | | | | |
| 7 | ethyl acetate | 7.52 | 0.96 | 613 | 615 | F/V/M | pineapple | 91 | 90 | 91 | 91 | | | | 91 | 92 | | 90 | 88 | 96 | 95 | 93 | 94 | | 91 | |
| 8 | hexane | 7.59 | 0.45 | 614 | 618 | | alkane | | | | | | | | 88 | 86 | | | | | | | | | | |
| 9 | acetic acid | 8.11 | 1.66 | 622 | 622 | F/V/R | sour, vinegary | 96 | 92 | 91 | 92 | 90 | 91 | 90 | 92 | 90 | 91 | 90 | 96 | 96 | | | | | | |
| 10 | 2-methyl-3-Buten-2-ol | 8.28 | 1.78 | 626 | 600 | | herb | 95 | 95 | 94 | 93 | 92 | | 91 | 93 | | 90 | | | | | | | 94 | | |
| 11 | 2-Butenal | 8.61 | 1.87 | 632 | 615 | | green, fruit | 96 | 97 | 96 | 97 | 96 | 92 | 93 | 97 | | 91 | | | | | | | | | |
| 12 | butan-1-ol | 9.47 | 2.92 | 648 | 662 | V/M * | Winey | 92 | 95 | | 91 | | 90 | 88 | 91 | 89 | 88 | | | | | | | 93 | | |
| 13 | 3-methylbutanal | 9.52 | 0.96 | 649 | 643 | F/V | Malty | 91 | 92 | | | | | | 94 | 92 | | | | | | | | | | |
| 14 | pent-1-en-3-one | 10.61 | 1.39 | 668 | 677 | M | Mustard | 95 | 95 | 95 | 95 | 96 | 94 | 95 | 95 | 90 | 95 | 87 | 90 | 98 | | 98 | 98 | 98 | 98 | |
| 15 | pent-1-en-3-ol | 10.71 | 2.77 | 670 | 671 | F/Mo | wet earth | 96 | 95 | 95 | 91 | 95 | 92 | 96 | 95 | 96 | 95 | 96 | 98 | 98 | 98 | 98 | 98 | 98 | 98 | |
| 16 | pentan-2-one | 10.77 | 1.18 | 671 | 671 | V * | ether | 90 | 87 | 89 | 96 | | | 88 | 87 | | 89 | 93 | | 93 | 90 | | 90 | 95 | 90 | |
| 17 | 2-methylhexane | 11.01 | 0.42 | 676 | 680 | | - | | | | 95 | | | | | | 91 | | | | | | | | | |
| 18 | methyl 2-methylpropanoate | 11.02 | 0.87 | 676 | 680 | | flower | | | 92 | | | | | | | | | | | | | | | | |
| 19 | 2-Ethylfuran | 11.02 | 0.99 | 676 | 689 | | sweet, ethereal | 93 | 94 | 90 | 93 | | | | 90 | | 91 | | | | | | | | | |
| 20 | 2-methylbutanal | 11.02 | 1.15 | 676 | 662 | V/M * | Malty | 89 | 95 | 90 | 95 | 90 | 91 | 90 | 94 | 92 | 90 | 91 | | | | 92 | 90 | 92 | 90 | |
| 21 | 3,4-dihydro-2H-pyran | 11.03 | 1.60 | 676 | 691 | F/V * | - | 90 | | 90 | 91 | | 90 | 89 | 91 | | 88 | | | | | | | | | |
| 22 | 3-methylhexane | 11.34 | 0.45 | 682 | 685 | | - | 90 | | | 93 | | | | | | | | | | | | | | | |
| 23 | pentan-3-ol | 11.70 | 1.87 | 688 | 681 | V | fruit | 96 | | | 96 | | 95 | | 96 | 96 | | 94 | 94 | 95 | 93 | | | | | |
| 24 | 2-butanol-3-one | 11.75 | 4.97 | 689 | 716 | V * | butter, cream | | | | | | | | | | | 95 | 93 | | | | | | | |
| 25 | ethyl propanoate | 12.10 | 0.93 | 696 | 708 | F | Fruit, strong | 95 | 94 | 96 | 96 | | | | | | | | | | | | | | | |
| 26 | propyl acetate | 12.27 | 0.96 | 699 | 686 | F * | sour | 82 | | | 94 | | | | | | | | | | | | | | | |
| 27 | heptane | 12.34 | 0.45 | 700 | 700 | | alkane | 95 | 94 | 90 | 96 | | | | 96 | | | | | | | | | 93 | | |
| 28 | methyl butanoate | 12.43 | 1.05 | 702 | 718 | F | ether, fruit, sweet | | 95 | | 95 | | | | 88 | | | | | | | | | | | |
| 29 | (2Z)-pent-2-enal | 13.11 | 1.75 | 714 | 727 | Mo | strawberry, fruit, tomato, green, pleasant | 93 | 92 | | 92 | 94 | 94 | 94 | 94 | | 94 | | | | | | | | | |
| 30 | (3Z)-pent-3-en-2-one | 13.12 | 1.93 | 714 | 733 | | fish, pungent | 94 | 93 | 95 | 91 | | | | 90 | 90 | | | | | | | 90 | | 90 | |
| 31 | 3-methylbutan-1-ol | 13.46 | 2.62 | 720 | 731 | F/M/Mo | whiskey, malt, burnt | 90 | 94 | 92 | 90 | | 87 | 96 | 96 | 91 | 95 | 93 | 91 | 90 | 94 | 95 | 96 | 91 | 96 | 90 |
| 32 | (2E)-pent-2-enal | 13.53 | 1.93 | 721 | 751 | V | green, apple, tomato, pungent | 95 | 95 | 95 | 96 | 96 | 95 | 95 | 95 | 96 | 95 | 94 | | 93 | 90 | 91 | 92 | 92 | 94 | 93 |
| 33 | 2-(chloromethyl)but-1-ene | 14.02 | 0.96 | 730 | 711 | | - | | | | | 86 | | | | | | | | | | | | | | |
| 34 | toluene | 14.27 | 1.18 | 735 | 715 | | paint | 89 | 97 | 90 | 96 | 93 | 95 | 97 | 95 | 95 | 97 | 94 | | | | | | 85 | | |
| 35 | 2-methylpentan-3-one | 14.35 | 0.96 | 736 | 690 | | mint | 80 | 91 | | 90 | | | | | | | | | | | | | | | |
| 36 | (2Z)-pent-2-en-1-ol | 14.57 | 3.89 | 740 | 767 | | butter, pungent | | | | 89 | 95 | 93 | 90 | 91 | 90 | | 90 | | | | 80 | | | 90 | |
| 37 | (2E)-pent-2-en-1-ol | 14.65 | 4.04 | 742 | 769 | | butter, pungent | 88 | 90 | | | 97 | 96 | 97 | 96 | | 96 | 94 | 90 | 91 | 95 | 90 | 90 | 91 | 91 | |
| 38 | ethyl 2-methylpropanoate | 14.85 | 0.81 | 745 | 752 | | rancid, butter, cheese | | | 95 | 92 | | | | | | | | | | | | | | | |
| 39 | pentan-1-ol | 14.88 | 2.92 | 746 | 752 | F/M/V | Fruity | 95 | 96 | 93 | 94 | 97 | 95 | 95 | 95 | 95 | 96 | 92 | 91 | 91 | 91 | 95 | 92 | 95 | 96 | |
| 40 | (2E)-2-methylbut-2-enal | 15.12 | 2.38 | 750 | 738 | Mo * | green, fruit, aromatic | 92 | 93 | 91 | 93 | 92 | 90 | 90 | 92 | 91 | 92 | | | | | | 90 | | | |
| 41 | methyl 3-methylbutanoate | 15.35 | 0.99 | 754 | 721 | | sweat, acid, rancid | | | 92 | | | | | | | | | | | | | | | | |
| 42 | 2-methylpropyl acetate | 15.43 | 0.93 | 756 | 768 | F | fruit, apple, banana | 91 | | | 93 | | | | | | | | | | | | | | | |
| 43 | hexan-2-one | 15.61 | 1.30 | 759 | 754 | | ether, grape | | | 94 | | | | | | 96 | 93 | | | | | | | | | |
| 44 | 2-methylpent-1-en-3-ol | 15.71 | 2.59 | 761 | 747 | | alcohol, pungent | 93 | 94 | | 94 | 95 | 88 | | 82 | 80 | | | | | | | | | | |
| 45 | (3Z)-hex-3-enal | 15.86 | 1.66 | 764 | 797 | Fr | green | 92 | 95 | 90 | 87 | 95 | 95 | 95 | 96 | | 96 | 92 | 90 | 93 | 90 | 91 | 92 | 93 | 93 | |
| 46 | Hexanal | 15.94 | 1.33 | 773 | 801 | F/Mo/V/R | green apple, grassy | 95 | 96 | 96 | 97 | 97 | 96 | 96 | 96 | 98 | 97 | 95 | 95 | 96 | 96 | 96 | 95 | 90 | 96 | 95 |
| 47 | 1-methoxybutan-2-ol | 15.99 | 4.64 | 766 | 756 | R | - | 91 | 91 | 90 | 91 | 92 | | | | | | | | | | | 92 | | | |
| 48 | ethyl butanoate | 17.10 | 0.93 | 786 | 803 | F | sweet, fruity | 95 | 85 | | 96 | | | | 96 | 97 | | | 90 | | | | | | | |
| 49 | oct-1-ene | 17.18 | 0.54 | 787 | 807 | M | - | 97 | 96 | 97 | 97 | | | 98 | 96 | 98 | 96 | | | | | | 89 | | | |
| 50 | propyl propanoate | 17.43 | 0.93 | 792 | 785 | Mo | pineapple | | | | 94 | | | | | | | | | | | | | | | |

Table 2. continued.

| ID# | Compound | ¹ D (min) | ² D (s) | I ^r _s | I ^r _s lib | Sensory defect | sensory descriptor | Vinegary | Fusty | Mould | Muddy | Rancids | EVO-1 | EVO-2 | EVO-3 | EVO-4 | EVO-5 | EVO-6 | EVO-7 | EVO-8 | EVO-9 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|-------------------------------------|-------------------------|-----------------------|-----------------------------|------------------------------------|-------------------|--------------------------|----------|-------|-------|-------|---------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------|--------|--------|--------|--------|
| 101 | ethyl (2E)-2-methylbut-2-enoate | 24.11 | 1.33 | 913 | 929 | | - | | | 95 | | | | | 94 | 96 | | | | | | | | | | |
| 102 | gamma-pentanolactone | 24.12 | 1.93 | 913 | 941 | | - | 92 | | 84 | 88 | | | 87 | | | | | | | | | | | | |
| 103 | (2Z)-hept-2-enal | 24.28 | 1.72 | 917 | 927 | R | oxidized, tallowy | 86 | 88 | 88 | 87 | 86 | | 88 | | | 87 | | | | | 92 | 88 | 94 | | |
| 104 | (3E)-nona-1,3-diene | 24.43 | 0.78 | 919 | 914 | R * | - | | | | 88 | 90 | | | | 94 | | | | | | | | | | |
| 105 | 3-(propan-2-yl)cyclohex-1-ene | 24.44 | 1.30 | 920 | 897 | | - | 84 | 81 | 81 | 81 | 90 | 83 | 83 | 84 | 83 | 35 | | | | | | | 85 | | |
| 106 | 2,6-dimethylocta-1,3,7-triene | 24.44 | 1.51 | 920 | 916 | | | | | | 83 | | | | | | | | | | | | | | | |
| 107 | Benzaldehyde | 24.66 | 4.49 | 924 | 952 | | almond, burnt sugar | 96 | 94 | 96 | 90 | 92 | 92 | 97 | 92 | | 96 | 93 | | 93 | 95 | 92 | 94 | 93 | 92 | 92 |
| 108 | 3,4-dihydro-2H-pyran-2-carbaldehyde | 24.68 | 0.93 | 924 | 940 | Mo | - | 88 | 90 | 90 | 91 | 90 | 92 | 90 | 91 | | 93 | 92 | | | | | | | 88 | |
| 109 | methyl 2,4-dimethylpent-4-enoate | 24.77 | 0.87 | 926 | 887 | | - | | | | 82 | | | | | | | | | | | | | | | |
| 110 | (2E)-hept-2-enal | 24.78 | 1.87 | 926 | 931 | Mo/R | - | 97 | 95 | 97 | 97 | 97 | 97 | 96 | 97 | 96 | 97 | 95 | 95 | 92 | 93 | 90 | | 90 | | 88 |
| 111 | 1-butylcyclopent-1-ene | 25.18 | 0.64 | 932 | 969 | R | - | | | | | | | | | 96 | | | | | | | | | | |
| 112 | 2,6-dimethylheptan-4-one | 24.82 | 4.16 | 926-8 | 945 | | - | | | | | | | | | | | | | | | | | 88 | | |
| 113 | 7-ethylcyclohepta-1,3,5-triene | 25.19 | 1.15 | 934 | 946 | | - | | | | | | | | | | 91 | | | | | | | | | |
| 114 | 3-ethylocta-1,5-diene | 25.59 | 0.60 | 942 | 949 | | - | 95 | 90 | | 96 | 96 | 96 | 94 | 94 | 94 | 84 | 95 | 90 | 90 | 90 | 94 | 90 | 90 | 90 | 90 |
| 115 | cumene | 25.77 | 1.30 | 945 | 928 | | - | 90 | | 88 | 90 | 88 | 93 | 90 | 89 | 90 | 91 | | | | | | 90 | | | |
| 116 | (2E)-hept-2-en-1-ol | 25.81 | 3.40 | 946 | 964 | | - | 91 | 96 | 93 | 92 | | | 94 | 92 | 96 | | | | | | 90 | | 94 | | |
| 117 | 3-ethylocta-1,5-diene isomer | 25.83 | 4.49 | 946 | 949 | | - | 96 | | 90 | 95 | 96 | 95 | 82 | 95 | 94 | 95 | 90 | | 90 | 90 | 90 | | 90 | 90 | 90 |
| 118 | Phenol | 25.99 | 4.67 | 949 | 901 | | Phenol | | | | 91 | | | | | | | | 95 | | | | | | | |
| 119 | 3-ethylocta-1,5-diene isomer | 26.09 | 0.60 | 951 | 949 | | - | | | | 92 | | | | | | | | | | | | | | | |
| 120 | heptan-1-ol | 26.13 | 2.53 | 952 | 960 | | herb | 90 | 91 | 92 | 90 | | | 97 | 96 | 97 | | | | | | 90 | 92 | 90 | 91 | |
| 121 | 3,3-dimethyloctane | 26.27 | 0.93 | 954 | 931 | | | | | | 84 | | | | | | | | | | | | | | | |
| 122 | (2Z,4E)-hepta-2,4-dienal | 26.30 | 2.98 | 955 | 974 | R/Mo/F* | Fatty, rancid | 84 | 85 | 84 | 90 | | | 86 | 88 | | | 88 | | | | | | 90 | | |
| 123 | oct-1-en-3-one | 26.44 | 1.42 | 958 | 972 | Mo | mushroom, mould | 95 | 95 | 95 | 95 | 92 | | | 94 | 93 | | | | | | | | | | |
| 124 | (methoxymethyl)benzene | 26.54 | 2.11 | 959 | 969 | V | - | 90 | 93 | | 94 | | | | 93 | 95 | 91 | | | | | | | | | |
| 125 | (2E,4Z)-hepta-2,4-dienal | 26.63 | 2.77 | 961 | 971 | R/Mo/F* | Fatty, rancid | 93 | 93 | 93 | 93 | 92 | 93 | 93 | 92 | 92 | 92 | 90 | 91 | 90 | 92 | 95 | 91 | 91 | | |
| 126 | 6-methylhept-5-en-2-one | 26.69 | 1.63 | 962 | 981 | Mo/F/R | pungent, green | 92 | 93 | 85 | 86 | 94 | 94 | 94 | 94 | 94 | 94 | 93 | 91 | | | | 92 | 95 | 91 | |
| 127 | octan-3-one | 26.94 | 1.18 | 967 | 952 | M | herb, butter, resin | 90 | 94 | 94 | 93 | | 92 | 95 | 93 | 95 | | | | | | | 91 | 90 | | |
| 128 | 1-ethyl-3-methylbenzene | 26.94 | 1.36 | 967 | 967 | | - | 90 | 91 | 88 | 93 | | | 90 | 93 | | | | | | | 90 | | | | |
| 129 | Sabinene | 27.01 | 0.66 | 968 | 972 | | pepper, turpentine, wood | 88 | 90 | 90 | 91 | | 90 | 95 | 88 | 90 | | | | | | | 85 | | | |
| 130 | 1,3,5-trimethylbenzene | 27.02 | 1.21 | 969 | 994 | | - | 88 | | | 92 | | 93 | 93 | | 93 | 94 | | | | | | | 83 | | |
| 131 | oct-1-en-3-ol | 27.04 | 2.26 | 969 | 969 | Mo | mould, earthy | 96 | 96 | 96 | 96 | 96 | 95 | 96 | 96 | 96 | | 92 | | 91 | | | 94 | 95 | | |
| 132 | hex-3-ene-2,5-diol | 27.07 | 3.92 | 969-5 | 950 | M | - | | | | | | | | | | | | | | | | 85 | | | |
| 133 | 2-pentylfuran | 27.35 | 1.02 | 975 | 984 | | green bean, butter | 92 | 91 | 90 | 92 | 94 | | | | | | | | | | | | | | |
| 134 | octan-2-one | 27.35 | 1.27 | 975 | 989 | V | mould, green | 96 | 96 | 95 | 95 | 91 | 92 | 93 | 95 | 96 | 96 | 95 | 95 | | 93 | 90 | 94 | 92 | 93 | |
| 135 | (2E,4E)-hepta-2,4-dienal | 27.38 | 2.98 | 976 | 971 | R/Mo/F* | Fatty, rancid | 95 | 96 | 95 | 96 | 96 | 95 | 95 | 95 | 95 | 95 | | 93 | | 90 | 91 | 90 | 91 | 90 | |
| 136 | 1-methoxy-2-methylbenzene | 27.45 | 1.99 | 977 | 983 | | - | | | | | 94 | | | | | 90 | | | | | | | | | |
| 137 | 6-methylhept-5-en-2-ol | 27.70 | 2.23 | 982 | 989 | Mo | herbaceous, pungent | 90 | | 91 | 93 | | | | | | | | | | | | | | | |
| 138 | octanal | 27.85 | 1.27 | 984 | 998 | Mo/R | fatty, sharp | 97 | 97 | 96 | 97 | 97 | 97 | 97 | 95 | 94 | 95 | 95 | 91 | 93 | 85 | 96 | 90 | 97 | 95 | |
| 139 | ethyl hexanoate | 27.93 | 0.99 | 986 | 997 | F | apple peel, fruit | 95 | 92 | 91 | 95 | | | | 94 | 93 | 92 | | | | | | | | | |
| 140 | octan-2-ol | 27.94 | 1.63 | 986 | 999 | Mo | moss, nut, mushroom | | 91 | 95 | 92 | | | | | 91 | | | | | | | | | | |
| 141 | phenylmethanol | 28.08 | 4.70 | 989 | 1026 | | - | | 86 | | | | | | | | | | | | | | | | | |
| 142 | 3-ethylocta-1,5-diene isomer | 28.09 | 0.51 | 989 | | | - | 95 | | | 92 | | 94 | 94 | 94 | | 92 | | | | | 90 | | | 90 | |
| 143 | (3Z)-hex-3-en-1-yl acetate | 28.19 | 1.30 | 991 | 984 | | green, banana | 94 | 95 | 94 | 91 | 96 | 94 | 96 | 92 | 95 | 94 | 93 | 91 | 95 | 95 | 95 | 91 | 96 | 95 | 96 |
| 144 | 3-ethylocta-1,5-diene isomer | 28.26 | 0.66 | 992 | | | - | 82 | | | | | | | | | | | | | | 90 | 90 | 90 | 90 | |
| 145 | 1,2,4-trimethylbenzene | 28.28 | 1.54 | 992 | 1020 | | - | | | | | | | | | | 93 | | | | | | | | | |
| 146 | dec-1-ene | 28.34 | 0.57 | 994 | 1005 | | - | | | | | | | | | 96 | | | | | | | | | | |
| 147 | 3-ethylocta-1,5-diene isomer | 28.34 | 0.63 | 994 | | | - | | | | 93 | 94 | 94 | | 94 | | | 90 | 90 | | | | 90 | 90 | 90 | |
| 148 | dimethyl-2,5-dihydrofuran-2,5-dione | 28.46 | 2.32 | 996 | 996 | | - | 93 | 91 | | 94 | | | | | 90 | 92 | | | | | | | | | |
| 149 | hexyl acetate | 28.52 | 1.08 | 997 | 984 | Fr | green, fruity, sweet | 95 | 95 | 95 | 96 | 96 | 95 | 96 | 95 | 94 | 94 | 97 | 96 | 97 | 94 | 97 | 95 | 96 | 97 | 96 |
| 150 | (2E)-hex-2-en-1-yl acetate | 28.52 | 1.18 | 997 | 993 | | - | 90 | | | 91 | 85 | 91 | 85 | 88 | 90 | 93 | | 91 | 94 | 93 | | 90 | 90 | 93 | |

Table 2. continued.

| ID# | Compound | ¹ D (min) | ² D (s) | I ^T _s | I ^T _{lib} | Sensory defect | sensory descriptor | Vinegary | Fusty | Mould | Muddy | Rancids | EVO-1 | EVO-2 | EVO-3 | EVO-4 | EVO-5 | EVO-6 | EVO-7 | EVO-8 | EVO-9 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|------------------------------------|-------------------------|-----------------------|-----------------------------|-------------------------------|-------------------|---------------------------|----------|-------|-------|-------|---------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------|--------|--------|--------|--------|
| 151 | alpha-phellandrene | 28.60 | 0.69 | 999 | 1007 | | dill | | | | 88 | | | 93 | 90 | 89 | 88 | | | | | | | | | |
| 152 | 4-hydroxyhex-2-enoic acid | 28.72 | 2.95 | 1001 | 968 | | spice | 90 | 87 | 88 | 87 | 90 | 88 | 91 | 90 | 89 | | | | | | | | | | |
| 153 | 5-ethyl-2,5-dihydrofuran-2-one | 28.79 | 2.20 | 1002 | 977 | | - | | | 88 | | 87 | | | | | | | | | | | | | | |
| 154 | decane | 28.84 | 0.51 | 1003 | 1000 | | alkane | | 90 | 92 | 96 | 93 | | 91 | 90 | 90 | | | | | | | | | | |
| 155 | 2-ethylbutanal | 28.97 | 3.25 | 1006 | 997 | | - | 88 | 81 | 85 | 83 | 85 | | 86 | 84 | | 83 | | | | | | | | | |
| 156 | alpha-terpinene | 29.01 | 0.72 | 1007 | 1018 | | lemon | | | | | | | 95 | | | | | | | | | | | | |
| 157 | (1E)-ethylidenecycloheptane | 29.03 | 1.99 | 1007 | 1023 | | - | 84 | 83 | 82 | 83 | 83 | | 84 | | | | | | | | | | | | |
| 158 | (3E)-oct-3-en-2-one | 29.45 | 1.72 | 1016 | 1036 | | nut, crushed bug | 95 | 95 | 95 | 95 | 95 | 92 | 96 | 95 | 93 | 91 | 95 | | | | | | 95 | 91 | |
| 159 | limonene | 29.52 | 0.87 | 1017 | 1025 | | citrus, mint | 92 | 86 | 85 | 96 | 94 | | 95 | 96 | 95 | 90 | 92 | 88 | 93 | | 92 | 88 | | | |
| 160 | gamma-hexalactone | 29.60 | 0.78 | 1019 | 1042 | | - | 95 | 87 | 88 | 90 | 89 | | 88 | | 89 | | | | | | | | | | |
| 161 | 5-ethenyl-5-methyloxolan-2-one | 29.64 | 3.49 | 1020 | 1035 | | - | | | 87 | | | | | | | | | | | | | | | | |
| 162 | (2Z)-oct-2-enal | 29.69 | 1.66 | 1021 | 1049 | | green leaf, walnut | 89 | 90 | 89 | | 90 | | 88 | | | | | | | | | | 90 | | |
| 163 | beta-phellandrene | 29.85 | 0.75 | 1024 | 1025 | | mint, turpentine | | | 89 | | | | 92 | | | | | | | | | | | | |
| 164 | ethyl (2E)-hex-2-enoate | 29.86 | 1.30 | 1024 | 1041 | | - | 92 | | 94 | | | | | | | | | | | | | | | | |
| 165 | methyl 2-oxohexanoate | 30.27 | 1.18 | 1033 | 1020 | | - | 87 | 85 | 86 | 88 | | | | | | | | | | | | | | | |
| 166 | acetophenone | 30.32 | 4.37 | 1035 | 1029 | | must, flower, almond | 90 | 85 | 92 | 92 | | | 90 | 91 | | | | | | | | | 85 | 90 | 90 |
| 167 | (Z)-beta-ocimene | 30.35 | 0.93 | 1034 | 1044 | Mo | citrus, herb, flower | | 87 | 94 | 93 | | 96 | 95 | 96 | 96 | | | | | | | | 93 | 95 | 90 |
| 168 | butan-2-ylbenzene | 30.35 | 1.18 | 1034 | 1028 | | - | | | | 85 | | | | | | | | | | | | | | | |
| 169 | (2E)-oct-2-enal | 30.36 | 1.78 | 1034 | 1049 | R | green, nut, fat | 97 | 97 | 97 | 97 | 97 | | 97 | 97 | 97 | | | | | | | | | 95 | |
| 170 | (E)-beta-ocimene | 30.43 | 0.90 | 1036 | 1044 | Mo | sweet, herb | | | | 94 | | | | | | 96 | | | | | | | | | |
| 171 | 2,2,7,7-tetramethyloctane | 30.59 | 0.43 | 1039 | 1045 | | - | | | | | | | | | 93 | | | | | | | | | | |
| 172 | p-cymene | 30.60 | 1.21 | 1039 | 1040 | | solvent, gasoline, citrus | | | | 82 | | | 95 | | | 93 | | | | | | | | | |
| 173 | camphenilone | 30.69 | 1.54 | 1041 | 1078 | | - | 84 | 84 | 84 | 84 | | 84 | 84 | 84 | 84 | | | | | | | | | | |
| 174 | gamma-terpinene | 30.76 | 0.81 | 1043 | 1018 | | gasoline, turpentine | | | | | | | 91 | | 91 | | | | | | | | | | |
| 175 | methyl 2-ethylhexanoate | 30.77 | 0.93 | 1043 | 1019 | | - | 91 | | | | | | | | | | | | | | | | | | |
| 176 | 1-(cyclohex-1-en-1-yl)ethan-1-one | 30.79 | 2.38 | 1043 | 1027 | | - | 84 | 81 | 81 | 82 | 82 | 83 | 82 | 81 | | 83 | | | | | | 84 | 83 | | |
| 177 | 1-chlorooctane | 31.02 | 0.90 | 1048 | 1042 | | - | 88 | 87 | 90 | 92 | | | 88 | 89 | | | | | | | | | | | |
| 178 | o-guaiacol | 31.07 | 4.19 | 1049 | 1090 | Mo | smoke, sweet, medicine | | | 95 | 94 | | | | | | | | | | | | | | | |
| 179 | 3-methylbutyl butanoate | 31.27 | 0.93 | 1053 | 1057 | F/Mo | - | | | | 94 | | | | 91 | 90 | | | | | | | | | | |
| 180 | octan-1-ol | 31.54 | 2.32 | 1058 | 1059 | | moss, nut, mushroom | 91 | 89 | 90 | 90 | | 92 | 98 | 96 | 96 | 97 | | | | | | | | 95 | |
| 181 | trans-dihydro-rose oxide | 31.63 | 2.89 | 1060 | 1068 | | - | 87 | 86 | 86 | 86 | | | 88 | | 86 | | | | | | | | | | |
| 182 | cyclopenta[b]pyridine | 31.80 | 3.07 | 1064 | 1040 | | - | | | 86 | 86 | | | | | 85 | | | | | | | | | | |
| 183 | methyl benzoate | 31.80 | 3.04 | 1064 | 1060 | | - | 90 | | | 85 | 95 | 93 | 91 | 91 | 88 | 93 | 88 | 90 | 97 | 88 | 94 | 94 | 90 | | |
| 184 | o-cymene | 31.86 | 1.30 | 1065 | 1042 | | | | | | | | | 90 | 88 | | 90 | | | | | | | | | |
| 185 | (3E)-octa-3,5-dien-2-one | 31.88 | 2.65 | 1065 | 1068 | V/Mo * | geranium-like | 93 | 81 | 81 | 93 | 87 | 85 | | 86 | | 85 | | | | | | | | 88 | |
| 186 | meta-cymenene | 32.02 | 1.45 | 1068 | 1082 | Me | - | | | | | | | 92 | | | | | | | | | | | | |
| 187 | 2-(oxolan-3-yl)oxolane | 32.04 | 2.29 | 1069 | 1079 | | - | | 87 | 86 | 86 | | | 87 | | | | | | | | | | | | |
| 188 | terpinolene | 32.10 | 0.84 | 1070 | 1086 | | - | | | | | | | 92 | | | | | | | | | | | | |
| 189 | undeca-1,10-diene | 32.18 | 0.64 | 1071 | 1095 | | - | | | | | | | | | 87 | | | | | | | | | | |
| 190 | 1-ethenyl-3-ethylbenzene | 32.19 | 1.66 | 1072 | 1096 | | - | 93 | 91 | 90 | | | 90 | 88 | 90 | 88 | | | | | | | | | | |
| 191 | (3E)-6-methylhepta-3,5-dien-2-one | 32.38 | 2.77 | 1076 | 1102 | V/Mo * | - | 93 | 93 | 93 | 94 | | | 91 | 88 | 90 | 91 | | | | | | | | | |
| 192 | nonan-2-one | 32.44 | 1.27 | 1077 | 1089 | | hot milk, soap, green | 96 | 95 | 96 | 96 | | | 94 | 96 | 91 | 92 | | | | | | | 90 | | |
| 193 | cis-sabinene hydrate | 32.44 | 1.42 | 1077 | 1069 | | balsamic | | | | | | | 93 | | | | | | | | | | | | |
| 194 | (4E)-non-4-enal | 32.44 | 1.54 | 1077 | 1112 | | - | 90 | 88 | 90 | 91 | | | | 90 | 88 | | | | | | | | | | |
| 195 | pinocampnone | 32.59 | 0.54 | 1080 | 1109 | | - | | | | | | | 80 | | | | | | | | | | | | |
| 196 | nonanal | 32.77 | 1.33 | 1084 | 1107 | R | fatty, waxy, pungent | 89 | 91 | 95 | 97 | 98 | 91 | 97 | 97 | 98 | 98 | 96 | 90 | 95 | 96 | 92 | 96 | 92 | 96 | 94 |
| 197 | p-ethyl-anisole | 32.95 | 2.08 | 1087 | 1082 | | - | | | 92 | 90 | | | | | 88 | | | | | | | | | | |
| 198 | Phenethyl alcohol | 32.97 | 3.13 | 1088 | 1113 | | - | 90 | 91 | 90 | 91 | 94 | 90 | 92 | 91 | 88 | 96 | 90 | 90 | | 91 | 90 | | | 90 | |
| 199 | (3Z)-4,8-dimethylnona-1,3,7-triene | 33.26 | 0.72 | 1094 | 1113 | | - | 88 | 90 | | | 91 | | 94 | 88 | | 90 | 91 | | | 90 | | | | | |
| 200 | hexyl propanoate | 33.27 | 1.02 | 1094 | 1106 | | - | 90 | | 94 | | | | | | | | | | | | | | | | |

Table 2. continued.


| ID# | Compound | ¹ D (min) | ² D (s) | I ^T _s | I ^T _s lib | Sensory defect | sensory descriptor | Vinegary | Fusty | Mould | Muddy | Rancids | EVO-1 | EVO-2 | EVO-3 | EVO-4 | EVO-5 | EVO-6 | EVO-7 | EVO-8 | EVO-9 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|---|-------------------------|-----------------------|-----------------------------|------------------------------------|-------------------|--------------------------|----------|-------|-------|-------|---------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------|--------|--------|--------|--------|
| 201 | undec-1-ene | 33.34 | 0.48 | 1095 | 1105 | | - | | | | | | | 97 | | 97 | | | | | | | | | | |
| 202 | 2-ethenyl-1,1-dimethyl-3-methylidenecyclohexane | 33.35 | 0.84 | 1095 | 1071 | | - | | | | | | | | | | 85 | | | | | | | | | |
| 203 | heptyl acetate | 33.60 | 1.08 | 1100 | 1114 | | - | | | 95 | | | | | | | 90 | | | | | | | | | |
| 204 | linalool | 33.70 | 1.84 | 1103 | 1101 | | flower, lavender | | | | 94 | | | | 96 | 90 | 91 | | | | | | | 95 | | |
| 205 | undecane | 33.84 | 0.54 | 1106 | 1100 | | alkane | | | | 96 | | | 90 | | | | | | | | | | | | |
| 206 | (2Z)-undec-2-ene | 33.93 | 0.57 | 1108 | 1123 | | - | | | | | | | | | 94 | | | | | | | | | | |
| 207 | allo-ocimene | 33.94 | 1.08 | 1108 | 1128 | | - | | | | 90 | | 91 | 90 | 92 | 90 | 90 | | | | | | | | | |
| 208 | 1,2-dimethoxybenzene | 33.98 | 3.80 | 1109 | 1147 | | - | | | 93 | 93 | | | | | | | | | | | | | | | |
| 209 | (3E)-4,8-dimethylnona-1,3,7-triene | 34.02 | 0.90 | 1109 | 1113 | | - | 93 | 95 | 93 | 95 | | 95 | 90 | 95 | | 91 | | 91 | 93 | | | 94 | 92 | 90 | 91 |
| 210 | p-Mentha-1,3,8-triene | 34.02 | 1.27 | 1110 | 1108 | | - | | | | | | 90 | 91 | 91 | 90 | | | | | | | | | | |
| 211 | methyl octanoate | 34.19 | 1.11 | 1113 | 1123 | | orange | 94 | 90 | 91 | 92 | | | | 90 | 88 | | | | | | | | | | |
| 212 | (2E)-undec-2-ene | 34.34 | 0.60 | 1117 | 1123 | | - | | | | | | | | | 95 | | | | | | | | | | |
| 213 | ethyl cyclohexanecarboxylate | 34.44 | 1.33 | 1119 | 1147 | Fr | aromatic, fruity | | | | 93 | | | | | | | | | | | | | | | |
| 214 | neo-allo-ocimene | 34.60 | 1.15 | 1122 | 1140 | | - | | | | 88 | | 90 | 91 | 91 | 91 | 91 | | | | | | | 91 | | |
| 215 | (3E)-non-3-en-2-one | 34.61 | 1.51 | 1122 | 1137 | | - | | | | | 95 | | | | | | | | | | | | | | |
| 216 | (4E)-undeca-1,4-diene | 34.68 | 0.75 | 1124 | 1113 | | - | | | | | | | | | 91 | | | | | | | | | | |
| 217 | (2Z)-non-2-enal | 35.44 | 1.57 | 1141 | 1144 | R | green, fatty | | | | | 96 | | | | | | | | | | | | | | |
| 218 | (2E)-non-2-enal | 35.45 | 1.75 | 1141 | 1163 | R | paper-like, fatty | 96 | 97 | 97 | 97 | 94 | | 97 | 95 | 92 | 94 | | | | | | 91 | 90 | | |
| 219 | 4-ethylphenol | 35.48 | 3.92 | 1141 | 1141 | | - | 93 | 90 | 93 | 95 | | | | 94 | | | | | | | | | | | |
| 220 | 2-ethylphenol | 35.65 | 3.90 | 1145 | 1148 | | - | 90 | 93 | 93 | 92 | | | | | 93 | | | | | | | | | | |
| 221 | ethyl benzoate | 35.71 | 2.77 | 1146 | 1160 | | camomile, flower, celery | 80 | | 88 | 88 | 90 | 93 | | 92 | 90 | 91 | | | | | | | | | |
| 222 | (2E,9E)-undeca-2,9-diene | 36.01 | 0.82 | 1153 | 1131 | | - | | | | | | | | | 86 | | | | | | | | | | |
| 223 | 2-methoxy-4-methylphenol | 36.29 | 2.41 | 1159 | 1188 | | - | | | | 91 | | | | | | | | | | | | | | | |
| 224 | methyl 2-hydroxybenzoate | 36.57 | 3.95 | 1165 | 1190 | | - | 92 | | | 92 | 96 | 94 | 93 | 91 | 94 | 94 | 90 | | 91 | | 92 | 91 | 93 | 90 | 92 |
| 225 | (2Z)-non-2-en-1-ol | 37.02 | 1.21 | 1175 | 1162 | | - | 85 | 82 | | | | | | | | | | | | | | | | | |
| 226 | (2E)-non-2-en-1-ol | 37.35 | 1.21 | 1182 | 1163 | | - | | 82 | | 81 | | | | | | | | | | | | | | | |
| 227 | decan-2-one | 37.44 | 1.27 | 1184 | 1196 | | - | 95 | | 90 | 90 | | | | | 88 | 90 | | | | | | | | | |
| 228 | hexyl butanoate | 37.52 | 0.87 | 1186 | 1195 | | apple peel | | | | | | | | 82 | | | | | | | | | | | |
| 229 | (Z)-linalool oxide | 37.52 | 1.05 | 1186 | 1164 | | flower | 85 | 81 | | | 83 | | | | | | | | | | | | 83 | | |
| 230 | ethyl octanoate | 37.68 | 1.02 | 1189 | 1196 | V* | fruit, fat | 96 | | 93 | 96 | | | | | 90 | 92 | | 88 | | 90 | | 90 | | | |
| 231 | alpha-terpineol | 37.71 | 2.41 | 1190 | 1186 | | oil, anise, mint | | | 94 | 93 | | | | | | | | | | | | | | | |
| 232 | (2E,4E)-nona-2,4-dienal | 37.71 | 2.59 | 1190 | 1210 | R | watermelon | | | 92 | 95 | 93 | 90 | | 88 | | | | | | | | | | | |
| 233 | dodec-1-ene | 37.93 | 0.60 | 1195 | 1191 | | - | 97 | 91 | 94 | 97 | | | 90 | | 97 | | | | | | | | | | |
| 234 | (E)-linalool oxide | 37.93 | 1.05 | 1195 | 1164 | | flower | 81 | 80 | 81 | 80 | 80 | 81 | 88 | | 80 | 85 | | | | | | | | | |
| 235 | Decanal | 37.94 | 1.27 | 1195 | 1201 | R | penetrating, sweet, waxy | 96 | 92 | 97 | 96 | 96 | | 92 | 94 | 92 | 92 | | | | | | 93 | 95 | | |
| 236 | octyl acetate | 38.19 | 1.11 | 1200 | 1211 | | fruit | | | 94 | 91 | | | | | | | | | | | | | | | |
| 237 | Dodecane | 38.34 | 0.57 | 1204 | 1200 | | alkane | 93 | 97 | | 91 | | 91 | 93 | | 97 | 90 | | | | | | | | | |
| 238 | (2Z)-dodec-2-ene | 39.18 | 0.60 | 1223 | 1194 | | - | 91 | | 93 | 92 | 92 | 92 | 92 | 92 | 93 | 91 | 93 | 92 | | | 90 | 91 | 90 | 90 | 90 |
| 239 | (8E)-dodeca-1,8-diene | 39.18 | 0.78 | 1223 | 1212 | | - | | | | | | | 86 | | 87 | | | | | | | | | | |
| 240 | alkene | 39.68 | 0.57 | 1235 | | | - | | 92 | 92 | 93 | | 91 | | 93 | 92 | 92 | | | 90 | | | | | | |
| 241 | 2-phenylethyl acetate | 39.80 | 3.16 | 1238 | 1254 | | - | 95 | 93 | | | | | | | | 90 | | | | | | | | | |
| 242 | dodeca-1,11-diene | 39.93 | 0.69 | 1241 | 1194 | | - | 91 | | | | | | | | 90 | | | | | | | | | | |
| 243 | Menth-1en-7-al | 39.97 | 3.01 | 1242 | 1196 | | - | 80 | | 80 | 80 | | | | | 81 | 80 | | | | | | | | | |
| 244 | alkene | 40.01 | 0.57 | 1243 | | | - | 91 | | 91 | | | | | | 92 | | | | | | | | | | |
| 245 | (2E)-dec-2-enal | 40.03 | 1.75 | 1243 | 1265 | R | painty, fishy, fatty | 97 | 97 | 96 | 97 | 97 | 97 | 97 | 97 | 97 | 96 | 95 | 94 | 95 | 93 | 90 | 92 | 96 | 92 | 94 |
| 246 | tridec-1-ene | 40.09 | 0.60 | 1245 | 1290 | | - | | | | | | | | | 94 | 92 | | | | | | | | | |
| 247 | alkene | 40.09 | 0.60 | 1245 | | | - | | | | 91 | | | | | 90 | | | | | | | | | | |
| 248 | 4-ethylguaiaacol | 40.27 | 1.21 | 1249 | 1275 | | spice, clove | 94 | 90 | 95 | 94 | | | | | | | | | | | | | | | |
| 249 | (2E, 4Z)-deca-2,4-dienal | 41.29 | 2.32 | 1273 | 1292 | R | deep-fried | | 94 | 93 | 94 | 95 | | | 90 | | | | | | | | 94 | | | |
| 250 | ethyl nonanoate | 42.10 | 1.05 | 1292 | 1297 | | - | 92 | | | 93 | | | | | | | | | | | | | | | |

Table 2. continued.

| ID# | Compound | ¹ D (min) | ² D (s) | I ^r _s | I ^r _s lib | Sensory defect | sensory descriptor | Vinegary | Fusty | Mould | Muddy | Rancids | EVO-1 | EVO-2 | EVO-3 | EVO-4 | EVO-5 | EVO-6 | EVO-7 | EVO-8 | EVO-9 | EVO-10 | EVO-11 | EVO-12 | EVO-13 | EVO-14 |
|-----|------------------------------|-------------------------|-----------------------|-----------------------------|------------------------------------|-------------------|----------------------|----------|-------|-------|-------|---------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--------|--------|--------|--------|--------|
| 251 | (2E,4E)-deca-2,4-dienal | 42.29 | 2.53 | 1296 | 1315 | R | deep-fried | | 96 | 94 | 95 | 96 | | 91 | 90 | 91 | | | | | | | | 92 | | |
| 252 | 4-ethyl-1,2-dimethoxybenzene | 42.30 | 2.98 | 1296 | 1271 | | - | | | 93 | 92 | | | | | | | | | | | | | | | |
| 253 | methyl 2-methoxybenzoate | 42.84 | 0.48 | 1310 | 1334 | | warm, flower, walnut | | | | | | 90 | | 93 | 91 | | | | | | | | | | |
| 254 | (2E)-2-undecen-1-ol | 43.91 | 4.70 | 1336 | 1365 | | - | 81 | 82 | 84 | 82 | | | | | | | | | | | | | | | |
| 255 | gamma-nonalactone | 44.16 | 4.58 | 1343 | 1358 | | coconut, peach | 94 | | 91 | 90 | | | | | | | | | | | | | | | |
| 256 | ethyl decanoate | 46.19 | 1.08 | 1393 | 1395 | V* | grape | 94 | | 90 | 94 | | | | | | | 91 | 90 | 91 | 92 | 92 | 93 | 91 | 91 | 90 |
| 257 | copaene | 46.35 | 0.87 | 1397 | 1374 | | wood, spice | 94 | 95 | 93 | 95 | 95 | 92 | 92 | 93 | 95 | 95 | 90 | | | | 91 | 92 | 91 | 90 | |
| 258 | tetradec-1-ene | 46.43 | 0.66 | 1399 | 1403 | | - | | 91 | 90 | 93 | | | 90 | | 90 | 91 | | | | | | | | | |
| 259 | tetradecane | 46.76 | 0.60 | 1408 | 1400 | | alkane | | 96 | | 94 | | | | 90 | | 90 | | | | | | | | | |
| 260 | geranyl acetone | 48.45 | 1.72 | 1451 | 1453 | | magnolia, green | 93 | 93 | 94 | 94 | | | | 90 | | 91 | | | | | | | | | |
| 261 | alpha-farnesene | 50.44 | 1.18 | 1502 | 1505 | | wood, sweet | 94 | 91 | 92 | 90 | 95 | 90 | 90 | 90 | 96 | 95 | 90 | 91 | 90 | | 90 | 91 | 90 | 92 | 93 |

Muddy
Vinegary
Mould
Fusty
Rancid
EVO-1
EVO-2
EVO-3
EVO-4
EVO-5
EVO-6
EVO-7
EVO-8
EVO-9
EVO-10
EVO-11
EVO-12
EVO-13
EVO-14

a



Muddy
Vinegary
Mould
Fusty
Rancid
EVO-1
EVO-2
EVO-3
EVO-4
EVO-5
EVO-6
EVO-7
EVO-8
EVO-9
EVO-10
EVO-11
EVO-12
EVO-13
EVO-14

b

