

AperTO - Archivio Istituzionale Open Access dell'Università di Torino

Wrinkles, folds and calcifications: reduced durability after transcatheter aortic valve-in-valve replacement

This is the author's manuscript

Original Citation:

Availability:

This version is available <http://hdl.handle.net/2318/1627575> since 2017-03-10T11:34:06Z

Published version:

DOI:<http://dx.doi.org/10.1016/j.jtcvs.2016.08.018>

Terms of use:

Open Access

Anyone can freely access the full text of works made available as "Open Access". Works made available under a Creative Commons license can be used according to the terms and conditions of said license. Use of all other works requires consent of the right holder (author or publisher) if not exempted from copyright protection by the applicable law.

(Article begins on next page)

This is the author's final version of the contribution published as:

Grubitzsch, H., Galloni, M., Falk, V. - Wrinkles, folds and calcifications: reduced durability after transcatheter aortic valve-in-valve replacement - *Journal of thoracic and cardiovascular surgery*, 153 (2), 2016, p. 266-268,
<http://dx.doi.org/10.1016/j.jtcvs.2016.08.018>

The publisher's version is available at:

<http://www.sciencedirect.com/science/article/pii/S0022522316310376>

When citing, please refer to the published version.

Link to this full text:

<http://hdl.handle.net/2318/1627575>

This full text was downloaded from iris-AperTO: <https://iris.unito.it/>

Wrinkles, folds and calcifications: reduced durability after transcatheter aortic valve-in-valve replacement

Herko Grubitzsch¹, MD, PhD, Karl Stangl², MD, PhD, Marco Galloni^{3,4}, MSc, Volkmar Falk¹, MD, PhD

¹Klinik für Kardiovaskuläre Chirurgie, Charité - Universitätsmedizin Berlin, Berlin, Germany,

²Medizinische Klinik mit Schwerpunkt Kardiologie und Angiologie, Charité – Universitätsmedizin, Berlin, Berlin, Germany

³Dipartimento di Scienze Veterinarie, Laboratorio di Biomateriali, Università degli Studi di Torino, Turin, Italy

⁴Life and Device S.r.l., Turin, Italy

Transcatheter valve-in-valve implantation is an acceptable alternative therapy for failed aortic or mitral bioprostheses in selected patients.¹⁻³ With improved survival, durability of transcatheter valves becomes increasingly important because early structural failure frequently requires complex surgery.⁴

Figure 1 shows a self-expandable transcatheter aortic valve (CoreValve, 23 mm; Medtronic, Minneapolis, Minn) implanted through transfemoral access in a stented bioprosthesis (Sorin Mitroflow, 21 mm; LivaNova PLC, London, United Kingdom) 3 years previously (**Figure 1, A**, and **Video 1**), leading to an effective aortic valve area of 1.5 cm² and a reduction of the mean transvalvular gradient from 50 to 31 mm Hg initially. This valve was explanted (**Figure 1, B**, CoreValve-Mitroflow-corpus) because of severe aortic valve stenosis.

Effective orifice area of the valve was 0.8 cm² (mean transvalvular gradient, 48mmHg), as determined by preoperative echocardiography. The following findings were revealed by structured examination of the explanted transcatheter valve.

Wrinkled pericardium at the outflow aspect (**Figure 1, B**) and inflow aspect (**Figure 1, C**) of the prosthesis, an ovalized, noncircular cross-sectional profile (**Figure 1, D**),

and deep folds at each leaflet (Figure 1, B, C, and D) were indicative of an incomplete prosthesis expansion. Radiographic examination demonstrated multiple calcifications at all leaflets (Figure 1, D). Apparently, deposition of calcium started within the leaflets (Figure 1, E); these deposits evolved into vegetating calcifications (Figure 1, B and C) and were present in all leaflets at the inflow side (Figure 1, C) and in 2 of 3 leaflets at the outflow side (Figure 1, B). Altered collagen bundles, a bubbly tissue texture, and pericardial delamination were associated with leaflet folds (Figure 1, F) and may represent early stages of bioprosthetic degeneration. Because of the risk of accelerated structural degeneration, incomplete prosthesis expansion as a result of oversizing has to be prevented when performing transcatheter valve-in-valve procedures, especially for failed stented bioprostheses in the aortic position. The true internal diameter should be used for selecting an appropriate transcatheter heart valve device.⁵

References

1. Bapat V, Attia R, Redwood S, Hancock J, Wilson K, Young C, et al. Use of transcatheter heart valves for a valve-in-valve implantation in patients with degenerated aortic bioprosthesis: technical considerations and results. *J Thorac Cardiovasc. Surg.* 2012;144:1372-9; discussion 1379-80.
2. Conradi L, Silaschi M, Seiffert M, Lubos E, Blankenberg S, Reichenspurner H, et al. Transcatheter valve-in-valve therapy using 6 different devices in 4 anatomic positions: clinical outcomes and technical considerations. *J Thorac Cardiovasc. Surg.* 2015;150:1557-65. 1567.e1-3; discussion 1565-7.
3. Ye J, Webb JG, Cheung A, Soon JL, Wood D, Thompson CR, et al. Transapical transcatheter aortic valve-in-valve implantation: clinical and hemodynamic outcomes beyond 2 years. *J Thorac Cardiovasc Surg.* 2013;145:1554-62.
4. Seeburger J, Weiss G, Borger MA, Mohr FW. Structural valve deterioration of a CoreValve prosthesis 9 months after implantation. *Eur Heart J.* 2013;34:1607.
5. Bapat V, Attia R, Thomas M. Effect of valve design on the stent internal diameter of a bioprosthetic valve: a concept of true internal diameter and its implications for the valve-in-valve procedure. *JACC Cardiovasc Interv.* 2014;7:115-27.




FIGURE 1. Early structural failure after transcatheter aortic valve-in-valve replacement. A, Section of patient's chest radiograph showing a 23-mm Core Valve transcatheter valve (Medtronic, Minneapolis, Minn) within a 21-mm Mitroflow conventional stented bioprosthesis (LivaNova PLC, London, United Kingdom). B, The explanted CoreValve-Mitroflow-corpus presenting pericardial folds and vegetating calcifications at the outflow aspects. C, Wrinkled pericardium and calcifications at the inflow aspect of the explanted CoreValve prosthesis. D, Radiographic examination demonstrates an ovalized, noncircular cross-sectional profile of the transcatheter valve within the stented bioprosthesis, as well as multiple calcifications at all leaflets of the CoreValve prosthesis.

E, Alizarin red staining of a leaflet fold with intrinsic calcium nodule (3-fold magnification). F, Hematoxylin-eosin staining of a leaflet fold showing homogenized and locally disrupted collagen bundles, a bubbly tissue texture, and pericardial delamination (20-fold magnification).

VIDEO 1. Fluoroscopy demonstrating the degenerated 21-mm Mitroflow prosthesis (LivaNova PLC, London, United Kingdom) and the valve-in-valve procedure with a 23-mm CoreValve prosthesis (Medtronic, Minneapolis, Minn). Video available at: [http://www.jtcvsonline.org/article/S0022-5223\(16\)31037-6/addons](http://www.jtcvsonline.org/article/S0022-5223(16)31037-6/addons).