

AperTO - Archivio Istituzionale Open Access dell'Università di Torino

Long-term culture mantains stemness in bovine mammary epithelial cells

This is the author's manuscript

Original Citation:

Availability:

This version is available <http://hdl.handle.net/2318/153666> since 2015-09-22T08:03:34Z

Published version:

DOI:10.1111/rda.12220

Terms of use:

Open Access

Anyone can freely access the full text of works made available as "Open Access". Works made available under a Creative Commons license can be used according to the terms and conditions of said license. Use of all other works requires consent of the right holder (author or publisher) if not exempted from copyright protection by the applicable law.

(Article begins on next page)



UNIVERSITÀ DEGLI STUDI DI TORINO

This is an author version of the contribution published on:

Questa è la versione dell'autore dell'opera:

Long-term culture mantains stemness in bovine mammary epithelial cells

D. Cravero, E. Martignani, S. Miretti, P. Accornero, M. Baratta

Reprod Dom Anim 48 (Suppl. 1), 72–124 (2013);

doi: 10.1111/rda.12220 ISSN 1439-0531

Long-term culture maintains stemness in bovine mammary epithelial cells

D Cravero, E Martignani, S Miretti, P Accornero, M Baratta

Department of Veterinary Science, University of Turin, Turin, Italy

The regenerative potential of mammary gland stem cells may be exploited to regulate milk production in dairy cattle. Our recent works demonstrated the intrinsic regenerative potential found in bovine mammary tissue, but a greater challenge is to purify and enrich stem cells that reside in the mammary gland. We describe here the phenotype of mammary epithelial cells emerging in two defined culture conditions and the assessment of their *in vivo* regenerative potential after long-term culture. We used Epicult B and SF7 media (both are used to culture mouse and human mammary epithelial cells) to expand bovine cells. We analyzed cell phenotype at different passages (P0, 1, 3 and 5; cells were passaged every 5 days). We also tested the regenerative potential of cells cultured in different media through Colony Forming Cell assay (CFC) and through xenotransplantation in NOD/SCID mice. Cytokeratin 18 (CK18) and cytokeratin 14 (K14) were observed in both media until P5. Significant differences in marker expression were found that were associated with the medium used. Vimentin (a stromal marker) was differently expressed in EpiB and SF7 at P5 (0% for EpiB and 2% for SF7, $p < 0.05$). K14 and K18, P63 and EpCam expression was also evaluated in regenerated polarized structures found in NOD/SCID mice. These results demonstrate that a long-term culture is influenced by the medium used but maintains a multi-potent cell subpopulation with intrinsic regenerative potential.

(This work was supported by PRIN 2010-2011 grant.)