

Potential of the Reversed-Inject Differential Flow Modulator for Comprehensive Two-dimensional Gas Chromatography in the Quantitative Profiling of Complex Natural Samples

Chiara Cordero¹, Carlo Bicchi¹, Luigi Cobelli², Gianluca Stani², Armando Miliazza³, Roger Firor⁴

¹ Dipartimento di Scienza e Tecnologia del Farmaco, Università di Torino, Turin, Italy

² SRA Instruments SpA, Cernusco sul Naviglio, Milan, Italy

³ Agilent Technologies Italia SpA, Cernusco sul Naviglio, Milan, Italy

⁴ Agilent Technologies, Wilmington DE, USA

Abstract

A “second generation” design of differential flow modulation for comprehensive two-dimensional gas chromatography (GC×GC) was recently introduced [1,2]. This new configuration, based on the original device from Seeley et al. [3], adopts a reverse fill/flush injection dynamic instead of the forward flush/fill implemented in the “first generation” devices. Advantages include: (a) higher efficiency of band re-injection with improved ²D peak-widths and symmetry, (b) adjustable collection channel volume, (c) better handling of the overloading phenomenon without dramatic loss of peak-capacity/ resolution, and operation with hydrogen and helium carrier.

In this study, the first Capillary Flow Technology (CFT) reverse-inject differential flow modulator was implemented with different column configurations (lengths, diameters and stationary phase coupling) and detector combinations (Mass Spectrometry -MS and Flame Ionization Detection - FID) to evaluate its potential in the quantitative profiling of medium-to-highly complex essential oils and food volatiles fractions. In particular, a parallel dual-secondary column dual-detection configuration that has shown to improve the information potential of thermally modulated GC×GC platforms in terms of MS identification, reliability and accurate FID quantitation [4], was used here for the detailed characterization of essential oils (lavender, mint and vetiver) and food volatiles of interest for aroma blueprinting (extra-virgin olive oil and roasted hazelnuts).

Experimental results demonstrate that by changing column dimensions and configuration, the system flexibility can be greatly improved affording: (a) to extend the modulation period to match for ²D selectivity exploitation, (b) to operate with carrier gas linear velocities close-to-optimal values in both dimensions and (c) to handle overloaded peaks without dramatic losses in ²D peak-capacity and resolution.

References

[1] Griffith JF, Winniford WL, Sun K, Edam R, Luong JC. A reversed-flow differential flow modulator for comprehensive two-dimensional gas chromatography. *J Chromatogr A* 2012;1226:116-123

- [2] Xu S, Luong J, Aybar J, Miliazza A, Firor R. Characterization of a New Reversed Flow Modulator for GC×GC. Proceedings of the 38th ISCC and 11th GC×GC Symposium May 2014 Riva del Garda Italy
- [3] Bueno PA, Seeley JV. Flow-switching device for comprehensive two-dimensional gas chromatography. J. Chromatogr A 2004;1027:3-10.
- [4] Nicolotti L, Cordero C, Bressanello D, Cagliero C, Liberto E, Magagna F, Rubiolo P, Sgorbini B, Bicchi C. Parallel dual secondary column-dual detection: A further way of enhancing the informative potential of two-dimensional comprehensive gas chromatography. J. Chromatogr A. 2014;1360:264-274